

SYNERGY™ XD

Everolimus-Eluting Platinum Chromium
Coronary Stent System



 **Delivering More.**
Xtra Deliverability. Optimal Healing.

SYNERGY™ XD

Everolimus-Eluting Platinum Chromium
Coronary Stent System




Introducing: the next generation SYNERGY™ Bioabsorbable Polymer (BP) Stent. The SYNERGY XD Stent System is equipped with a delivery system designed to make optimal healing **even more deliverable**.

Innovative Delivery System for Xtra Deliverability

- Increased Pushability*
- Increased Trackability*

Fast Polymer Absorption for Optimal Healing

- Proven SYNERGY BP Stent
- Trusted Drug and Abluminal Polymer



* Based on bench testing of SYNERGY and SYNERGY XD BP Stents. Radial Guide Tracking and Radial Push Transmission; 3.00 mm stent systems tested, n=15 units. Bench testing performed by Boston Scientific. Results not necessarily indicative of clinical performance.

Xtra Deliverability

Innovative Delivery System



The **SYNERGY XD™ Stent System** combines added **trackability** and **pushability** with a **low-profile** stent platform to enable better navigation of challenging anatomy

Increased Trackability
Extended Lubricious
Coating*

Increased Pushability
Enhanced Laser-cut
Hypotube Technology*



Designed for Flexibility
Short Red Tip and 2-Connector Stent Design




Low Crossing Profile
Thin, Rounded Strut Design



Exceptional Overexpansion
Up to 5.75 mm†

Best-in-Class Visibility
Proprietary PtCr Alloy




* Based on bench testing of SYNERGY and SYNERGY XD BP Stents. Radial Guide Tracking and Radial Push Transmission; 3.00 mm stent systems tested, n=15 units.
Bench testing performed by Boston Scientific. Results not necessarily indicative of clinical performance.
† Based on the 4.00 - 5.00 mm stent models labeled post-dilatation limits. Please refer to the directions for use for more detail on safety information.
Please refer to the directions for use for the SYNERGY XD BP Stent.

Optimal Healing

Fast Absorbing Bioabsorbable Polymer Advantage

Unlike permanent polymer DES technology, the SYNERGY XD Bioabsorbable Polymer is gone shortly after the drug is completely eluted¹


Optimized Everolimus Drug Delivery
Synchronous drug elution & polymer absorption



Fast Polymer Absorption
The 4-month absorption time is up to 6x faster than other BP Stent[‡]



Abluminal Polymer Application
Suppresses neointimal growth while promoting early healing



Excellent Safety & Long Term Outcomes

The SYNERGY BP Stent has been studied in over 35,000 patients across various patient and lesion complexities



The HBR indication is supported by the data from the EVOLVE Short DAPT Trial.

0.2% ST

After patients stopped DAPT at 3-months through 15-months
EVOLVE Short DAPT Trial[§]

ZERO ST
After patients stopped DAPT at 1-month through 12-months
SENIOR Trial[§]

RANKED #1
Lowest relative risk of Def/Prob ST
Kang Network Meta-Analysis[§]

Excellent Results In
COMPLEX PATIENTS
Lowest ST rate in statistically more complex patients
Real-World SCAAR Registry[§]

§ The HBR indication excludes the SYNERGY XD 48 mm stent

** EVOLVE Short DAPT is a prospective, multicenter, single-arm trial defining the safety of 3-month DAPT in subjects at high risk for bleeding undergoing PCI with the SYNERGY BP Stent. Approximately 74% of patients enrolled discontinued DAPT at 3-months. N = 1,396 (patients with respective event or sufficient follow-up). Co-primary endpoints: ARC Def/Prob ST and Death/MI from 3-15 months.

3. Varenne O et al. Drug-eluting stents in elderly patients with coronary artery disease (SENIOR): a randomised single-blind trial. *Lancet*. 2018;391(10115):45-50. doi:10.1016/S0140-6736(17)32713-7.

4. Kang S et al. *J Am Coll Cardiol Interv*. doi:10.1016/j.jcin.2016.03.038.

5. EVOLVE Short DAPT Subgroup Analyses. Presented by Robert Stoler, MD at ACC 2020.

6. Sarno G et al. *Cathet Cardiovasc Intervent*. doi:10.1002/ccd.27030.

‡ Compared to Orsiro BP Stent.

1. Wilson GJ et al. *Catheter Cardiovasc Interv*.

2. Soucy N, Feygin J et al. *EuroIntervention*. 2010 Nov;6(5):630-7.

ø (mm)	STENT LENGTH (mm)										Expansion Limit
	8	12	16	20	24	28	32	38	48		
2.25	H7493941808220	H7493941812220	H7493941816220	H7493941820220	H7493941824220	H7493941828220	H7493941832220	H7493941838220			3.50
2.50	H7493941808250	H7493941812250	H7493941816250	H7493941820250	H7493941824250	H7493941828250	H7493941832250	H7493941838250	H7493941848250		3.50
2.75	H7493941808270	H7493941812270	H7493941816270	H7493941820270	H7493941824270	H7493941828270	H7493941832270	H7493941838270	H7493941848270		3.50
3.00	H7493941808300	H7493941812300	H7493941816300	H7493941820300	H7493941824300	H7493941828300	H7493941832300	H7493941838300	H7493941848300		4.25
3.50	H7493941808350	H7493941812350	H7493941816350	H7493941820350	H7493941824350	H7493941828350	H7493941832350	H7493941838350	H7493941848350		4.25
4.00	H7493941808400	H7493941812400	H7493941816400	H7493941820400	H7493941824400	H7493941828400	H7493941832400	H7493941838400	H7493941848400		5.75
4.50	H7493941812450	H7493941816450	H7493941820450	H7493941824450	H7493941828450	H7493941832450					5.75
5.00	H7493941821500	H7493941816500	H7493941820500	H7493941824500	H7493941828500	H7493941832500					5.75

SYNERGY™ XD Monorail Coronary Stent System

INDICATIONS FOR USE: The SYNERGY™ XD Everolimus-Eluting Platinum Chromium Coronary Stent System is indicated for improving luminal diameter in patients, including those with diabetes mellitus, with symptomatic heart disease, stable angina, unstable angina, non-ST elevation MI or documented silent ischemia due to atherosclerotic lesions in native coronary arteries <2.25 mm in diameter in lesions <44 mm in length and for high bleeding risk patients with coronary arteries <2.25 mm to <5.00 mm in diameter in lesions <34 mm in length. **CONTRAINDICATIONS:** Use of the SYNERGY XD Everolimus-Eluting Platinum Chromium Stent System is contraindicated in patients with known hypersensitivity to: • 316L stainless steel, platinum, chromium, nickel or molybdenum • Everolimus or structurally-related compounds • The polymer or their individual components (see Section 2.4.2 Polymer Carrier in the eUfU) Coronary Artery Stenting is contraindicated for use in: • Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or delivery device. • Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy (see Sections 6.2 Pre- and Post-Procedure Antiplatelet Regimen for more information in the eUfU). **WARNINGS:** • To maintain sterility, the package should not be opened or damaged prior to use. • The use of this product carries the risks associated with coronary artery stenting, including stent thrombosis, vascular complications, and/or bleeding events. • This product should not be used in patients who are not likely to comply with recommended antiplatelet therapy. **GENERAL PRECAUTIONS:** • Careful consideration should be given to the risks and benefits of use in patients with history of severe reaction to contrast agents. • Stent thrombosis is a rare event and is frequently associated with myocardial infarction (MI) or death. In the clinical trials analyzed to date, differences in the incidence of stent thrombosis have not been associated with an increased risk of cardiac death, MI, or device mortality. • Patients with a history of bleeding should be carefully monitored for bleeding events. • The following outcomes may be observed during the EVOLVE clinical trials: • Compared to use with the specified Indications for Use, the use of DES in patients with lesions outside of the labeled indications may have an increased risk of adverse events, including stent thrombosis, stent embolization, MI or death. When treating such patients, physicians should be aware of this increased risk and consider available data and the limitations of such data. • Orally-administered everolimus combined with cyclosporine is associated with increased serum cholesterol and triglyceride levels. **Pre- and Post-Procedure Antiplatelet Regimen** The optimal duration of antiplatelet therapy, specifically P2Y12 inhibitor therapy is unknown and DES thrombosis may still occur despite continued therapy beyond current professional society guidelines. **Oral Antiplatelet Therapy:** Continuation of combination treatment with aspirin and a P2Y12 inhibitor after PCI appears to reduce major adverse cardiac events. On the basis of randomized clinical trials the 2016 ACC/AHA guidelines recommend aspirin 81 mg daily should be given indefinitely after PCI in patients who are not at high risk of bleeding. A P2Y12 inhibitor should be given daily for at least 6 months in stable ischemic heart disease patients and for at least 12 months in acute coronary syndrome (ACS) patients. Full guidelines are provided at the following website: <http://www.onlineacc.org>. Based upon the results of the EVOLVE Short-DAPT Study the SYNERGY XD stent can be safely used in conjunction with shortened DAPT in patients at high risk for bleeding. In the EVOLVE Short-DAPT Study, high bleeding risks were defined as meeting one or more of the following criteria: ≥75 years of age and, in the opinion of the investigator, the risk of bleeding associated with >3 months of DAPT outweighs the benefit; need for chronic or lifelong anticoagulation; history of major bleeding (severe IFRG) or Stent thrombosis; or conditional contraindications up to 12 months of the index procedure; history of stroke (ischemic or hemorrhagic); renal insufficiency (creatinine >2.0 mg/dL) or failure (dialysis dependent); platelet count <100,000/μL. It is very important that the patient is compliant with the post-procedure antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, MI or death. Prior to PCI, if a surgical or dental procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventional cardiologist and patient should carefully consider whether a DES and its associated recommended antiplatelet therapy is the appropriate PCI choice. **Pediatric Use:** The safety and effectiveness of the SYNERGY Stent System in pediatric patients have not been established. **Lesion/Vessel Characteristics:** The safety and effectiveness of the SYNERGY Stent System have not been established in the cerebral, carotid, or peripheral vasculature or in the following patient populations: • Patients with vessel thrombus at the lesion site. • Patients with coronary artery ostial/lesion complex diameters <2.25 or >5.00 mm. • Patients with coronary artery lesions longer than 44mm or requiring more than one SYNERGY Stent. • Patients with lesions located in saphenous vein grafts, in the left main coronary artery, aortic valve/corona, or complex bifurcation (i.e. bifurcation lesion requiring treatment with more than one stent). • Patients with diffuse disease or reduced blood flow distal to the identified lesions. • Patients with a recent acute ST elevation myocardial infarction where there is evidence of thrombosis or poor flow. • Patients with in-stent restenosis. • Patients with a coronary total occlusion. • Patients with 3 vessel disease. **Magnetic Resonance Imaging (MRI) Safety Information:** Non-clinical data have demonstrated safety of the SYNERGY XD Stent System in conditional contraindications up to 94 mm. A patient with this device can be safely scanned in a Magnetic Resonance system meeting the following conditions: • Static magnetic field of 3.0 and 1.5 Tesla only • Maximum spatial gradient magnetic field of 2000 gauss/cm (23.7 mT/m) • Maximum Magnetic Resonance system reported, whole body averaged specific absorption rate (SAR) of ≤2 W/kg (Normal Operating Mode) Under the scan conditions defined above, the SYNERGY XD Stent is expected to produce a maximum temperature rise of < 0.5 °C less after 15 minutes of continuous scanning. MRI image quality may be compromised if the area of interest is within the lumen or relatively near the stent. Therefore, it may be necessary to optimize MRI imaging parameters for the presence of the stent. The image artifact extends approximately 1 cm from the stent when scanned in non-clinical MRI testing specified in ASTM F2119-07. The artifact does obscure the device lumen. Image artifact was minimized using the spin echo sequence versus gradient echo. **POTENTIAL ADVERSE EVENTS:** Potential adverse events (in alphabetical order) which may be associated with the use of a coronary stent in native coronary arteries include but are not limited to: • Abrupt stent dislodge • Allergic reaction to anti-coagulant and/or antiplatelet therapy, contrast medium, or stent materials • Angina • Arrhythmias, including ventricular fibrillation, ventricular tachycardia and heart block • Cardiogenic shock/pulmonary edema • Death • Embolization (air, tissue or thrombotic material or material from device) (s) used in the procedure), including stent embolization or migration • Heart failure • Hemorrhage, which may require transfusion, including bleeding and hematoma • Hypertension/hypertensive • Infection, local or systemic, including fever and pyrogen reaction • Myocardial ischemia or infarction • Pain, chest or access site • Pericardial effusion or cardiac tamponade • Renal insufficiency or failure • Respiratory failure • Restenosis or aneurysm of stented segment • Stent deformation, collapse or fracture • Stent thrombosis/occlusion • Stroke/cerebrovascular accident/transient ischemic attack • Vessel trauma requiring surgical repair or reintervention; including coronary, femoral or radial artery spasm, dissection; occlusion; perforation, rupture, or pseudoaneurysm • Zoster™, the oral formulation of everolimus developed by Novartis Pharmaceuticals Corporation, has been evaluated in clinical trials and is approved in the United States for the prevention of organ rejection in adult kidney transplant recipients at the dose of 15 mg/day. Outside the U.S., Zoster is sold under the brand name, Certican™, in more than 70 countries. Everolimus is also approved in the United States under the name of Afinitor™ for patients with advanced renal cell carcinoma (cancer), after failure of treatment with sunitinib or sorafenib, at doses of 5 to 20 mg/day when taken by mouth. The following list includes the known risks of everolimus at the oral doses listed above. The amount of drug that circulates in the bloodstream following implantation of a SYNERGY™ Stent is several folds lower than that obtained with oral doses (15 mg to 20 mg/day, see Section 7.2, Pharmacokinetics of the eUfU). • Abdominal pain • Anemia • Angioedema • Anorexia • Asthenia • Constipation • Cough • Delayed wound healing • Fluid accumulation • Diarrhea • Dyslipidemia (including hyperlipidemia and hypercholesterolemia) • Dyspnea • Dyspnea • Dyspnea • Dry skin • Edema (peripheral) • Epistaxis • Fatigue • Headache • Hematuria • Hyperkalemia • Hyperglycemia (may include new onset of diabetes) • Hypokalemia • Hyperlipidemia • Hypertension • Hypokalemia • Hypomagnesaemia • Hypophosphatemia • Increased serum creatinine • Infections and serious infections: bacterial, viral, fungal, and protozoal infections (may include herpes virus infection, polyoma virus infection which may be associated with BK virus associated nephropathy, and/or other opportunistic infections) • Ischemia • Interaction with strong inhibitors and inducers of CYP3A4 • Leukopenia • Lymphoma and other malignancies (including skin cancer) • Male infertility (azoospermia and/or oligospermia) • Mucosal inflammation (including oral ulceration and oral mucositis) • Nausea • Neurotension • Non-infectious pneumonitis • Pain; extremity; injection site and procedural, back, chest, musculoskeletal • Proteinuria • Pruritus • Rash • Stomatitis • Thrombocytopenia • Thrombotic microangiopathy (TMA)/Thrombotic thrombocytopenia purpura (TTP)/Hemolytic uremic syndrome (HUS) • Tremor • Upper respiratory tract infection • Urinary tract infection • Vomiting • Skin vaccines should be avoided and dose contact with those that have had live vaccines should be avoided. Fetal harm can occur when administered to a pregnant woman. There may be other potential adverse events that are unforeseen at this time. **CAUTION:** Federal law (FDA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete Instructions for Use for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

9257883B

Boston Scientific
Advancing science for life™

Interventional Cardiology
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

*To order product or for more information
contact customer service at 1.888.272.1001.*

© 2020 Boston Scientific Corporation
or its affiliates. All rights reserved.

IC-751701-AB