

Peripheral Interventions

Product Catalogue - Produktkatalog
Catalogue Produits - Catálogo de Productos

**Boston
Scientific**
Advancing science for life™



1. Guidewires

- 1.1 0.014" (0.36mm) Guidewires**
 - Journey™
 - Thruway™
 - Platinum™ Plus
 - V-14™ ControlWire®
 - Victory™
- 1.2 0.018" (0.46mm) Guidewires**
 - Thruway™
 - Platinum™ Plus
 - V-18™ ControlWire™
 - Victory™
 - ZIPWire™
- 1.3 0.035" (0.89 mm) +/- Guidewires**
 - PTFE coated - Starter™
 - Platinum™ Plus
 - ZIPWire™
- 1.4 Specialty Guidewires**
 - Magic Torque™
 - Amplatz Super Stiff™
 - Back-up Meier™

2. Balloon Dilatation Catheters

- 2.1 0.014" (0.36mm) Guidewire Compatible**
 - Coyote™
 - Coyote™ ES
 - Ultra-Soft™ SV
- 2.2 0.018" (0.46mm) Guidewire Compatible**
 - Sterling™
 - Sterling SL™
 - Symmetry™
 - Ultra-Soft™ SV
- 2.3 0.035" (0.89 mm) Guidewire Compatible**
 - Mustang™
 - Wanda™
 - XXL™

3. Specialty Balloons

Peripheral Cutting Balloon™
Equalizer™ Occlusion Balloon™

4. CTO Devices

OffRoad™ Re-Entry Catheter System
TruePath™ CTO Device

5. Catheters

- 5.1 Diagnostic Catheter - Imager II™**
- 5.2 Guiding Catheter - Mach1™**
- 5.3 Rubicon™ Support Catheter**

6. Accessories

Medikit Supersheath™ Introducer Sheath
Medikit Supersheath™ RO Introducer Sheath
Medikit Supersheath™ XL Introducer Sheath
Encore™ 26 Inflation Device
Gateway™ Plus Y-Adaptor
Vascular Dilators
Femoral Entry Needles
Torque Visers
Connecting Tubes
Floswitch™ HP High Pressure Flow Control Device

7. Stents

- 7.1 Balloon Expandable Stenting**
 - Express™ Vascular LD
 - Express™ Vascular SD
 - PROMUS Element™ Plus DES BTK
- 7.2 Self-Expandable Stenting**
 - WALLSTENT-Uni™ Endoprosthesis
 - WALLSTENT™ RP
 - Epic™ Nitinol Vascular
 - Innova™
 - 7F (2.37 mm) Placehit WALLSTENT™
 - WALLGRAFT™ Covered Stent
- 7.3 Carotid Stenting**
 - Adapt™ Carotid Stent
 - Carotid WALLSTENT™ Monorail™ Endoprosthesis

8. Embolic Protection

- 8.1 FilterWire EZ™ Embolic Protection System**

9. Peripheral Embolisations

- 9.1 Guidewires & Microcatheters**
 - FATHOM™-16 Steerable Guidewires
 - Transend™ Steerable Guidewires
 - Renegade™ Fiber Braided Microcatheters
 - Renegade™ HI-FLO™ Fathom™
 - FasTracker™ 325 Microcatheters
- 9.2 Embolic Agents & Occlusion Balloons**
 - Interlock™ Fibered IDC™ Occlusion System Family
 - IDC™ Interlocking Detachable Coils
 - Complex Helical -18 Fibered Platinum Coils
 - Pushable Coils
 - Occlusion Balloons
 - Contour™ SE Microspheres
 - Contour™ Embolisation Particles

10. Non-Vascular Interventions

- 10.1 All-Purpose Drainage Catheters**
 - Flexima™
 - Percuflex™
- 10.2 vanSonnenberg Drainage Catheters**
 - Percuflex™
 - Flexima™
- 10.3 Biliary Drainage Catheters**
 - Flexima™
 - Carey-Coons™
 - Amplatz Anchor™
- 10.4 Urinary Drainage Catheters**
 - Flexima™
 - Percuflex™
- 10.5 Ureteral Drainage Catheters**
 - Flexima™
 - Percuflex™
 - Amplatz Anchor™
 - Malecot Nephrostomy Catheters
- 10.6 Introducers and Accessories**

11. Radiofrequency Ablation

- 11.1 RF 3000™ Generator**
- 11.2 Needle Electrodes and Accessories**
 - LeVeen™ Needle Electrode
 - LeVeen™ CoAccess Electrode System
 - LeVeen™ SuperSlim™ Needle Electrode
 - Soloist™ Single Needle Electrode

12. Miscellaneous

- 12.1 Thrombolysis**
 - Katzen™ Infusion Wire

1. Führungsdrähte

- 1.1 **0.014" (0.36mm) Führungsdrähte**
Journey™
Thruway™
Platinum™ Plus
V-14™ ControlWire®
Victory™
- 1.2 **0.018" (0.46mm) Führungsdrähte**
Thruway™
Platinum™ Plus
V-18™ ControlWire™
Victory™
ZIPWire™
- 1.3 **0.035" (0.89 mm) +/- Führungsdrähte**
PTFE coated - Starter™
Platinum™ Plus
ZIPWire™
- 1.4 **Spezial-Führungsdrähte**
Magic Torque™
Amplatz Super Stiff™
Back-up Meier™

2. Ballondilatationskatheter

- 2.1 **0.014" (0.36mm) Führungsdraht-Kompatibilität**
Coyote™
Coyote™ ES
Ultra-Soft™ SV
- 2.2 **0.018" (0.46mm) Führungsdraht-Kompatibilität**
Sterling™
Sterling SL™
Symmetry™
Ultra-Soft™ SV
- 2.3 **0.035" (0.89 mm) Führungsdraht-Kompatibilität**
Mustang™
Wanda™
XXL™

3. Spezialballons

Peripheral Cutting Balloon™
Equalizer™ Occlusions-Ballon

4. CTO Devices

OffRoad™ Re-Entry Katheter-System
TruePath™ CTO Device

5. Katheter

- 5.1 Angiographie-Katheter - Imager II™
- 5.2 Führungskatheter - Mach1™
- 5.3 Rubicon™ Supportkatheter

6. Zubehör

Medikit Supersheath™ Einführschleuse
Medikit Supersheath™ RO Einführschleuse
Medikit Supersheath™ XL Einführschleuse
Inflatoren - Encore™ 26
Y-Adaptor - Gateway™ Plus
Gefäß-Dilatatoren
Femorale Einführnadeln
Torque Visers
Connecting Tubes
Floswitch™ HP Hochdruck-Durchflussregler

7. Stents

- 7.1 **Balloon Expandable Stenting**
Express™ Vascular LD
Express™ Vascular SD
PROMUS Element™ Plus DES BTK
- 7.2 **Self-Expandable Stenting**
WALLSTENT-Uni™ Endoprothesis
WALLSTENT™ RP
Epic™ Nitinol Vascular
Innova™
7F (2.37 mm) Placehit WALLSTENT™
WALLGRAFT™ Covered Stent
- 7.3 **Carotid Stenting**
Adapt™ Carotid Stent
Carotid WALLSTENT™ Monorail™ Endoprosthesis

8. Embolieschutz

- 8.1 **FilterWire EZ™ Embolieprotektionssystem**

9. Periphere Embolisationen

- 9.1 **Führungsdrähte & Mikrokatheter**
FATHOM™-16 Steerable Guidewires
Transend™ Steerable Guidewires
Renegade™ Fiber Braided Microcatheters
Renegade™ HI-FLO™ Fathom™
FasTracker™ 325 Microcatheters
- 9.2 **Embolisate & Okklusions-Ballons**
Interlock™ Fibered IDC™ Occlusion System Family
IDC™ Interlocking Detachable Coils
Complex Helical -18 Fibered Platinum Coils
Pushable Coils
Occlusion Balloons
Contour™ SE Microspheres
Contour™ Embolisation Particles

10. Nicht-Vaskuläre Interventionen

- 10.1 **Universal-Drainage-Katheter**
Flexima™
Percuflex™
- 10.2 **vanSonnenberg-Drainagekatheter**
Percuflex™
Flexima™
- 10.3 **Gallendrainage-Katheter**
Flexima™
Carey-Coons™
Amplatz Anchor™
- 10.4 **Harndrainage-Katheter**
Flexima™
Percuflex™
- 10.5 **Ureteral-Drainage-Katheter**
Flexima™
Percuflex™
Amplatz Anchor™
Malecot Nephrostomy Catheters
- 10.6 **Einführschleusen und Zubehör**

11. Radiofrequenz-Ablation

- 11.1 **RF 3000™ Generator**
- 11.2 **Nadelelektroden und Zubehör**
LeVeen™ Nadelelektrode
LeVeen™ CoAccess Elektrodensystem
LeVeen™ SuperSlim™ Nadelelektrode
Soloist™ Single Needle Elektrode

12. Sonstiges

- 12.1 **Thrombolyse**
Katzen™ Infusionsdraht

1. Guides

- 1.1 **0.014" (0.36mm) Guides**
Journey™
Thruway™
Platinum™ Plus
V-14™ ControlWire®
Victory™
- 1.2 **0.018" (0.46mm) Guides**
Thruway™
Platinum™ Plus
V-18™ ControlWire™
Victory™
ZIPWire™
- 1.3 **0.035" (0.89 mm) +/- Guides**
PTFE coated - Starter™
Platinum™ Plus
ZIPWire™
- 1.4 **Guides Spéciaux**
Magic Torque™
Amplatz Super Stiff™
Back-up Meier™

2. Cathéters de dilatation à ballonnet

- 2.1 **0.014" (0.36mm) Guide Compatible**
Coyote™
Coyote™ ES
Ultra-Soft™ SV
- 2.2 **0.018" (0.46mm) Guide Compatible**
Sterling™
Sterling SL™
Symmetry™
Ultra-Soft™ SV
- 2.3 **0.035" (0.89 mm) Guide Compatible**
Mustang™
Wanda™
XXL™

3. Ballonnets Spéciaux

Peripheral Cutting Balloon™
Equalizer™ Occlusion Balloon™

4. CTO Devices

OffRoad™ Re-Entry Catheter System
TruePath™ CTO Device

5. Cathéters

- 5.1 Diagnostic Catheter - Imager II™
- 5.2 Guiding Catheter - Mach1™
- 5.3 Rubicon™ Support Catheter

6. Accessoires

Medikit Supersheath™ Introducteur
Medikit Supersheath™ RO Introducteur
Medikit Supersheath™ XL Introducteur
Inflateur - Encore™ 26
Adaptateur en Y Gateway™ Plus
Dilatateurs vasculaires
Aiguilles de ponction fémorale
Torqueurs
Raccords pour injection de contraste
Robinet 1 voie Floswitch™ HP

7. Stents

- 7.1 **Balloon Expandable Stenting**
Express™ Vascular LD
Express™ Vascular SD
PROMUS Element™ Plus DES BTK
- 7.2 **Self-Expandable Stenting**
WALLSTENT-Uni™ Endoprosthesis
WALLSTENT™ RP
Epic™ Nitinol Vascular
Innova™
7F (2.37 mm) Placehit WALLSTENT™
WALLGRAFT™ Covered Stent
- 7.3 **Carotid Stenting**
Adapt™ Carotid Stent
Carotid WALLSTENT™ Monorail™ Endoprosthesis

8. Protection Embolique

- 8.1 **FilterWire EZ™ Embolic Protection System**

9. Embolisations Périphériques

- 9.1 **Guides et microcathéters**
FATHOM™-16 Steerable Guidewires
Transend™ Steerable Guidewires
Renegade™ Fiber Braided Microcatheters
Renegade™ HI-FLO™ Fathom™
FasTracker™ 325 Microcatheters
- 9.2 **Agents emboliques et ballons d'occlusion**
Interlock™ Fibered IDC™ Occlusion System Family
IDC™ Interlocking Detachable Coils
Complex Helical -18 Fibered Platinum Coils
Pushable Coils
Occlusion Balloons
Contour™ SE Microspheres
Contour™ Embolisation Particles

10. Interventions non Vasculaires

- 10.1 **Cathéters de drainage universels**
Flexima™
Percuflex™
- 10.2 **Cathéters de drainage vanSonnenberg**
Percuflex™
Flexima™
- 10.3 **Cathéters de drainage biliaire**
Flexima™
Carey-Coons™
Amplatz Anchor™
- 10.4 **Cathéters de drainage urinaire**
Flexima™
Percuflex™
- 10.5 **Cathéters de drainage urétéral**
Flexima™
Percuflex™
Amplatz Anchor™
Malecot Nephrostomy Catheters
- 10.6 **Introducteurs et Accessoires**

11. Ablation par Radiofréquence

- 11.1 **Générateur RF 3000™**
- 11.2 **Électrodes et Accessoires**
Électrode LeVeen™
Système d'électrode LeVeen™ CoAccess
Électrode LeVeen™ SuperSlim™
Électrode à aiguille droite Soloist™

12. Divers

- 12.1 **Thrombolyse**
Guide d'infusion Katzen™

1. Guías

- 1.1 **0.014" (0.36mm) Guías**
Journey™
Thruway™
Platinum™ Plus
V-14™ ControlWire®
Victory™
- 1.2 **0.018" (0.46mm) Guías**
Thruway™
Platinum™ Plus
V-18™ ControlWire™
Victory™
ZIPWire™
- 1.3 **0.035" (0.89 mm) +/- Guías**
PTFE coated - Starter™
Platinum™ Plus
ZIPWire™
- 1.4 **Specialty Guidewires**
Magic Torque™
Amplatz Super Stiff™
Back-up Meier™

2. Catéteres balón de dilatación

- 2.1 **0.014" (0.36mm) Guidewire Compatible**
Coyote™
Coyote™ ES
Ultra-Soft™ SV
- 2.2 **0.018" (0.46mm) Guidewire Compatible**
Sterling™
Sterling SL™
Symmetry™
Ultra-Soft™ SV
- 2.3 **0.035" (0.89 mm) Guidewire Compatible**
Mustang™
Wanda™
XXL™

3. Balones Especiales

Cutting balón periférico
Balón de oclusión Equalizer™

4. CTO Devices

OffRoad™ Re-Entry Catheter System
TruePath™ CTO Device

5. Catéteres

- 5.1 **Diagnostic Catheter - Imager II™**
- 5.2 **Guiding Catheter - Mach1™**
- 5.3 **Rubicon™ Support Catheter**

6. Accesorios

Medikit Supersheath™ Introducer Sheath
Medikit Supersheath™ RO Introducer Sheath
Medikit Supersheath™ XL Introducer Sheath
Dispositivo de inflado - Encore™ 26
Adaptador Y - Gateway™ Plus
Dilatadores Vasculares
Agujas de acceso femoral
Torque
Tubos de conexión
Dispositivos de control de flujo de alta presión
Floswitch™ HP

7. Stents

- 7.1 **Balloon Expandable Stenting**
Express™ Vascular LD
Express™ Vascular SD
PROMUS Element™ Plus DES BTK
- 7.2 **Self-Expandable Stenting**
WALLSTENT-Uni™ Endoprosthesis
WALLSTENT™ RP
Epic™ Nitinol Vascular
Innova™
7F (2.37 mm) Placehit WALLSTENT™
WALLGRAFT™ Covered Stent
- 7.3 **Carotid Stenting**
Adapt™ Carotid Stent
Carotid WALLSTENT™ Monorail™ Endoprosthesis

8. Protección Embólica

- 8.1 **Embólico FilterWire EZ™**

9. Embolización Periférica

- 9.1 **Guías y microcatéteres**
FATHOM™-16 Steerable Guidewires
Transend™ Steerable Guidewires
Renegade™ Fiber Braided Microcatheters
Renegade™ HI-FLO™ Fathom™
FasTracker™ 325 Microcatheters
- 9.2 **Agentes embólicos y balones de oclusión**
Interlock™ Fibered IDC™ Occlusion System Family
IDC™ Interlocking Detachable Coils
Complex Helical -18 Fibered Platinum Coils
Pushable Coils
Occlusion Balloons
Contour™ SE Microspheres
Contour™ Embolisation Particles

10. Intervenciones no Vasculares

- 10.1 **Catéteres de drenaje multiusos**
Flexima™
Percuflex™
- 10.2 **Catéteres de drenaje vanSonnenberg**
Percuflex™
Flexima™
- 10.3 **Catéteres de drenaje biliar**
Flexima™
Carey-Coons™
Amplatz Anchor™
- 10.4 **Catéteres de drenaje urinario**
Flexima™
Percuflex™
- 10.5 **Catéteres de drenaje uretral**
Flexima™
Percuflex™
Amplatz Anchor™
Malecot Nephrostomy Catheters
- 10.6 **Introdutores y accesorios**

11. Ablación por Radiofrecuencia

- 11.1 **RF 3000™ Generator**
- 11.2 **Electrodos de agujas y accesorios**
Electrodo de agujas LeVeen™
Sistema de electrodo LeVeen™ CoAccess
Electrodo de agujas LeVeen™ SuperSlim™
Electrodo de aguja única Soloist™

12. Otros

- 12.1 **Trombólisis**
Guía de infusión Katzen™

1 Guidewires

1.1 0.014" (0.36mm) Guidewires

Journey™

Thruway™

Platinum™ Plus

V-14™ ControlWire™

Victory™




1.1 0.014" (0.36mm) Guidewires

Journey™

Peripheral Guidewire / Peripherer Führungsdraht / Guide périphérique / Guía periférica

- Designed to exceed conventional guidewire performance by providing far superior torque control without compromising flexibility
- Bessere Eigenschaften als konventionelle Führungsdrähte durch deutlich bessere Drehmomentkontrolle bei gleichzeitig guter Flexibilität
- Conçu pour surpasser les performances des guides conventionnels en offrant un contrôle du torque bien supérieur sans compromis sur la flexibilité
- Diseñadas para superar el rendimiento convencional de la guía proporcionando un control de torsión superior sin comprometer la flexibilidad




Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Tip Shape	Tip Length (cm)
M001391260	0.014	0.36	185	Straight	2
M001391270	0.014	0.36	185	Angled	2
M001391280	0.014	0.36	300	Straight	2
M001391290	0.014	0.36	300	Angled	2

Thruway™

Peripheral Guidewire / Peripherer Führungsdraht / Guide périphérique / Guía periférica

- Indicated for diagnostic and interventional intravascular procedures, including renal intervention; combining precision, strength and balance; featuring radiopaque distal markers
- Indiqué pour les procédures diagnostiques et les interventions intravasculaires, dont les procédures rénales. Associe précision, résistance et équilibre ; avec marqueurs radio-opaques distaux
- Indiziert für diagnostische und interventionelle intravaskuläre Eingriffe, beispielsweise an der Niere; kombiniert Präzision und Festigkeit mit guter Sichtbarkeit durch distale Röntgenmarker
- Indicada para procedimientos de diagnóstico e intervención intravascular, incluidas intervenciones renales. Combina precisión, resistencia y equilibrio, y ofrece marcadores distales radiopacos



Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Tip Taper Length	Tip Shape	Distal Coil (cm)
M001492871	0.014	0.36	130	Short	Straight	3
M001492881	0.014	0.36	130	Short	J	3
M001492891	0.014	0.36	130	Long	Straight	3
M001492901	0.014	0.36	130	Long	J	3
M001492911	0.014	0.36	130	Extra Long	Straight	5
M001492921	0.014	0.36	190	Short	Straight	3
M001492931	0.014	0.36	190	Short	J	3
M001492941	0.014	0.36	190	Long	Straight	3
M001492951	0.014	0.36	190	Long	J	3
M001492961	0.014	0.36	190	Extra Long	Straight	5
M001492971	0.014	0.36	300	Short	Straight	3
M001492981	0.014	0.36	300	Long	Straight	3
M001492991	0.014	0.36	300	Extra Long	Straight	5


Packaged 1/box

1.1 0.014" (0.36mm) Guidewires

Platinum™ Plus

Guidewire ST / Führungsdraht ST / Guide ST / Guía ST

- For access in anatomy with short distal purchase. Distal platinum coils for good visualization
- Für den Zugang zu Anatomien mit kurzem, distalen Ende. Distale Platinspitze für gute Sichtbarkeit.
- Pour accéder aux anatomies à courtes prises distales. Coils distaux en platine pour une bonne visualisation
- Para intervenciones en anatomías con base distal corta. Coils distales de platino para una buena visualización



Platinum Plus™ 0.014" (0.36 mm) with Glidex™ Coating

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Radiopaque Tip Length (cm)	Floppy Tip Length (cm)
M001467400	0.014	0.36	180	5	3
M001467120	0.014	0.36	300	3	3

Packaged 1/box

V-14™ ControlWire®

Polymer sleeve with ICE™ Hydrophilic Coating for excellent tracking / Polymer-Ummantelung mit hydrophiler ICE™ Beschichtung für exzellente Steuerung / Polymère avec revêtement hydrophile ICE™ pour une navigabilité excellente / Revestimiento de polímero con recubrimiento hidrofílico ICE™ para una navegación excelente

- V-14™ ControlWire® Guidewire offers the choice of two tip load options with excellent rail support
- V-14™ ControlWire® Guidewire bietet zwei Spitzen zur Auswahl und exzellente Führungsunterstützung
- V-14™ ControlWire® Guidewire permet de choisir entre deux extrémités de rigidité différentes
- V-14™ ControlWire® Guidewire ofrece dos opciones de cargas de punta con un excelente soporte de riel

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Tip Taper Length	Tip Shape
H74939216718210	0.014	0.36	182	Short	Straight
H74939216718220	0.014	0.36	182	Short	Angled
H74939216730010	0.014	0.36	300	Short	Straight
H74939216730020	0.014	0.36	300	Short	Angled
H74939216918210	0.014	0.36	182	Long	Straight
H74939216918220	0.014	0.36	182	Long	Angled
H74939216930010	0.014	0.36	300	Long	Straight
H74939216930020	0.014	0.36	300	Long	Angled

Packaged 1/box

Victory™

Peripheral Guidewire / Peripherer Führungsdraht / Guide périphérique / Guía periférica

- Victory™ is designed to break through and cross resistant lesions, offering increasing tip loads from 12g to 30g for control and a range of pushability
- Le Guide Périphérique Victory™ est conçu pour franchir et traverser des lésions résistantes ; il offre un bon contrôle et garantit de multiples niveaux de pression
- Victory™ ermöglicht leichtes Passieren resistenter Läsionen und bietet zudem mit 4 unterschiedlichen Spitzengewichten von 12g bis zu 30g stets höchste Kontrolle und präzises Vorschieben
- Las guías Victory™ están diseñadas para avanzar y cruzar lesiones resistentes mediante puntas que van de los 12g a los 30g para proporcionar el control y la capacidad de empuje necesarios

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Radiopaque Tip Length (cm)	Tip Load (g)
H74939230195120	0.014	0.36	195	2.5	12
H74939230300120	0.014	0.36	300	2.5	12
H74939230195180	0.014	0.36	195	2.5	18
H74939230300180	0.014	0.36	300	2.5	18
H74939231195250	0.014	0.36	195	2	25
H74939231300250	0.014	0.36	300	2	25
H74939231195300	0.014	0.36	195	2	30
H74939231300300	0.014	0.36	300	2	30

Ordering Information - 8
Peripheral Interventions

Boston
Scientific



1 Guidewires

1.2 0.018" (0.46mm) Guidewires


Thruway™

Platinum™ Plus

V-18™ ControlWire™

Victory™

ZIPWire™



1.2 0.018" (0.46mm) Guidewires

Thruway™

Peripheral Guidewire / Peripherer Führungsdraht / Guide périphérique / Guía periférica

- Indicated for diagnostic and interventional intravascular procedures, including renal intervention; combining precision, strength and balance; featuring radiopaque distal markers
- Indiziert für diagnostische und interventionelle intravaskuläre Eingriffe, beispielsweise an der Niere; kombiniert Präzision und Festigkeit mit guter Sichtbarkeit durch distale Röntgenmarker
- Indiqué pour les procédures diagnostiques et les interventions intravasculaires, dont les procédures rénales. Associe précision, résistance et équilibre ; avec marqueurs radio-opaques distaux
- Indicada para procedimientos de diagnóstico e intervención intravascular, incluidas intervenciones renales. Combina precisión, resistencia y equilibrio, y ofrece marcadores distales radiopacos


Order Number	Diameter* (in) (mm)		Length (cm)	Tip Length	Tip Shape
M001492751	0.018	0.46	130	Short	Straight / shapeable
M001492761	0.018	0.46	130	Short	J
M001492771	0.018	0.46	130	Long	Straight / shapeable
M001492781	0.018	0.46	130	Long	J
M001492791	0.018	0.46	190	Short	Straight / shapeable
M001492801	0.018	0.46	190	Short	J
M001492811	0.018	0.46	190	Long	Straight / shapeable
M001492821	0.018	0.46	190	Long	J
M001492831	0.018	0.46	300	Short	Straight / shapeable
M001492841	0.018	0.46	300	Short	J
M001492851	0.018	0.46	300	Long	Straight / shapeable
M001492861	0.018	0.46	300	Long	J
H74922150-010	0.014	0.36	145	Trooper Patriot™ Extension Wire*	

*Compatible with 130cm and 190cm Thruway™ Guidewires

Platinum™ Plus

Guidewire LT / Führungsdraht LT / Guide LT / Guía LT

- For negotiation of tortuous anatomy and contralateral approaches. Distal platinum coils for good visualization
- Für gewundene Anatomien und kontralaterale Einsätze. Distale Platinspitze für gute Sichtbarkeit
- Pour négocier les anatomies sinueuses et pour les approches contralatérales. Coils distaux en platine pour une bonne visualisation
- Para el tratamiento de anatomías tortuosas y para procedimientos contralaterales. Coils distales de platino para una buena visualización



Long Tip

Order Number	Diameter (in) (mm)		Length (cm)	Radiopaque Tip Length (cm)	Floppy Tip Length (cm)
M001466000	0.018	0.46	145	15	10
M001466010	0.018	0.46	180	15	10
M001466020	0.018	0.46	260	15	10

Packaged 1/box

Short Tip

Order Number	Diameter (in) (mm)		Length (cm)	Radiopaque Tip Length (cm)	Floppy Tip Length (cm)
M001466040	0.018	0.46	180	7	3
M001466050	0.018	0.46	260	7	3

Packaged 1/box


V-18™ ControlWire™

Peripheral Guidewire / Peripherer Führungsdraht / Guide périphérique / Guía periférica

- Scitanium stainless steel alloy core for torque control and manoeuvrability; ICE™ hydrophilic coating; soft 2 cm shapeable tip
- Drehmomentkontrolle und Lenkbarkeit durch Kerndraht aus Scitanium-Edelstahllegierung; hydrophile ICE™-Beschichtung; weiche, formbare Spitze (2 cm lang)
- Noyau en alliage d'acier inoxydable Scitanium pour un bon contrôle du torque et une bonne maniabilité ; revêtement hydrophile ICE™ ; extrémité souple malléable de 2 cm
- Núcleo con aleación de acero inoxidable Scitanium para un control de torsión y maniobrabilidad; recubrimiento hidrofílico ICE™; punta flexible de 2 cm

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Flexible Tip Length (cm)
M001468480	0.018	0.46	110	8
M001468500	0.018	0.46	150	8
M001468560	0.018	0.46	150	12
M001468520	0.018	0.46	200	8
M001468580	0.018	0.46	200	12
M001468540	0.018	0.46	300	8
M001468600	0.018	0.46	300	12

Packaged 1/box



Victory™

Peripheral Guidewire / Peripherer Führungsdraht / Guide périphérique / Guía periférica

- Victory™ is designed to break through and cross resistant lesions, offering increasing tip loads from 12g to 30g for control and a range of pushability
- Victory™ ermöglicht leichtes Passieren resistenter Läsionen und bietet zudem mit 4 unterschiedlichen Spitzengewichten von 12g bis zu 30g stets höchste Kontrolle und präzises Vorschieben
- Le Guide Périphérique Victory™ est conçu pour franchir et traverser des lésions résistantes ; il offre un bon contrôle et garantit de multiples niveaux de pression
- Las guías Victory™ están diseñadas para avanzar y cruzar lesiones resistentes mediante puntas que van de los 12g a los 30g para proporcionar el control y la capacidad de empuje necesarios

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Radiopaque Tip Length (cm)	Tip Load (g)
H74939232195120	0.018	0.46	195	2.5	12
H74939232300120	0.018	0.46	300	2.5	12
H74939232195180	0.018	0.46	195	2	18
H74939232300180	0.018	0.46	300	2	18
H74939233195250	0.018	0.46	195	2	25
H74939233300250	0.018	0.46	300	2	25
H74939233195300	0.018	0.46	195	2	30
H74939233300300	0.018	0.46	300	2	30

ZIPWire™

Hydrophilic Guidewire - Standard / Hydrophiler Führungsdraht - Standard / Guide hydrophile - standard / Guía hidrofílica estándar

- The ZIPWire™ Hydrophilic Guidewire has focused design improvements to enhance steerability over the leading guide wire technology for better control. Now you have a choice when selecting your guide wire, without compromising performance*
- Der hydrophile Führungsdraht ZIPWire™ wurde dahingehend verbessert, dass die Steuerbarkeit die führende Steuerungstechnologie von Führungsdrähten übertrifft. Bei der Bestimmung des Führungsdrachts haben Sie nun die Auswahl, ohne Leistungseinschränkungen*
- Le guide hydrophile ZIPWire™ bénéficie d'une conception de pointe axée sur l'amélioration de la manœuvrabilité, ce qui permet un contrôle précis. Désormais, vous pouvez choisir votre guide, sans pour autant compromettre la performance*
- La guía hidrofílica ZIPWire™ presenta mejoras de diseño que aportan una mejor navegabilidad frente a la tecnología líder en guías para proporcionar un mayor control. Ahora dispone de un amplio abanico de guías, sin tener que renunciar a un buen rendimiento*

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Core Taper (cm)	Tip Shape
M00146231B1	0.018	0.46	150	3	Angled
M00146232B1	0.018	0.46	180	3	Angled
M00146235B1	0.018	0.46	150	3	Straight
M00146236B1	0.018	0.46	180	3	Straight

Packaged 5/box



*Testing completed by Boston Scientific Corporation. Data on file. Bench test results may not necessarily be indicative of clinical performance (n=15).

Banc d'essai effectué par Boston Scientific Corporation. Données sur fichier. Les résultats du banc d'essai ne reflètent pas nécessairement les résultats cliniques. (n=15).

Tests durchgeführt von Boston Scientific Corporation. Daten im Archiv. Testergebnisse stimmen ggf. nicht mit klinischen Ergebnissen überein. (n=15).

Pruebas efectuadas por Boston Scientific Corporation. Datos disponibles. Los resultados de las pruebas comparativas no son necesariamente indicativos de los resultados clínicos. (n=15).

1 Guidewires

1.3 0.035" (0.89 mm) +/- Guidewires

PTFE coated - Starter™ Guidewire

Platinum™ Plus Guidewire

ZIPWire™ Hydrophilic Guidewire




1.3 0.035" (0.89 mm) +/- Guidewires

PTFE coated - Starter™

Guidewire - FIXED CORE / Führungsdraht - FESTER KERNDRAHT / Guide - AME FIXE / Guía de núcleo FIJO

- For selective catheterization and all diagnostic and interventional procedures
- Zur selektiven Katheterisierung sowie zur allgemeinen Diagnose und für chirurgische Eingriffe und Prozeduren
- Pour la cathétérisation sélective et tous types de procédures diagnostiques et d'interventions
- Para cateterización selectiva y todos los procedimientos de diagnóstico e intervención



Straight Flexible Tip (3 cm Flexible Tip)

Order Number	Diameter (in)	Diameter (mm)	Length (cm)
M001492001	0.018	0.46	150
M001492011	0.025	0.64	40
M001492021	0.025	0.64	100
M001492031	0.025	0.64	125
M001491001	0.025	0.64	150
M001492051	0.028	0.71	40
M001492061	0.028	0.71	150
M001492071	0.032	0.81	40
M001492081	0.032	0.81	150
M001492091	0.035	0.89	40
M001492311	0.035	0.89	80
M001492101	0.035	0.89	100
M001492111	0.035	0.89	125
M001491011	0.035	0.89	150
M001492121	0.035	0.89	180
M001491031	0.035	0.89	260
M001492141	0.038	0.97	100
M001491021	0.038	0.97	150
M001492131	0.038	0.97	260

Packaged 10/box for lengths up to 150 cm included
Packaged 5/box for lengths superior to 150 cm

J-Curved Flexible Tip (3 cm Flexible Tip)


Order Number	Diameter (in)	Diameter (mm)	Length (cm)	J-Curved Tip (mm)
M001491161	0.025	0.64	150	3
M001492211	0.032	0.81	150	3
M001491151	0.035	0.89	150	1.5
M001492221	0.035	0.89	40	3
M001492301	0.035	0.89	80	3
M001492231	0.035	0.89	100	3
M001491171	0.035	0.89	125	3
M001491181	0.035	0.89	150	3
M001491201	0.035	0.89	180	3
M001491211	0.035	0.89	260	3
M001492251	0.035	0.89	150	6
M001492261	0.035	0.89	180	6
M001492271	0.035	0.89	150	12
M001491191	0.038	0.97	150	3
M001492281	0.038	0.97	150	12

Packaged 10/box for lengths up to 150 cm included
Packaged 5/box for lengths superior to 150 cm

PTFE coated - Starter™

Guidewire - MOVEABLE CORE / Führungsdraht - BEWEGLICHER KERNDRAHT / Guide - AME MOBILE / Guía de núcleo MÓVIL

- Offering varying tip flexibility for tortuous vessels and all diagnostic and interventional procedures
- Flexible Spitze für gewundene Gefäße sowie alle diagnostischen und chirurgischen Eingriffe
- Flexibilité variable de l'extrémité pour les vaisseaux d'anatomie sinueuse et pour tout type de procédures diagnostiques et d'interventions
- Ofrece una variedad de puntas flexibles para vasos tortuosos en todos los procedimientos de diagnóstico e intervención



Straight Tip

Order Number	Diameter (in)	Diameter (mm)	Length (cm)
M001491101	0.035	0.89	150
M001491111	0.038	0.97	150

Packaged 10/box

J-Curved

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	J-Curved Tip (mm)
M001491301	0.035	0.89	150	1.5
M001491321	0.035	0.89	150	3
M001491331	0.038	0.97	150	3
M001491351	0.035	0.89	150	6
M001491341	0.035	0.89	260	3

Packaged 10/box for lengths up to 150 cm included
Packaged 5/box for lengths superior to 150 cm

1.3 0.035" (0.89 mm) +/- Guidewires

PTFE coated - Starter™

Guidewire - DOUBLE ENDED / Führungsdraht - MIT 2 ENDEN / Guide - À EXTRÉMITÉ DOUBLE / Guía DE DOBLE PUNTA


Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Straight Tip (cm)	J-Curved Tip (mm)
M001492201	0.035	0.89	150	3	1.5
M001492241	0.035	0.89	150	3	3

Packaged 10/box

PTFE coated - Starter™

Guidewire - NEWTON / Führungsdraht - NEWTON / Guide - NEWTON / Guía DE NEWTON

- For access in tortuous and highly diseased vessels, and smooth advancement during peripheral interventions*
- Für den Zugang in stark geschädigte und gewundene Gefäße und zum reibungslosen Vorschub bei peripheren Eingriffen*
- Pour un accès aux vaisseaux sinueux et sténosés, et une progression aisée au cours d'interventions périphériques*
- Para un acceso en vasos tortuosos y altamente dañados, y un suave avance durante las intervenciones periférico*




Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Flexible Tip (mm)
M001491501	0.035	0.89	150	10
M001491511	0.035	0.89	150	15

Packaged 10/box

PTFE coated - Starter™

Guidewire - NEW BENTSON™ / Führungsdraht - NEW BENTSON™ / Guide - NEW BENTSON™ / Guía DE NEW BENTSON™

- For smooth advancement during peripheral intervention and for all extremely delicate diagnostic and interventional procedures*
- Zum sanften Vorschub bei peripheren Eingriffen sowie für alle extrem empfindlichen Diagnoseprozeduren und chirurgischen Eingriffe*
- Pour une progression aisée au cours d'interventions périphériques et pour toutes les procédures diagnostiques et toutes les interventions délicates*
- Para un avance suave durante la intervención periférico y para todos los procedimientos de diagnóstico e intervención extremadamente delicados*



Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Flexible Tip (mm)
M001491451	0.035	0.89	150	15
M001491471	0.035	0.89	180	15

Packaged 10/box for lengths up to 150 cm included
Packaged 5/box for lengths superior to 150 cm

* Data on file at Boston Scientific.
* Daten in den Firmenunterlagen verfügbar.
* Données sur fichier Boston Scientific.
* Datos disponibles en los archivos de Boston Scientific.

1.3 0.035" (0.89 mm) +/- Guidewires

PTFE coated - Starter™

Guidewire - HEAVY DUTY FIXED CORE / Führungsdraht - FESTER KERNDRAHT FÜR HOHE BEANSPRUCHUNG
Guide - À ÂME FIXE ET RENFORCÉ / Guía DE NÚCLEO FIJO MUY RESISTENTE

- PTFE Coated - Tip Length 3 cm
- PTFE-beschichtet - Länge der Spitze 3 cm
- Revêtement PTFE - Longueur de l'extrémité 3 cm
- Recubrimiento PTFE - Longitud de punta 3 cm

Order Number	Diameter (in)	Diameter (mm)	Length (cm)
M001492501	0.035	0.89	150

Packaged 10/box


PTFE coated - Starter™

Guidewire - ROSEN HEAVY DUTY / Führungsdraht - ROSEN FÜR HOHE BEANSPRUCHUNG / Guide - ROSEN RENFORCÉ
Guía ROSEN HEAVY DUTY

- For additional support during procedures in highly diseased vessels, renal procedures and interventions in anatomy with short distal purchase
- Pour un support supplémentaire au cours des procédures dans des vaisseaux très sténosés, des procédures rénales et d'interventions dans les anatomies à accroche distale courte
- Zur zusätzlichen Unterstützung bei Eingriffen an stark geschädigten Gefäßen, Nierengefäßen sowie chirurgischen Eingriffen bei Anatomien mit kurzem distalem Zugang
- Para obtener soporte adicional durante los procedimientos en vasos altamente dañados, procedimientos renales e intervenciones en anatomía con base distal corta

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Flexible Tip (mm)
M001491551	0.035	0.89	150	1.5


Packaged 10/box



Platinum™ Plus

Guidewire ST / Führungsdraht ST / Guide ST / Guía ST

- For access in anatomy with short distal purchase. Distal platinum coils for good visualization
- Für den Zugang zu Anatomien mit kurzem, distalen Ende. Distale Platinspitze für gute Sichtbarkeit
- Pour accéder aux anatomies à courtes prises distales. Coils distaux en platine pour une bonne visualisation
- Para intervenciones en anatomías con base distal corta. Coils distales de platino para una buena visualización



Short Tip

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Radiopaque Tip Length (cm)	Floppy Tip Length (cm)
M001466100	0.025	0.64	180	9	3
M001466110	0.025	0.64	260	9	3


Packaged 1/box

1.3 0.035" (0.89 mm) +/- Guidewires

ZIPWire™

Hydrophilic Guidewire / Hydrophiler Führungsdraht / Guide hydrophile / Guía hidrofílica

- The ZIPWire™ Hydrophilic Guidewire has focused design improvements to enhance steerability over the leading guide wire technology for better control. Now you have a choice when selecting your guide wire, without compromising performance*
- Le guide hydrophile ZIPWire™ bénéficie d'une conception de pointe axée sur l'amélioration de la manœuvrabilité, ce qui permet un contrôle précis. Désormais, vous pouvez choisir votre guide, sans pour autant compromettre la performance*
- Der hydrophile Führungsdraht ZIPWire™ wurde dahingehend verbessert, dass die Steuerbarkeit die führende Steuerungstechnologie von Führungsdrähten übertrifft. Bei der Bestimmung des Führungsdrahts haben Sie nun die Auswahl, ohne Leistungseinschränkungen*
- La guía hidrofílica ZIPWire™ presenta mejoras de diseño que aportan una mejor navegabilidad frente a la tecnología líder en guías para proporcionar un mayor control. Ahora dispone de un amplio abanico de guías, sin tener que renunciar a un buen rendimiento*



Standard

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Core Taper (cm)	Tip Shape
M00146151B1	0.035	0.89	150	3	Angled
M00146152B1	0.035	0.89	180	3	Angled
M00146154B1	0.035	0.89	260	3	Angled
M00146155B1	0.035	0.89	150	3	Straight
M00146156B1	0.035	0.89	180	3	Straight
M00146159B1	0.035	0.89	260	3	Straight
M00146161B1	0.038	0.97	150	3	Angled
M00146162B1	0.038	0.97	180	3	Angled
M00146164B1	0.038	0.97	260	3	Angled
M00146165B1	0.038	0.97	150	3	Straight
M00146166B1	0.038	0.97	180	3	Straight
M00146169B1	0.038	0.97	260	3	Straight

Packaged 5/box

Stiff


Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Core Taper (cm)	Tip Shape
M00146300B1	0.035	0.89	80	3	Angled
M00146305B1	0.035	0.89	150	3	Angled
M00146306B1	0.035	0.89	150	3	Straight
M00146308B1	0.035	0.89	180	3	Angled
M00146309B1	0.035	0.89	180	3	Straight
M00146316B1	0.035	0.89	260	3	Angled
M00146317B1	0.035	0.89	260	3	Straight
M00146301B1	0.038	0.97	80	3	Angled
M00146311B1	0.038	0.97	150	3	Angled
M00146312B1	0.038	0.97	150	3	Straight
M00146314B1	0.038	0.97	180	3	Angled
M00146315B1	0.038	0.97	180	3	Straight

Packaged 5/box

Standard Long Taper

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Core Taper (cm)	Tip Shape
M00146325B1	0.035	0.89	150	8	Angled
M00146331B1	0.038	0.97	150	8	Angled


Packaged 5/box



Stiff Long Taper

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Core Taper (cm)	Tip Shape
M00146318B1	0.035	0.89	150	5	Angled
M00146319B1	0.035	0.89	260	5	Angled
M00146320B1	0.035	0.89	150	5	Straight
M00146321B1	0.038	0.97	150	5	Angled
M00146322B1	0.038	0.97	150	5	Straight

Packaged 5/box



J-Tip

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Core Taper (cm)	J-Curved Tip (mm)
M00146350B1	0.035	0.89	150	3	3
M00146355B1	0.038	0.97	150	3	3

Packaged 5/box

*Testing completed by Boston Scientific Corporation. Data on file. Bench test results may not necessarily be indicative of clinical performance (n=15).
 Banc d'essai effectué par Boston Scientific Corporation. Données sur fichier. Les résultats du banc d'essai ne reflètent pas nécessairement les résultats cliniques. (n=15).
 Tests durchgeführt von Boston Scientific Corporation. Daten im Archiv. Testergebnisse stimmen ggf. nicht mit klinischen Ergebnissen überein. (n=15).
 Pruebas efectuadas por Boston Scientific Corporation. Datos disponibles. Los resultados de las pruebas comparativas no son necesariamente indicativos de los resultados clínicos. (n=15).

1 Guidewires

1.4 Speciality Guidewires

Magic Torque™ *Glidex™ Coated Hydrophilic Guidewire*

Amplatz Super Stiff™ *Guidewire*


Back-up Meier™ *Steerable Guidewire*



Magic Torque™

Glidex™ Coated Hydrophilic Guidewire / Glidex™ beschichteter hydrophiler Führungsdraht / Guide hydrophile à revêtement Glidex™
Guía hidrofílica con recubrimiento Glidex™

- 4 distal platinum Magic™ markers spaced at 1 cm increment for good visualization
- 4 distale Platin-Magic™ Marker im Abstand von 1 cm zur besseren Visualisierung
- 4 marqueurs distaux en platine Magic™ espacés de 1 cm améliorent la visualisation
- 4 marcas de platino distales Magic™ espaciadas un centímetro para una visualización mejorada



Tip Length 5 cm


Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Tip Length (cm)
M001465911	0.035	0.89	180	Shapeable
M001465921	0.035	0.89	260	Shapeable

Packaged 3/box

Amplatz Super Stiff™

Guidewire / Führungsdraht / Guide / Guía

- Combining great stiffness, strength and stability with safety of a soft, atraumatic floppy tip
- Kombination aus höherer Steifigkeit, Widerstandskraft und Stabilität sowie der Sicherheit durch eine weiche, atraumatische Spitze
- Associe un niveau élevé de rigidité, de résistance et de stabilité à la sécurité d'une extrémité souple atraumatique
- Combinación de gran rigidez, resistencia y estabilidad con la seguridad de una punta flexible y atraumática



Straight (6 cm Flexible Tip)

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Tip Length (cm)
M001465631	0.035	0.89	75	7
M001465231	0.035	0.89	145	7
M001465251	0.035	0.89	180	7
M001465261	0.035	0.89	260	7
M001465641	0.038	0.97	75	6
M001465171	0.038	0.97	145	6
M001465191	0.038	0.97	180	6
M001465201	0.038	0.97	260	6

Packaged 5/box

Straight (ST Short Tip)

Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Tip Length (cm)
M001465091	0.035	0.89	260	1
M001465241	0.035	0.89	145	4
M001465101	0.038	0.97	260	1
M001465181	0.038	0.97	145	4

Packaged 5/box

3 mm J-tip


Order Number	Diameter (in)	Diameter (mm)	Length (cm)	Tip Length (cm)
M001465001	0.035	0.89	145	6
M001465011	0.035	0.89	180	6
M001465021	0.035	0.89	260	6
M001465031	0.038	0.97	145	6
M001465041	0.038	0.97	180	6

Packaged 5/box

Back-up Meier™

Steerable Guidewire / Lenkbarer Führungsdraht / Guide orientable / Guía orientable

- Designed to balance the support, flexibility and radiopacity required for efficient catheter placement
- Ausgewogene Kombination von Eigenschaften wie Unterstützung, Flexibilität und Röntgensichtbarkeit, die für die effiziente Katheterplatzierung erforderlich sind
- Offre le niveau équilibré de support, de flexibilité et de radio-opacité nécessaire pour une mise en place efficace du cathéter
- Diseñada para equilibrar el soporte, la flexibilidad y radiopacidad necesaria para una colocación eficaz del catéter



Flexible Tip

Order Number	Diameter (in) (mm)		Length (cm)	Tip Configuration	Tip Length (cm)
H965SCH306001	0.035	0.89	185	J	16
H965SCH306011	0.035	0.89	300	C-Shaped	11

Packaged 5/box

2 Balloon Dilatation Catheters

2.1 0.014" (0.36mm) Guidewire Compatible

Coyote™ Balloon Dilatation Catheter

Coyote™ ES Balloon Dilatation Catheter

Ultra-Soft™ SV Monorail™ Balloon Dilatation Catheter



2.1 0.014" (0.36mm) Guidewire Compatible

Coyote™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter balón de dilatación

- GET TO - with four catheter configurations optimized for maximum push and excellent track
GET THROUGH - with an ultra-low lesion entry profile (0.017"/0.43 mm) and best-in-class crossing profile (0.030"/0.76 mm)
GET GOING - with deflation times of less than 10 seconds* and balloon lengths up to 220 mm
- GET TO - Mit vier Katheteroptionen für maximalen Vorschub und ausgezeichnete Steuerbarkeit
GET THROUGH - Mit einem extrem niedrigen Läsionseingangsprofil (0,017"/0,43 mm) und dem besten Crossing-Profil seiner Klasse (0,030"/0,76 mm)
GET GOING - Mit Deflationszeiten von weniger als 10 Sekunden* und Ballonlängen bis zu 220 mm
- ACCÉDER - Quatre configurations optimales pour un maximum de push et une excellente navigabilité.
FRANCHIR - Un profil d'entrée ultra-fin de 0.017" (0.43mm) et un profil de franchissement incomparable (0.76mm)
AGIR - Un temps de déflation inférieur à 10 secondes* et des longueurs de ballon jusqu'à 220mm
- LLEGA - Con cuatro configuraciones de catéter optimizadas para lograr el máximo empuje y una excelente navegabilidad
CRUZA - Con un perfil de entrada en lesión ultrabajo (0,017 pulg./0,43 mm) y el mejor perfil de cruce de su clase (0,030 pulg./0,76 mm)
ACTÚA - Con tiempos de desinflado inferiores a 10 segundos* y longitudes de balón de hasta 220mm

Coyote™ Monorail™ Balloon Dilatation Catheter

Order Number	Order Number	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
Shaft Length 90 cm	Shaft Length 150 cm						
H74939185150490	H74939185150410	1.5	40	4	1.33	14	1418
H74939185200690	H74939185200610	2.0	60	4	1.33	14	1418
H74939185200890	H74939185200810	2.0	80	4	1.33	14	1418
H74939185201090	H74939185201010	2.0	100	4	1.33	14	1418
H74939185201290	H74939185201210	2.0	120	4	1.33	14	1418
H74939185201590	H74939185201510	2.0	150	4	1.33	14	1418
H74939185202290	H74939185202210	2.0	220	4	1.33	14	1418
H74939185250690	H74939185250610	2.5	60	4	1.33	14	1418
H74939185250890	H74939185250810	2.5	80	4	1.33	14	1418
H74939185251090	H74939185251010	2.5	100	4	1.33	14	1418
H74939185251290	H74939185251210	2.5	120	4	1.33	14	1418
H74939185251590	H74939185251510	2.5	150	4	1.33	14	1418
H74939185252290	H74939185252210	2.5	220	4	1.33	14	1418
H74939185300690	H74939185300610	3.0	60	4	1.33	14	1418
H74939185300890	H74939185300810	3.0	80	4	1.33	14	1418
H74939185301090	H74939185301010	3.0	100	4	1.33	14	1418
H74939185301290	H74939185301210	3.0	120	4	1.33	14	1418
H74939185301590	H74939185301510	3.0	150	4	1.33	14	1418
H74939185302290	H74939185302210	3.0	220	4	1.33	14	1418
H74939185350690	H74939185350610	3.5	60	4	1.33	14	1418
H74939185350890	H74939185350810	3.5	80	4	1.33	14	1418
H74939185351090	H74939185351010	3.5	100	4	1.33	14	1418
H74939185351290	H74939185351210	3.5	120	4	1.33	14	1418
H74939185351590	H74939185351510	3.5	150	4	1.33	14	1418
H74939185352290	H74939185352210	3.5	220	4	1.33	14	1418
H74939185400690	H74939185400610	4.0	60	4	1.33	14	1418
H74939185400890	H74939185400810	4.0	80	4	1.33	14	1418
H74939185401090	H74939185401010	4.0	100	4	1.33	14	1418
H74939185401290	H74939185401210	4.0	120	4	1.33	14	1418
H74939185401590	H74939185401510	4.0	150	4	1.33	14	1418
H74939185402290	H74939185402210	4.0	220	4	1.33	14	1418

* Average measurements taken by Boston Scientific (n=3, 2 mm x 120 mm balloon). Date on file. Bench testing may not be representative of clinical performance

2.1 0.014" (0.36mm) Guidewire Compatible

Coyote™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter balón de dilatación

Coyote™ Over-the-Wire Balloon Dilatation Catheter

Order Number Catheter Length 90 cm	Order Number Catheter Length 150 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
H74939186150490	H74939186150410	1.5	40	4	1.33	14	1418
H74939186200690	H74939186200610	2.0	60	4	1.33	14	1418
H74939186200890	H74939186200810	2.0	80	4	1.33	14	1418
H74939186201090	H74939186201010	2.0	100	4	1.33	14	1418
H74939186201290	H74939186201210	2.0	120	4	1.33	14	1418
H74939186201590	H74939186201510	2.0	150	4	1.33	14	1418
H74939186202290	H74939186202210	2.0	220	4	1.33	14	1418
H74939186250690	H74939186250610	2.5	60	4	1.33	14	1418
H74939186250890	H74939186250810	2.5	80	4	1.33	14	1418
H74939186251090	H74939186251010	2.5	100	4	1.33	14	1418
H74939186251290	H74939186251210	2.5	120	4	1.33	14	1418
H74939186251590	H74939186251510	2.5	150	4	1.33	14	1418
H74939186252290	H74939186252210	2.5	220	4	1.33	14	1418
H74939186300690	H74939186300610	3.0	60	4	1.33	14	1418
H74939186300890	H74939186300810	3.0	80	4	1.33	14	1418
H74939186301090	H74939186301010	3.0	100	4	1.33	14	1418
H74939186301290	H74939186301210	3.0	120	4	1.33	14	1418
H74939186301590	H74939186301510	3.0	150	4	1.33	14	1418
H74939186302290	H74939186302210	3.0	220	4	1.33	14	1418
H74939186350690	H74939186350610	3.5	60	4	1.33	14	1418
H74939186350890	H74939186350810	3.5	80	4	1.33	14	1418
H74939186351090	H74939186351010	3.5	100	4	1.33	14	1418
H74939186351290	H74939186351210	3.5	120	4	1.33	14	1418
H74939186351590	H74939186351510	3.5	150	4	1.33	14	1418
H74939186352290	H74939186352210	3.5	220	4	1.33	14	1418
H74939186400690	H74939186400610	4.0	60	4	1.33	14	1418
H74939186400890	H74939186400810	4.0	80	4	1.33	14	1418
H74939186401090	H74939186401010	4.0	100	4	1.33	14	1418
H74939186401290	H74939186401210	4.0	120	4	1.33	14	1418
H74939186401590	H74939186401510	4.0	150	4	1.33	14	1418
H74939186402290	H74939186402210	4.0	220	4	1.33	14	1418



2.1 0.014" (0.36mm) Guidewire Compatible

Coyote™ ES

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter balón de dilatación

- This ultra-low profile balloon is available in both Over-The-Wire (OTW) and Monorail™ platforms and is designed for exceptional crossing
- Dieser extrem flache Ballonkatheter ist sowohl für Over the Wire (OTW) als auch Monorail™ Systeme verfügbar und für hervorragende Passageeigenschaften konzipiert
- Ce ballonnet à profil ultra bas existe en version Coaxiale (OTW) et Monorail™ et est conçu pour assurer un franchissement exceptionnel
- Esta balón de perfil ultrabajo está disponible tanto para la plataforma Over-the-Wire (OTW) como para la plataforma Monorail™ y se ha diseñado con una capacidad de cruce excepcional

Coyote™ ES Over-the-Wire Balloon Dilatation Catheter

Catheter Length 142 cm	Order Number Catheter Length 143 cm	Catheter Length 144 cm	Catheter Length 145 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
H74939134152010				1.5	20	4	1.33	14	1418
H74939134202010				2.0	20	4	1.33	14	1418
	H74939134203010			2.0	30	4	1.33	14	1418
		H74939134204010		2.0	40	4	1.33	14	1418
	H74939134252010			2.5	20	4	1.33	14	1418
		H74939134253010		2.5	30	4	1.33	14	1418
			H74939134254010	2.5	40	4	1.33	14	1418
	H74939134302010			3.0	20	4	1.33	14	1418
		H74939134303010		3.0	30	4	1.33	14	1418
			H74939134304010	3.0	40	4	1.33	14	1418
	H74939134352010			3.5	20	4	1.33	12	1216
		H74939134353010		3.5	30	4	1.33	12	1216
			H74939134354010	3.5	40	4	1.33	12	1216
	H74939134402010			4.0	20	4	1.33	12	1216
		H74939134403010		4.0	30	4	1.33	12	1216
			H74939134404010	4.0	40	4	1.33	12	1216

Coyote™ ES Monorail™ Balloon Dilatation Catheter

Shaft Length 143 cm	Order Number Shaft Length 144 cm	Shaft Length 145 cm	Shaft Length 146 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
H74939135152010				1.5	20	4	1.33	14	1418
H74939135202010				2.0	20	4	1.33	14	1418
	H74939135203010			2.0	30	4	1.33	14	1418
		H74939135204010		2.0	40	4	1.33	14	1418
	H74939135252010			2.5	20	4	1.33	14	1418
		H74939135253010		2.5	30	4	1.33	14	1418
			H74939135254010	2.5	40	4	1.33	14	1418
	H74939135302010			3.0	20	4	1.33	14	1418
		H74939135303010		3.0	30	4	1.33	14	1418
			H74939135304010	3.0	40	4	1.33	14	1418
	H74939135352010			3.5	20	4	1.33	12	1216
		H74939135353010		3.5	30	4	1.33	12	1216
			H74939135354010	3.5	40	4	1.33	12	1216
	H74939135402010			4.0	20	4	1.33	12	1216
		H74939135403010		4.0	30	4	1.33	12	1216
			H74939135404010	4.0	40	4	1.33	12	1216

2.1 0.014" (0.36mm) Guidewire Compatible

Ultra-Soft™ SV

Monorail™ Balloon Dilatation Catheter / Monorail™ Ballondilatationskatheter / Cathéter de dilatation à ballonnet Monorail™ / Catéter balón de dilatación Monorail™

- The Ultra-Soft™ SV is a combination of proximal stiffness for pushability and distal shaft flexibility for improved tracking
- L'Ultra-Soft™ SV est à la fois rigide en proximal pour donner du push avec un shaft souple en distalité pour une excellente navigabilité
- Ultra-Soft™ SV USSV ist eine ausgezeichnete Kombination von proximaler Festigkeit für verbesserten Vorschub und höherer Flexibilität des distalen Schafts für optimierte Steuerbarkeit
- El Ultra-Soft™ SV es una combinación de rigidez en su parte proximal para empuje y flexibilidad en su vaina distal para mejorar la Navegabilidad

Order Number Shaft Length 140 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
M0013894015200	1.5	20	4	1.33	12	1216
M0013894020200	2	20	4	1.33	12	1216
M0013894020300	2	30	4	1.33	12	1216
M0013894025200	2.5	20	4	1.33	14	1419
M0013894025300	2.5	30	4	1.33	14	1419
M0013894030200	3	20	4	1.33	14	1419
M0013894030300	3	30	4	1.33	14	1419
M0013894035200	3.5	20	4	1.33	12	1216
M0013894035300	3.5	30	4	1.33	12	1216
M0013894040200	4	20	4	1.33	12	1216
M0013894040300	4	30	4	1.33	12	1216

2 Balloon Dilatation Catheters

2.2 0.018" (0.46mm) Guidewire Compatible

Sterling™
Sterling SL™
Symmetry™
Ultra-Soft™ SV



2.2 0.018" (0.46mm) Guidewire Compatible

Sterling™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter con balón de dilatación

- The Sterling™ Balloon Dilatation Catheter meets the challenges of lower extremity, carotid* and renal anatomy with a powerful combination of deliverability and trackability** resulting in clinical versatility. Excellent trackability and a broad size offering provide multiple solutions for all of your peripheral intervention needs
- Der Sterling™ Ballondilatationskatheter erfüllt auch die Anforderungen beim Einsatz in den unteren Extremitäten, in der Karotis sowie der Nierenarterie und ist als leistungsfähige Kombination aus Platzierbarkeit und Steuerbarkeit klinisch vielseitig einsetzbar. Durch seine hervorragende Steuerbarkeit und eine große Auswahl an Größen bietet er zahlreiche Lösungen für alle Anforderungen peripherer Eingriffe
- Le cathéter de dilatation à ballonnet Sterling™ répond aux difficultés anatomiques des membres inférieurs, de la carotide et de la rénale grâce à ses qualités de mise en place et de navigabilité, répondant ainsi aux besoins cliniques. L'excellente navigabilité et l'étendue de la gamme offrent des solutions multiples répondant à l'ensemble de vos procédures périphériques interventionnelles
- El balón de dilatación Sterling™ se adecua a los retos que suponen las extremidades inferiores, la anatomía carotídea y la renal con una potente combinación de liberación y navegabilidad, que da como resultado versatilidad clínica. La excelente navegabilidad y la amplitud de la oferta proporciona soluciones múltiples para todas las necesidades de las intervenciones periféricas.

* Sterling Monorail indication only; Sterling OTW is not indicated for use in carotid arteries.

**Data on File Boston Scientific 2010

Over-the-Wire - Shaft Length 40 cm / 80 cm / 135 cm

Shaft Length 40 cm	Order Number Shaft Length 80 cm	Shaft Length 135 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
H74939032402040	H74939032402080	H74939032402010	4	20	4	1.33	14	1418
—	H74939032403080	H74939032403010	4	30	4	1.33	14	1418
H74939032404040	H74939032404080	H74939032404010	4	40	4	1.33	14	1418
—	H74939032406080	H74939032406010	4	60	4	1.33	14	1418
—	H74939032408080	H74939032408010	4	80	4	1.33	14	1418
—	H74939032401080	H74939032401010	4	100	4	1.33	14	1418
H74939032502040	H74939032502080	H74939032502010	5	20	4	1.33	14	1418
—	H74939032503080	H74939032503010	5	30	4	1.33	14	1418
H74939032504040	H74939032504080	H74939032504010	5	40	4	1.33	14	1418
—	H74939032506080	H74939032506010	5	60	4	1.33	14	1418
—	H74939032508080	H74939032508010	5	80	4	1.33	14	1418
—	H74939032501080	H74939032501010	5	100	4	1.33	14	1418
H74939032602040	H74939032602080	H74939032602010	6	20	4	1.33	14	1418
—	H74939032603080	H74939032603010	6	30	4	1.33	14	1418
H74939032604040	H74939032604080	H74939032604010	6	40	4	1.33	14	1418
—	H74939032606080	H74939032606010	6	60	4	1.33	14	1418
—	H74939032608080	H74939032608010	6	80	5	1.7	14	1418
—	H74939032601080	H74939032601010	6	100	5	1.7	14	1418
H74939032702040	H74939032702080	H74939032702010	7	20	4	1.33	14	1418
—	H74939032703080	H74939032703010	7	30	4	1.33	14	1418
H74939032704040	H74939032704080	H74939032704010	7	40	4	1.33	14	1418
—	H74939032706080	H74939032706010	7	60	5	1.7	14	1418
—	H74939032708080	H74939032708010	7	80	6	2.0	14	1418
—	H74939032701080	H74939032701010	7	100	6	2.0	14	1418
H74939032802040	H74939032802080	H74939032802010	8	20	4	1.33	12	1216
—	H74939032803080	H74939032803010	8	30	4	1.33	12	1216
H74939032804040	H74939032804080	H74939032804010	8	40	4	1.33	12	1216
—	H74939032806080	H74939032806010	8	60	5	1.7	12	1216
—	H74939032808080	H74939032808010	8	80	6	2.0	12	1216
—	H74939032902080	H74939032902010	9	20	5	1.7	10	1013
—	H74939032903080	H74939032903010	9	30	5	1.7	10	1013
—	H74939032904080	H74939032904010	9	40	5	1.7	10	1013
—	H74939032906080	H74939032906010	9	60	5	1.7	10	1013
—	H74939032908080	H74939032908010	9	80	6	2.0	10	1013
—	H74939032102080	H74939032102010	10	20	5	1.7	10	1013
—	H74939032103080	H74939032103010	10	30	5	1.7	10	1013
—	H74939032104080	H74939032104010	10	40	5	1.7	10	1013
—	H74939032106080	H74939032106010	10	60	6	2.0	10	1013
—	H74939032108080	H74939032108010	10	80	6	2.0	10	1013

Specifications for all Sterling™ balloon dilatation catheters
Recommended guide wire: 0.014" / 0.018" (0.36 mm / 0.46 mm)
Recommended guide catheter size: 2F (0.67 mm) sizes higher than sheath compatibility
Nominal pressure: 6 ATM (608 kPa) across the range

Boston
Scientific



2.2 0.018" (0.46mm) Guidewire Compatible

Sterling™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter con balón de dilatación

Monorail™ - Shaft Length 80 cm / 135 cm

Order Number	Shaft Length 80 cm	Shaft Length 135 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
—	—	H74939031302010	3	20	4	1.33	14	1418
—	—	H74939031303010	3	30	4	1.33	14	1418
—	—	H74939031304010	3	40	4	1.33	14	1418
—	—	H74939031306010	3	60	4	1.33	14	1418
—	—	H74939031352010	3.5	20	4	1.33	14	1418
—	—	H74939031353010	3.5	30	4	1.33	14	1418
—	—	H74939031354010	3.5	40	4	1.33	14	1418
H74939031401080	—	—	4	10	4	1.33	14	1418
H74939031401580	—	—	4	15	4	1.33	14	1418
H74939031402080	—	H74939031402010	4	20	4	1.33	14	1418
H74939031403080	—	H74939031403010	4	30	4	1.33	14	1418
H74939031404080	—	H74939031404010	4	40	4	1.33	14	1418
—	—	H74939031406010	4	60	4	1.33	14	1418
H74939031451580	—	—	4.5	15	4	1.33	14	1418
H74939031452080	—	H74939031452010	4.5	20	4	1.33	14	1418
H74939031453080	—	H74939031453010	4.5	30	4	1.33	14	1418
H74939031454080	—	H74939031454010	4.5	40	4	1.33	14	1418
H74939031501080	—	—	5	10	4	1.33	14	1418
H74939031501580	—	—	5	15	4	1.33	14	1418
H74939031502080	—	H74939031502010	5	20	4	1.33	14	1418
H74939031503080	—	H74939031503010	5	30	4	1.33	14	1418
H74939031504080	—	H74939031504010	5	40	4	1.33	14	1418
—	—	H74939031506010	5	60	4	1.33	14	1418
H74939031551580	—	—	5.5	15	4	1.33	14	1418
H74939031552080	—	H74939031552010	5.5	20	4	1.33	14	1418
H74939031553080	—	H74939031553010	5.5	30	4	1.33	14	1418
H74939031554080	—	H74939031554010	5.5	40	4	1.33	14	1418
H74939031601080	—	—	6	10	4	1.33	14	1418
H74939031601580	—	—	6	15	4	1.33	14	1418
H74939031602080	—	H74939031602010	6	20	4	1.33	14	1418
H74939031603080	—	H74939031603010	6	30	4	1.33	14	1418
H74939031604080	—	H74939031604010	6	40	4	1.33	14	1418
—	—	H74939031606010	6	60	4	1.33	14	1418
H74939031651580	—	—	6.5	15	4	1.33	14	1418
H74939031652080	—	H74939031652010	6.5	20	4	1.33	14	1418
H74939031653080	—	H74939031653010	6.5	30	4	1.33	14	1418
H74939031654080	—	H74939031654010	6.5	40	4	1.33	14	1418
H74939031701080	—	—	7	10	4	1.33	14	1418
H74939031701580	—	—	7	15	4	1.33	14	1418
H74939031702080	—	H74939031702010	7	20	4	1.33	14	1418
H74939031703080	—	H74939031703010	7	30	4	1.33	14	1418
H74939031704080	—	H74939031704010	7	40	4	1.33	14	1418
—	—	H74939031706010	7	60	5	1.7	14	1418
H74939031801080	—	—	8	10	4	1.33	12	1216
H74939031801580	—	—	8	15	4	1.33	12	1216
H74939031802080	—	H74939031802010	8	20	4	1.33	12	1216
H74939031803080	—	H74939031803010	8	30	4	1.33	12	1216
H74939031804080	—	H74939031804010	8	40	4	1.33	12	1216
—	—	H74939031806010	8	60	5	1.7	12	1216

Specifications for all Sterling™ balloon dilatation catheters
 Recommended guide wire: 0.014" / 0.018" (0.36 mm / 0.46 mm)
 Recommended guide catheter size: 2F (0.67 mm) sizes higher than sheath compatibility
 Lesion Entry Profile: 0.022" (0.56 mm)
 Nominal pressure: 6 ATM (608 kPa) across the range

Ordering Information - 27
Peripheral Interventions

Boston
Scientific



2.2 0.018" (0.46mm) Guidewire Compatible

Sterling SL™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter con balón de dilatación

Over-the-Wire

Order Number	Shaft Length 90 cm	Shaft Length 150 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
H74939148208090		H74939148208010	2	80	4	1.33	14	1418
H74939148201090		H74939148201010	2	100	4	1.33	14	1418
H74939148201290		H74939148201210	2	120	4	1.33	14	1418
H74939148201590		H74939148201510	2	150	4	1.33	14	1418
H74939148258090		H74939148258010	2.5	80	4	1.33	14	1418
H74939148251090		H74939148251010	2.5	100	4	1.33	14	1418
H74939148251290		H74939148251210	2.5	120	4	1.33	14	1418
H74939148251590		H74939148251510	2.5	150	4	1.33	14	1418
H74939148308090		H74939148308010	3	80	4	1.33	14	1418
H74939148301090		H74939148301010	3	100	4	1.33	14	1418
H74939148301290		H74939148301210	3	120	4	1.33	14	1418
H74939148301590		H74939148301510	3	150	4	1.33	14	1418
H74939148358090		H74939148358010	3.5	80	4	1.33	14	1418
H74939148351090		H74939148351010	3.5	100	4	1.33	14	1418
H74939148351290		H74939148351210	3.5	120	4	1.33	14	1418
H74939148351590		H74939148351510	3.5	150	4	1.33	14	1418
—		H74939148408010	4	80	4	1.33	14	1418
—		H74939148401010	4	100	4	1.33	14	1418
H74939148401290		H74939148401210	4	120	4	1.33	14	1418
H74939148401590		H74939148401510	4	150	4	1.33	14	1418

Monorail

Order Number	Shaft Length 90 cm	Shaft Length 150 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
H74939147208090		H74939147208010	2	80	4	1.33	14	1418
H74939147201090		H74939147201010	2	100	4	1.33	14	1418
H74939147201290		H74939147201210	2	120	4	1.33	14	1418
H74939147201590		H74939147201510	2	150	4	1.33	14	1418
H74939147258090		H74939147258010	2.5	80	4	1.33	14	1418
H74939147251090		H74939147251010	2.5	100	4	1.33	14	1418
H74939147251290		H74939147251210	2.5	120	4	1.33	14	1418
H74939147251590		H74939147251510	2.5	150	4	1.33	14	1418
H74939147308090		H74939147308010	3	80	4	1.33	14	1418
H74939147301090		H74939147301010	3	100	4	1.33	14	1418
H74939147301290		H74939147301210	3	120	4	1.33	14	1418
H74939147301590		H74939147301510	3	150	4	1.33	14	1418
H74939147358090		H74939147358010	3.5	80	4	1.33	14	1418
H74939147351090		H74939147351010	3.5	100	4	1.33	14	1418
H74939147351290		H74939147351210	3.5	120	4	1.33	14	1418
H74939147351590		H74939147351510	3.5	150	4	1.33	14	1418
H74939147408090		H74939147408010	4	80	4	1.33	14	1418
H74939147401090		H74939147401010	4	100	4	1.33	14	1418
H74939147401290		H74939147401210	4	120	4	1.33	14	1418
H74939147401590		H74939147401510	4	150	4	1.33	14	1418

Specifications for all Sterling™ balloon dilatation catheters
 Recommended guide wire: 0.014" / 0.018" (0.36 mm / 0.46 mm)
 Recommended guide catheter size: 2F (0.67 mm) sizes higher than sheath compatibility
 Lesion Entry Profile: 0.022" (0.56 mm)
 Nominal pressure: 6 ATM (608 kPa) across the range

2.2 0.018" (0.46mm) Guidewire Compatible

Symmetry™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter con balón de dilatación

Order Number		Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F) (mm)		Rated Burst Pressure (ATM) (kPa)	
Shaft Length 90 cm	Shaft Length 150 cm						
M001103110	M001103350	1.5	20	4	1.33	15	1520
M001103120	M001103360	1.5	40	4	1.33	15	1520
—	M001103370	1.5	100	4	1.33	15	1520
M001103140	M001103400	2	20	4	1.33	15	1520
M001103150	M001103410	2	40	4	1.33	15	1520
M001103160	M001103420	2	100	4	1.33	15	1520
M001103170	M001103430	2.5	20	4	1.33	15	1520
M001103180	M001103440	2.5	40	4	1.33	15	1520
M001103190	M001103450	2.5	100	4	1.33	15	1520
M001103200	M001103460	3	20	4	1.33	15	1520
M001103210	M001103470	3	40	4	1.33	15	1520
M001103220	M001103480	3	100	4	1.33	15	1520
M001103230	M001103490	3.5	20	4	1.33	15	1520
M001103240	M001103500	3.5	40	4	1.33	15	1520
M001103250	M001103510	3.5	100	4	1.33	15	1520
M001103260	M001103520	4	20	4	1.33	15	1520
M001103270	M001103530	4	40	4	1.33	15	1520
M001103280	M001103540	4	100	5	1.7	15	1520
M001103290	M001103550	4.5	20	5	1.7	15	1520
M001103300	M001103560	4.5	40	5	1.7	15	1520
M001103310	M001103570	5	20	5	1.7	15	1520
M001103320	M001103580	5	40	5	1.7	15	1520
M001103850	M001103590	5.5	20	5	1.7	15	1520
M001103860	M001103600	5.5	40	5	1.7	15	1520
M001103330	M001103610	6	20	5	1.7	15	1520
M001103340	M001103620	6	40	5	1.7	15	1520



2.2 0.018" (0.46mm) Guidewire Compatible

Symmetry™

Balloon Dilatation Catheter - STIFF SHAFT / Ballondilatationskatheter - STEIFER SCHAFT

Cathéter de dilatation à ballonnet - SHAFT RIGIDE / Catéter con balón de dilatación DE EXTREMO RÍGIDO

Shaft Length 90 cm	Order Number Shaft Length 135 cm	Shaft Length 150 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
M001105030	M001105270	M001105510	1.5	20	4	1.33	15	1520
M001105040	M001105280	—	1.5	40	4	1.33	15	1520
M001105050	M001105290	M001105530	1.5	100	4	1.33	15	1520
M001105060	M001105300	M001105540	2	20	4	1.33	15	1520
M001105070	M001105310	M001105550	2	40	4	1.33	15	1520
M001105080	M001105320	M001105560	2	100	4	1.33	15	1520
M001105090	M001105330	M001105570	2.5	20	4	1.33	15	1520
M001105100	M001105340	M001105580	2.5	40	4	1.33	15	1520
M001105110	M001105350	M001105590	2.5	100	4	1.33	15	1520
M001105120	M001105360	M001105600	3	20	4	1.33	15	1520
M001105130	M001105370	M001105610	3	40	4	1.33	15	1520
M001105140	M001105380	M001105620	3	100	4	1.33	15	1520
M001105150	M001105390	—	3.5	20	4	1.33	15	1520
M001105160	M001105400	M001105640	3.5	40	4	1.33	15	1520
M001105170	M001105410	M001105650	3.5	100	4	1.33	15	1520
M001105180	M001105420	M001105660	4	20	4	1.33	15	1520
M001105190	M001105430	M001105670	4	40	4	1.33	15	1520
M001105200	M001105440	M001105680	4	100	5	1.7	15	1520
M001105210	M001105450	M001105690	4.5	20	5	1.7	15	1520
M001105220	M001105460	M001105700	4.5	40	5	1.7	15	1520
M001105230	M001105470	M001105710	5	20	5	1.7	15	1520
M001105240	M001105480	M001105720	5	40	5	1.7	15	1520
M001105250	M001105490	M001105730	6	20	5	1.7	15	1520
M001105260	M001105500	M001105740	6	40	5	1.7	15	1520



2.2 0.018" (0.46mm) Guidewire Compatible

Ultra-Soft™ SV

Monorail™ Balloon Dilatation Catheter / Monorail™ Ballondilationskatheter

Cathéter de dilatation à ballonnet Monorail™ / Catéter balón de dilatación Monorail™

Order Number	Order Number	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
Shaft Length 90 cm	Shaft Length 150 cm						
M0013899040150	M0013895040150	4.0	15	4	1.33	14	1419
M0013899040200	M0013895040200	4.0	20	4	1.33	14	1419
M0013899045150	M0013895045150	4.5	15	4	1.33	14	1419
M0013899045200	M0013895045200	4.5	20	4	1.33	14	1419
M0013899050150	M0013895050150	5.0	15	4	1.33	14	1419
M0013899050200	M0013895050200	5.0	20	4	1.33	14	1419
M0013899055150	M0013895055150	5.5	15	5	1.33	14	1419
M0013899055200	M0013895055200	5.5	20	5	1.33	14	1419
M0013899060150	M0013895060150	6.0	15	5	1.33	14	1419
M0013899060200	M0013895060200	6.0	20	5	1.33	14	1419
M0013899065150	M0013895065150	6.5	15	5	1.33	14	1419
M0013899065200	M0013895065200	6.5	20	5	1.33	14	1419
M0013899070150	M0013895070150	7.0	15	5	1.33	14	1419
M0013899070200	M0013895070200	7.0	20	5	1.33	14	1419



2 Balloon Dilatation Catheters

2.3 0.035" (0.89 mm) Guidewire Compatible

Mustang™ Balloon Dilatation Catheter

Wanda™ Balloon Dilatation Catheter

XXL™ Balloon Dilatation Catheter



2.3 0.035" (0.89 mm) Guidewire Compatible

Mustang™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter con balón de dilatación

- With excellent cross and track, powerful dilatation, longer lengths, and smaller sheath sizes, Mustang is the 0.035" (0.89mm) balloon that offers excellence across the board
- Mustang steht für Exzellenz in jeder Hinsicht: der 0,035" (0,89mm) große Ballon bietet mit größeren Längen und kleineren Schleusen hervorragende "cross and track"-Eigenschaften und eine kraftvolle Dilatation
- Grâce à ses capacités de franchissement et de navigabilité supérieures, sa force de dilatation, ses grandes longueurs, et sa compatibilité avec des introducteurs de plus petite taille, Mustang est le seul ballon 0,035" (0,89 mm) offrant l'excellence à tous les niveaux
- Con una navegabilidad y cruce superiores, una gran fuerza de dilatación, unas longitudes mayores y unos tamaños de vaina más reducidos, Mustang es el balón de 0,035" (0,89 mm) que ofrece una gran calidad en todos los niveles

Catheter Length 40 cm	Order Number Catheter Length 75 cm	Catheter Length 135 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
H74939171030240	H74939171030270	H74939171030210	3	20	5	1.7	24	2431
H74939171030340	H74939171030370	H74939171030310	3	30	5	1.7	24	2431
H74939171030440	H74939171030470	H74939171030410	3	40	5	1.7	24	2431
H74939171030640	H74939171030670	H74939171030610	3	60	5	1.7	24	2431
H74939171030840	H74939171030870	H74939171030810	3	80	5	1.7	24	2431
H74939171031040	H74939171031070	H74939171031010	3	100	5	1.7	24	2431
—	H74939171031270	H74939171031210	3	120	5	1.7	24	2431
—	H74939171031570	H74939171031510	3	150	5	1.7	24	2431
—	H74939171031870	H74939171031810	3	180	5	1.7	24	2431
—	H74939171032070	H74939171032010	3	200	5	1.7	24	2431
H74939171040240	H74939171040270	H74939171040210	4	20	5	1.7	24	2431
H74939171040340	H74939171040370	H74939171040310	4	30	5	1.7	24	2431
H74939171040440	H74939171040470	H74939171040410	4	40	5	1.7	24	2431
H74939171040640	H74939171040670	H74939171040610	4	60	5	1.7	24	2431
H74939171040840	H74939171040870	H74939171040810	4	80	5	1.7	24	2431
H74939171041040	H74939171041070	H74939171041010	4	100	5	1.7	24	2431
—	H74939171041270	H74939171041210	4	120	5	1.7	24	2431
—	H74939171041570	H74939171041510	4	150	5	1.7	24	2431
—	H74939171041870	H74939171041810	4	180	5	1.7	24	2431
—	H74939171042070	H74939171042010	4	200	5	1.7	24	2431
H74939171050240	H74939171050270	H74939171050210	5	20	5	1.7	24	2431
H74939171050340	H74939171050370	H74939171050310	5	30	5	1.7	24	2431
H74939171050440	H74939171050470	H74939171050410	5	40	5	1.7	24	2431
H74939171050640	H74939171050670	H74939171050610	5	60	5	1.7	24	2431
H74939171050840	H74939171050870	H74939171050810	5	80	5	1.7	24	2431
H74939171051040	H74939171051070	H74939171051010	5	100	5	1.7	24	2431
H74939171051240	H74939171051270	H74939171051210	5	120	5	1.7	24	2431
—	H74939171051570	H74939171051510	5	150	5	1.7	24	2431
—	H74939171051870	H74939171051810	5	180	5	1.7	24	2431
—	H74939171052070	H74939171052010	5	200	5	1.7	24	2431
H74939171060240	H74939171060270	H74939171060210	6	20	5	1.7	24	2431
H74939171060340	H74939171060370	H74939171060310	6	30	5	1.7	24	2431
H74939171060440	H74939171060470	H74939171060410	6	40	5	1.7	24	2431
H74939171060640	H74939171060670	H74939171060610	6	60	5	1.7	24	2431
H74939171060840	H74939171060870	H74939171060810	6	80	5	1.7	24	2431
H74939171061040	H74939171061070	H74939171061010	6	100	5	1.7	24	2431
H74939171061240	H74939171061270	H74939171061210	6	120	5	1.7	22	2229
—	H74939171061570	H74939171061510	6	150	5	1.7	22	2229
—	H74939171061870	H74939171061810	6	180	5	1.7	22	2229
—	H74939171062070	H74939171062010	6	200	5	1.7	22	2229
H74939171070240	H74939171070270	H74939171070210	7	20	5	1.7	20	2027
H74939171070340	H74939171070370	H74939171070310	7	30	5	1.7	20	2027
H74939171070440	H74939171070470	H74939171070410	7	40	5	1.7	20	2027
H74939171070640	H74939171070670	H74939171070610	7	60	5	1.7	20	2027
H74939171070840	H74939171070870	H74939171070810	7	80	5	1.7	20	2027
H74939171071040	H74939171071070	H74939171071010	7	100	5	1.7	20	2027
—	H74939171071270	H74939171071210	7	120	5	1.7	18	1823
—	H74939171071570	H74939171071510	7	150	5	1.7	18	1823
—	H74939171071870	H74939171071810	7	180	5	1.7	18	1823
—	H74939171072070	H74939171072010	7	200	5	1.7	18	1823
H74939171080240	H74939171080270	H74939171080210	8	20	6	2.0	20	2027
H74939171080340	H74939171080370	H74939171080310	8	30	6	2.0	20	2027
H74939171080440	H74939171080470	H74939171080410	8	40	6	2.0	20	2027
H74939171080640	H74939171080670	H74939171080610	8	60	6	2.0	20	2027
H74939171080840	H74939171080870	H74939171080810	8	80	6	2.0	20	2027
H74939171081040	H74939171081070	H74939171081010	8	100	6	2.0	20	2027
—	H74939171081570	H74939171081510	8	150	6	2.0	18	1823
—	H74939171081870	H74939171081810	8	180	6	2.0	18	1823
—	H74939171082070	H74939171082010	8	200	6	2.0	18	1823
H74939171090240	H74939171090270	H74939171090210	9	20	6	2.0	18	1823
H74939171090340	H74939171090370	H74939171090310	9	30	6	2.0	18	1823
H74939171090440	H74939171090470	H74939171090410	9	40	6	2.0	18	1823
H74939171090640	H74939171090670	H74939171090610	9	60	6	2.0	18	1823
H74939171090840	H74939171090870	H74939171090810	9	80	6	2.0	18	1823
H74939171100240	H74939171100270	H74939171100210	10	20	6	2.0	14	1419
H74939171100340	H74939171100370	H74939171100310	10	30	6	2.0	14	1419
H74939171100440	H74939171100470	H74939171100410	10	40	6	2.0	14	1419
H74939171100640	H74939171100670	H74939171100610	10	60	6	2.0	14	1419
H74939171100840	H74939171100870	H74939171100810	10	80	6	2.0	14	1419
H74939171120240	H74939171120270	H74939171120210	12	20	7	2.3	14	1419
H74939171120340	H74939171120370	H74939171120310	12	30	7	2.3	14	1419
H74939171120440	H74939171120470	H74939171120410	12	40	7	2.3	14	1419
H74939171120640	H74939171120670	H74939171120610	12	60	7	2.3	14	1419
H74939171120840	H74939171120870	H74939171120810	12	80	7	2.3	14	1419



2.3 0.035" (0.89 mm) Guidewire Compatible

Wanda™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter con balón de dilatación

Shaft Length 60 cm	Order Number Shaft Length 80 cm	Shaft Length 120 cm	Shaft Length 135 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
—	H965SCH505000	—	H965SCH505240	3	20	5	1.67	16	1621
—	H965SCH505010	—	H965SCH505250	3	40	5	1.67	16	1621
—	H965SCH505020	H965SCH505260	—	3	80	5	1.67	16	1621
—	H965SCH505030	—	H965SCH505270	4	20	5	1.67	16	1621
—	H965SCH505040	—	H965SCH505280	4	40	5	1.67	16	1621
—	H965SCH505050	H965SCH505290	—	4	80	5	1.67	16	1621
—	H965SCH505060	—	H965SCH505300	5	20	5	1.67	16	1621
—	H965SCH505070	—	H965SCH505310	5	40	5	1.67	16	1621
—	H965SCH505080	H965SCH505320	—	5	80	5	1.67	16	1621
—	H965SCH505090	—	H965SCH505330	6	20	5	1.67	15	1520
—	H965SCH505100	—	H965SCH505340	6	40	5	1.67	15	1520
—	H965SCH505110	H965SCH505350	—	6	80	5	1.67	15	1520
—	H965SCH505120	—	H965SCH505360	7	20	6	2.00	15	1520
—	H965SCH505130	—	H965SCH505370	7	40	6	2.00	15	1520
—	H965SCH505140	H965SCH505380	—	7	80	6	2.00	15	1520
—	H965SCH505150	—	H965SCH505390	8	20	6	2.00	15	1520
—	H965SCH505160	—	H965SCH505400	8	40	6	2.00	15	1520
—	H965SCH505170	H965SCH505410	—	8	80	6	2.00	15	1520
H965SCH505420	H965SCH505180	—	—	9	20	7	2.00	13	1317
H965SCH505430	H965SCH505190	—	—	9	40	7	2.00	13	1317
H965SCH505440	H965SCH505200	—	—	10	20	7	2.3	14	1419
H965SCH505450	H965SCH505210	—	—	10	40	7	2.3	14	1419
H965SCH505460	H965SCH505220	—	—	12	20	7	2.3	11	1115
H965SCH505470	H965SCH505230	—	—	12	40	7	2.3	11	1115

Recommended Guidewire 0.035" (0.89 mm), PTFE coated

XXL™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéters de dilatation à ballonnet / Catéter con balón de dilatación

Order Number Shaft Length 75 cm	Shaft Length 120 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath Compatibility (mm)	Rated Burst Pressure (ATM)	Rated Burst Pressure (kPa)
M001145080*	M001145090*	12	20	7	2.33	8	810
M001145100*	M001145110*	12	40	7	2.33	8	810
M001145300*	—	12	60	7	2.33	8	810
M001145120*	M001145130*	14	20	7	2.33	8	810
M001145140*	M001145150*	14	40	7	2.33	8	810
M001145320*	—	14	60	7	2.33	8	810
M001145500**	M001145510**	16	20	N/A	2.67	5	507
M001145520**	M001145530**	16	40	N/A	2.67	5	507
M001145540**	—	16	60	N/A	2.67	5	507
M001145550**	M001145560**	18	20	N/A	2.67	5	507
M001145570**	M001145580**	18	40	N/A	2.67	5	507
M001145590**	—	18	60	N/A	2.67	5	507

Recommended Introducer Sheaths 7F (2.33 mm) for 12 and 14 mm diameter balloons, 8F (2.67 mm) for 16 and 18 mm diameter balloons

Specifications for all XXL™ balloon dilatation catheters: Maximum guidewire 0.035" (0.89 mm) / Catheter Size 5.8F (1.93 mm) / Tip ID 0.037" (0.94 mm)

*Indicated for peripheral PTA

**Indicated for esophageal dilatation

3 Specialty Balloons

Peripheral Cutting Balloon™ Microsurgical Dilatation Device
Equalizer™ Balloon Dilatation Catheter



Peripheral Cutting Balloon™

Microsurgical Dilatation Device / Mikrochirurgisches Dilatationsgerät / Dispositif de dilatation microchirurgicale
Dispositivo de dilatación microquirúrgica

0.014" (0.36 mm), 1.5 cm Balloon Length - **Monorail™ System Flextome™***

Order Number Shaft Length 140 cm	Balloon Diameter (mm)	Balloon Length (cm)	Sheath Compatibility (F)	Sheath (mm)	Atherotomes	Rated Burst Pressure (ATM) (kPa)	
M001BPM2015140FO	2	15	6	2.00	3	12	1216
M001BPM2515140FO	2.5	15	6	2.00	3	12	1216
M001BPM3015140FO	3	15	6	2.00	3	12	1216
M001BPM3515140FO	3.5	15	6	2.00	4	12	1216
M001BPM4015140FO	4	15	6	2.00	4	12	1216

* Consult your local representative for availability in your country

0.018" (0.46 mm), 2.0 cm Balloon Length - **Over-The-Wire System**

Shaft Length 50 cm	Order Number Shaft Length 90 cm	Shaft Length 135 cm	Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Sheath (mm)	Rated Burst Pressure (ATM) (kPa)	
M001BP50520B0	M001BP90520B0	M001BP35520B0	5	20	7	2.33	10	1013
M001BP50620B0	M001BP90620B0	M001BP35620B0	6	20	7	2.33	10	1013
M001BP50720B0	M001BP90720B0	M001BP35720B0	7	20	7	2.33	10	1013
M001BP50820B0	M001BP90820B0	M001BP35820B0	8	20	7	2.33	10	1013

Maximum compatible guidewire 0.020" (0.50 mm)

Equalizer™

Balloon Dilatation Catheter / Ballondilatationskatheter / Cathéter de dilatation à ballonnet / Catéter balón de dilatación

Order Number	Inflated Balloon Diameter (mm)	Usable Length (cm)	Recommended Balloon Inflation Volume (ml)	Introducer Size (F)	Introducer Size (mm)
M001171050	20	65	4.8	14	4.0
M001171060	20	100	4.8	14	4.0
M001171070	27	65	10.0	14	4.7
M001171080	27	100	10.0	14	4.7
M001171090	33	65	19.0	16	5.3
M001171100	33	100	19.0	16	5.3
M001171110	40	65	34.5	N/A	N/A
M001171120	40	100	34.5	N/A	N/A

Specifications for all Equalizer™ Balloon Catheters:
Catheter Size 7F (2.3 mm); Number of Lumen 2; Minimum Tip ID 0.039" (0.99 mm); Recommended Guidewire 0.038" (0.97 mm)

4 CTO Devices

OffRoad™ Re-Entry Catheter System
TruePath™ CTO Device




OffRoad™

Re-entry Catheter System / Re-entry-Kathetersystem / Système de cathéter de ré-entrée / Sistema de catéter de reentrada

- OffRoad™: the industry's most elegant dual-component solution for challenging lesions. Re-entry is intuitive as the 5.4 mm articulating Positioning Balloon directs the system, creating a path for the Micro-Catheter Lancet to re-enter the vessel with precision
- OffRoad™: die elegante Zwei-Komponenten-Lösung für anspruchsvolle Läsionen. Der 5,4 mm große, abwinkelbare Platzierungsballon steuert das System und bildet einen Kanal für die Mikrokatheter-Nadel: für einfaches und präzises Re-entry
- OffRoad™: la solution à deux composants la plus élégante pour les lésions complexes. La ré-entrée est intuitive comme le ballonnet de positionnement de 5,4 mm est articulé et dirige le système, créant alors un chemin d'accès pour que le microcathéter biseauté puisse ré-entrer dans le vaisseau avec précision et facilité

- OffRoad™: la solución de doble componente más elegante de la industria para lesiones complejas. La re-entrada es intuitiva debido a que el balón de posicionamiento articulado de 5,4 mm dirige el sistema y crea una ruta para que el microcathéter Lanceta reentre en el vaso con precisión y facilidad



Compatibility Chart for OffRoad™ Re-entry Catheter–Short (70cm)

Order Number		Length (cm)	Catheter Shaft OD (F)	Catheter Shaft OD (mm)	Nominal Balloon Diameter (mm)	Nominal Inflation Pressure (ATM)	Nominal Inflation Pressure (kPa)	Maximum Inflation Pressure (ATM)	Maximum Inflation Pressure (kPa)	Sheath Diameter (F)	Sheath Diameter (mm)	Recommended Guidewire OD (in)	Recommended Guidewire OD (mm)
H74939202070540	Positioning Balloon Catheter	70	5	1.75	5.4	2	203	3.25	329	6	2	0.035	0.89
	Micro- Catheter Lancet	100	2.4	0.85	N/A	N/A	N/A	N/A	N/A	N/A	N/A	0.014*	0.36*

* Non-Coated. Do not use a 0.014 in (0.36 mm) guidewire that has a polymer (plastic) cover/sleeve with the Micro-Catheter Lancet. If using a 0.014 in (0.36 mm) PTFE coated guidewire, do not manipulate the PTFE coated part of the wire (i.e. back-and-forth movements) beyond the tip of the Micro-Catheter Lancet. Using an inappropriate guidewire may result in damage to the guidewire such as abrasion of the coating, release of polymer fragments or separation of the guidewire.

Compatibility Chart for OffRoad™ Re-entry Catheter–Long (100cm)

Order Number		Length (cm)	Catheter Shaft OD (F)	Catheter Shaft OD (mm)	Nominal Balloon Diameter (mm)	Nominal Inflation Pressure (ATM)	Nominal Inflation Pressure (kPa)	Maximum Inflation Pressure (ATM)	Maximum Inflation Pressure (kPa)	Sheath Diameter (F)	Sheath Diameter (mm)	Recommended Guidewire OD (in)	Recommended Guidewire OD (mm)
H74939202100540	Positioning Balloon Catheter	100	5	1.75	5.4	2	203	3.25	329	6	2	0.035	0.89
	Micro- Catheter Lancet	109	2.4	0.85	N/A	N/A	N/A	N/A	N/A	N/A	N/A	0.014*	0.36*

* Non-Coated. Do not use a 0.014 in (0.36 mm) guidewire that has a polymer (plastic) cover/sleeve with the Micro-Catheter Lancet. If using a 0.014 in (0.36 mm) PTFE coated guidewire, do not manipulate the PTFE coated part of the wire (i.e. back-and-forth movements) beyond the tip of the Micro-Catheter Lancet. Using an inappropriate guidewire may result in damage to the guidewire such as abrasion of the coating, release of polymer fragments or separation of the guidewire.


Recommended Guidewire: Platinum PLUS™ Guidewire

Order Number	Length (cm)	Diameter (in)	Diameter (mm)
M001467400	180	0.014	0.36
M001467120	300	0.014	0.36

TruePath™

CTO Device / CTO Device / CTO Device / CTO Device

- TruePath™: a powerful new solution for intraluminal treatment of chronic total occlusions. At 0.018"/0.46 mm in diameter, the guidewire-like profile is nearly half the size of competitive offerings, designed for optimal steering and crossing. Once positioned, the self-rotating tip works effortlessly
 - TruePath™: eine leistungsstarke neue Lösung für die intraluminale Behandlung chronischer Totalverschlüsse. Mit einem Durchmesser von 0,018" / 0,46 ist das Führungsdraht-ähnliche Profil nur halb so groß wie von Wettbewerbsprodukten und sorgt für optimale Steuerung und Passierbarkeit. Die Spitze dreht sich nach der Positionierung selbst – ohne weiteren Kraftaufwand
- TruePath™ : une nouvelle solution, optimale pour le traitement intraluminal des occlusions totales chroniques. Avec un diamètre de 0,018" / 0,46 mm, le profil de ce pseudo-guide est quasiment 50 % plus petit que les offres de la concurrence ; il est conçu pour une orientation et un franchissement optimaux. Une fois en place, l'extrémité autorotative travaille sans effort
 - TruePath™: una nueva y poderosa solución para el tratamiento intraluminal de oclusiones totales crónicas. Con un diámetro de 0,018" / 0,46 mm, el perfil similar a una guía tiene casi la mitad de tamaño que sus competidores y está diseñado para ofrecer un direccionamiento y cruce óptimos. Una vez posicionada, la punta autogiratoria funciona sin esfuerzos



TruePath™ CTO Device

Order Number	Description
H74939208181650	TruePath CTO Device
H74939208001601	TruePath Extension Wires (5/box)

TruePath™ Compatibility

Size OD (in)	Length (mm)	Length (cm)	Sheath Compatibility (in) (mm)		Min. Guide Catheter ID (in)	Catheter ID (mm)	Tip Profile (in) (mm)	
0.018	0.46	165	0.018 compatible catheters*		0.014	0.36	0.017	0.43

* With usable lengths of ≤ 135 cm

5 Catheters

5.1 Diagnostic Catheter

***Imager II™** Angiographic Catheter*




5.1 Diagnostic Catheter

Imager™ II


Angiographic Catheter Tip Shapes/ Angiographie-Katheter - Spitzenkonfigurationen

Formes des extrémités des cathéters d'angiographie / Formas de punta de los catéteres angiográficos


Flush Catheters



Selective Catheters Cerebral



Selective Catheters Visceral



5.1 Diagnostic Catheter

Imager™ II

Angiographic Catheter / Angiographie-Katheter / Cathéter d'angiographie / Catéter angiográfico

FLUSH CATHETERS - Tennis Racquet Catheters

Order Number	French Size (F) (mm)		Tip Shape	Length (cm)	Compatible Guidewire (in)	Guidewire (mm)	Maximum Flow Rate (ml/sec)
M001315001	4	1.38	Tennis Racquet	65	0.035	0.89	20
M001315011	4	1.38	Tennis Racquet	90	0.035	0.89	17
M001315021	4	1.38	Tennis Racquet	100	0.035	0.89	17
M001315031	5	1.70	Tennis Racquet	65	0.035	0.89	32
M001315041	5	1.70	Tennis Racquet	90	0.035	0.89	29
M001315051	5	1.70	Tennis Racquet	100	0.035	0.89	26
M001315061	5	1.70	Tennis Racquet	65	0.038	0.97	32
M001315071	5	1.70	Tennis Racquet	90	0.038	0.97	29
M001315081	5	1.70	Tennis Racquet	100	0.038	0.97	26

Packaged 5/box

FLUSH CATHETERS - Pigtail / Straight High Flow Catheters

Order Number	French Size (F) (mm)		Tip Shape	Length (cm)	Compatible Guidewire (in)	Guidewire (mm)	Maximum Flow Rate (ml/sec)
M001315091	4	1.38	Pigtail	65	0.035	0.89	20
M001315101	4	1.38	Straight	65	0.035	0.89	20
M001315111	4	1.38	Pigtail	90	0.035	0.89	17
M001315121	4	1.38	Straight	90	0.035	0.89	17
M001315131	5	1.70	Pigtail	65	0.038	0.97	32
M001315141	5	1.70	Pigtail	90	0.038	0.97	29
M001315151	5	1.70	Pigtail	100	0.038	0.97	26
M001315161	5	1.70	Straight	65	0.038	0.97	32
M001315171	5	1.70	Straight	90	0.038	0.97	29
M001315181	5	1.70	Straight	100	0.038	0.97	26
M001315191	5	1.70	Pigtail	65	0.035	0.89	32
M001315201	5	1.70	Straight	65	0.035	0.89	32
M001315211	5	1.70	Pigtail	90	0.035	0.89	29
M001315221	5	1.70	Straight	90	0.035	0.89	29
M001315231	5	1.70	Pigtail	100	0.035	0.89	26
M001315241	5	1.70	Straight	100	0.035	0.89	26

Packaged 5/box

FLUSH CATHETERS - Contralateral High Flow Catheters

Order Number	French Size (F) (mm)		Tip Shape	Length (cm)	Compatible Guidewire (in)	Guidewire (mm)	Maximum Flow Rate (ml/sec)
M001315251	4	1.38	Contralateral	65	0.035	0.89	20
M001315261	4	1.38	Contralateral	90	0.035	0.89	17
M001315271	4	1.38	Contralateral	100	0.035	0.89	17
M001315281	5	1.70	Contralateral	65	0.035	0.89	32
M001315291	5	1.70	Contralateral	90	0.035	0.89	29
M001315301	5	1.70	Contralateral	100	0.035	0.89	26
M001315311	5	1.70	Contralateral	65	0.038	0.97	32
M001315321	5	1.70	Contralateral	90	0.038	0.97	29
M001315331	5	1.70	Contralateral	100	0.038	0.97	26

Packaged 5/box



5.1 Diagnostic Catheter

SELECTIVE CATHETERS CEREBRAL - Braided Torque Cerebral Selective Catheters

Order Number	French Size (F)	Size (mm)	Tip Shape	Length (cm)	Compatible Guidewire (in)	Guidewire (mm)	Maximum Flow Rate (ml/sec)
M001316041	4	1.38	BARN	100	0.035	0.89	10
M001316061	4	1.38	BERN	100	0.035	0.89	10
M001316071	4	1.38	BERN	40	0.035	0.89	10
M001316091	4	1.38	BERN	65	0.035	0.89	10
M001316121	4	1.38	WEIN	100	0.035	0.89	10
M001316131	4	1.38	H1	65	0.035	0.89	10
M001316151	4	1.38	H1	100	0.035	0.89	10
M001316191	4	1.38	H1H	100	0.035	0.89	10
M001316211	4	1.38	JB1	65	0.035	0.89	10
M001316231	4	1.38	JB1	100	0.035	0.89	10
M001316251	4	1.38	JB2	100	0.035	0.89	10
M001316311	4	1.38	SIM1	65	0.035	0.89	12
M001316331	4	1.38	SIM1	100	0.035	0.89	12
M001316351	4	1.38	SIM2	100	0.035	0.89	12
M001316371	4	1.38	SIM3	100	0.035	0.89	12
M001314041	5	1.70	BARN	100	0.038	0.97	18
M001314051	5	1.70	BERN	100	0.038	0.97	18
M001314061	5	1.70	BERN	100	0.035	0.89	18
M001314071	5	1.70	BERN	40	0.035	0.89	20
M001314081	5	1.70	BERN	40	0.038	0.97	20
M001314091	5	1.70	BERN	65	0.035	0.89	20
M001314101	5	1.70	BERN	65	0.038	0.97	20
M001314111	5	1.70	WEIN	100	0.038	0.97	12
M001314121	5	1.70	WEIN	100	0.035	0.89	12
M001314131	5	1.70	H1	65	0.038	0.97	20
M001314141	5	1.70	H1	100	0.038	0.97	12
M001314151	5	1.70	H1	100	0.035	0.89	12
M001314161	5	1.70	H3	100	0.038	0.97	12
M001314171	5	1.70	H3	100	0.035	0.89	12
M001314181	5	1.70	H1H	100	0.038	0.97	12
M001314191	5	1.70	H1H	100	0.035	0.89	12
M001314211	5	1.70	JB1	65	0.038	0.97	20
M001314221	5	1.70	JB1	100	0.038	0.97	12
M001314231	5	1.70	JB1	100	0.035	0.89	12
M001314241	5	1.70	JB2	100	0.038	0.97	12
M001314251	5	1.70	JB2	100	0.035	0.89	12
M001314261	5	1.70	JB3	100	0.038	0.97	12
M001314271	5	1.70	JB3	100	0.035	0.89	12
M001314291	5	1.70	MANI	100	0.035	0.89	10

Packaged 5/box

5.1 Diagnostic Catheter

SELECTIVE CATHETERS CEREBRAL - Braided Torque Cerebral Selective Catheters (continuation)

Order Number	French Size (F) (mm)		Tip Shape	Length (cm)	Compatible Guidewire (in) (mm)		Maximum Flow Rate (ml/sec)
M001314301	5	1.70	SIM1	65	0.038	0.97	20
M001314311	5	1.70	SIM1	65	0.035	0.89	20
M001314321	5	1.70	SIM1	100	0.038	0.97	20
M001314331	5	1.70	SIM1	100	0.035	0.89	20
M001314341	5	1.70	SIM2	100	0.038	0.97	20
M001314351	5	1.70	SIM2	100	0.035	0.89	20
M001314361	5	1.70	SIM3	100	0.038	0.97	20
M001314371	5	1.70	SIM3	100	0.035	0.89	20
M001314381	5	1.70	SIM4	100	0.038	0.97	20
M001314391	5	1.70	HN1	100	0.038	0.97	10
M001314401	5	1.70	HN1	100	0.035	0.89	10
M001314411	5	1.70	HN2	100	0.038	0.97	10
M001314421	5	1.70	HN2	100	0.035	0.89	10
M001314431	5	1.70	HN3	100	0.038	0.97	10
M001314441	5	1.70	HN3	100	0.035	0.89	10
M001314451	5	1.70	HN4	100	0.038	0.97	10
M001314461	5	1.70	HN4	100	0.035	0.89	10
M001314471	5	1.70	HN5	100	0.038	0.97	10
M001314481	5	1.70	HN5	100	0.035	0.89	10

Packaged 5/box

SELECTIVE CATHETERS CEREBRAL - Braided Torque Multi-Purpose Selective Catheters

Order Number	French Size (F) (mm)		Tip Shape	Length (cm)	Compatible Guidewire (in) (mm)		Maximum Flow Rate (ml/sec)
M001316021	4	1.38	ST	65	0.035	0.89	10
M001316031	4	1.38	ST	100	0.035	0.89	10
M001316511	4	1.38	GEN*	100	0.035	0.89	10
M001314001	5	1.70	ST	65	0.038	0.97	20
M001314011	5	1.70	ST	100	0.038	0.97	18
M001314021	5	1.70	ST	65	0.035	0.89	20
M001314031	5	1.70	ST	100	0.035	0.89	18
M001314501	5	1.70	GEN	100	0.035	0.89	18

Packaged 5/box

*Denotes two side holes; all others no side holes



5.1 Diagnostic Catheter

SELECTIVE CATHETERS VISCERAL - Braided Torque Visceral Selective Catheters

Order Number	French Size (F) (mm)		Tip Shape	Length (cm)	Compatible Guidewire (in) (mm)		Maximum Flow Rate (ml/sec)
M001316581	4	1.38	CONTRA 2	65	0.035	0.89	6
M001316631	4	1.38	C1	65	0.035	0.89	10
M001316651	4	1.38	C1*	65	0.035	0.89	10
M001316671	4	1.38	C2	65	0.035	0.89	10
M001316691	4	1.38	C2*	65	0.035	0.89	10
M001316711	4	1.38	C3	65	0.035	0.89	10
M001316751	4	1.38	TRAIN	65	0.035	0.89	10
M001316771	4	1.38	MIK	65	0.035	0.89	6
M001316791	4	1.38	HK1.0	65	0.035	0.89	6
M001316811	4	1.38	SHK0.8	65	0.035	0.89	6
M001316831	4	1.38	SHK1.0	65	0.035	0.89	6
M001316851	4	1.38	RC1	65	0.035	0.89	6
M001316871	4	1.38	RC2	65	0.035	0.89	6
M001316891	4	1.38	RDC	65	0.035	0.89	10
M001314521	5	1.70	CHGB	65	0.038	0.97	10
M001314531	5	1.70	CHGB	65	0.035	0.89	10
M001314541	5	1.70	CHGC	65	0.038	0.97	10
M001314551	5	1.70	CHGC	65	0.035	0.89	10
M001314561	5	1.70	CHG2.5	65	0.038	0.97	10
M001314571	5	1.70	CHG2.5	65	0.035	0.89	10
M001314581	5	1.70	CONTRA 2	65	0.038	0.97	10
M001314591	5	1.70	CONTRA 2	65	0.035	0.89	10
M001314621	5	1.70	C1	65	0.038	0.97	20
M001314631	5	1.70	C1	65	0.035	0.89	20
M001314641	5	1.70	C1*	65	0.038	0.97	20
M001314651	5	1.70	C1*	65	0.035	0.89	20
M001314661	5	1.70	C2	65	0.038	0.97	20
M001314671	5	1.70	C2	65	0.035	0.89	20
M001314681	5	1.70	C2*	65	0.038	0.97	20
M001314691	5	1.70	C2*	65	0.035	0.89	20
M001314701	5	1.70	C3	65	0.038	0.97	20
M001314711	5	1.70	C3	65	0.035	0.89	20
M001314721	5	1.70	C3*	65	0.038	0.97	20
M001314731	5	1.70	C3*	65	0.035	0.89	20
M001314741	5	1.70	VAN*	65	0.038	0.97	10
M001314751	5	1.70	TRAIN	65	0.038	0.97	20
M001314761	5	1.70	MIK	65	0.038	0.97	10
M001314771	5	1.70	MIK	65	0.035	0.89	10
M001314781	5	1.70	HK1.0	65	0.038	0.97	10
M001314791	5	1.70	HK1.0	65	0.035	0.89	10

Packaged 5/box

*Denotes two side holes; all others no side holes



5.1 Diagnostic Catheter

SELECTIVE CATHETERS VISCERAL - Braided Torque Visceral Selective Catheters (continuation)

Order Number	French Size (F) (mm)		Tip Shape	Length (cm)	Compatible Guidewire (in) (mm)		Maximum Flow Rate (ml/sec)
M001314801	5	1.70	SHK0.8	65	0.038	0.97	10
M001314811	5	1.70	SHK0.8	65	0.035	0.89	10
M001314821	5	1.70	SHK1.0	65	0.038	0.97	10
M001314831	5	1.70	SHK1.0	65	0.035	0.89	10
M001314841	5	1.70	RC1	65	0.038	0.97	10
M001314851	5	1.70	RC1	65	0.035	0.89	10
M001314861	5	1.70	RC2	65	0.038	0.97	10
M001314871	5	1.70	RC2	65	0.035	0.89	10
M001314881	5	1.70	RDC	65	0.038	0.97	20
M001314891	5	1.70	RDC	65	0.035	0.89	20
M001314901	5	1.70	RDC*	65	0.038	0.97	20
M001314911	5	1.70	RDC*	65	0.035	0.89	20

Packaged 5/box

*Denotes two side holes; all others no side holes



5 Catheters

5.2 Guiding Catheter

Mach1™ Peripheral Guide Catheter







5.2 Guiding Catheter

Mach1TM

Peripheral Guide Catheter Tip Shapes / Peripherer Führungskatheter - Spitzenkonfigurationen

Formes des extrémités des cathéters guides périphériques / Formas de punta de los catéteres guía periféricos



5.2 Guiding Catheter

Mach1™

Peripheral Guide Catheter / Peripherer Führungskatheter / Cathéter guide périphérique / Catéter guía periférico

OD/ID
6F (2.1 mm) / 0.070" (1.7 mm)

Order Number	Tip Shape	Length (cm)
M001196000	HS	55
M001196660	RDC	55
M001196960	RDC 1	55
M001196040	RE-S	55
M001196080	RE-SS	55
M001196120	RE-L	55
M001196160	C1	55
M001196200	C2	55
M001196230	MP1	55
M001196260	MP1	90
M001196300	MP2	55
M001196330	MP2	90
M001196450	ST	55
M001196480	ST	90
M001196530	CROSSOVER 1	55
M001196560	CROSSOVER 2	55
M001196620	LIMA	55
M001196700	BATES 1	55
M001196800	BATES 2	55
M001196850	BATES 3	55
M001196900	BATES 4	55
M001196380	40°	90

OD/ID
7F (2.4 mm) / 0.081" (2.0 mm)

Order Number	Tip Shape	Length (cm)
M001197000	HS	55
M001197660	RDC	55
M001197960	RDC 1	55
M001197040	RE-S	55
M001197080	RE-SS	55
M001197120	RE-L	55
M001197160	C1	55
M001197200	C2	55
M001197230	MP1	55
M001197260	MP1	90
M001197300	MP2	55
M001197330	MP2	90
M001197450	ST	55
M001197480	ST	90
M001197530	CROSSOVER 1	55
M001197560	CROSSOVER 2	55
M001197620	LIMA	55
M001197700	BATES 1	55
M001197800	BATES 2	55
M001197850	BATES 3	55
M001197900	BATES 4	55
M001197380	40°	90

OD/ID
8F (2.7 mm) / 0.091" (2.3 mm)

Order Number	Tip Shape	Length (cm)
M001198000	HS	55
M001198660	RDC	55
M001198960	RDC 1	55
M001198040	RE-S	55
M001198080	RE-SS	55
M001198120	RE-L	55
M001198160	C1	55
M001198200	C2	55
M001198230	MP1	55
M001198260	MP1	90
M001198300	MP2	55
M001198330	MP2	90
M001198450	ST	55
M001198480	ST	90
M001198530	CROSSOVER 1	55
M001198560	CROSSOVER 2	55
M001198620	LIMA	55
M001198700	BATES 1	55
M001198800	BATES 2	55
M001198380	40°	90

Packaged 1/box

5 Catheters

5.3 Support Catheter

Rubicon™ Support Catheter




Rubicon™

Support Catheter / Supportkatheter / Cathéter de support / Catéter de soporte

- Enhances lesion access and entry
Facilitates crossing of tight occlusions & challenging vasculature
Enables delivery with and exchange of 0.014" (0.36mm), 0.018" (0.46mm) & 0.035" (0.89mm) guidewires
- Améliore l'accès et l'entrée à la lésion
Facilite le franchissement d'occlusions étroites et de vaisseaux tortueux
Permet la mise en place et l'échange de guides 0.014" (0.36mm), 0.018" (0.46mm) & 0.035" (0.89mm)
- Besserer Läsionszu- und eingang
Für die Passage enger Verschlüsse und anspruchsvoller Gefäße
Ermöglicht das Platzieren und Wechseln von 0.014" (0.36mm), 0.018" (0.46mm) & 0.035" (0.89mm) Führungsdrähten

- Mejora el acceso y la entrada en la lesión
Facilita el cruce de oclusiones estrechas y la vasculatura compleja
Permite la liberación a través de guías de 0.014" (0.36mm), 0.018" (0.46mm) & 0.035" (0.89mm) y su intercambio



Order Number	Description	Length (cm)
H74939212014131	Rubicon 14	135
H74939212014151	Rubicon 14	150

Packaged 5/box

Order Number	Description	Length (cm)
H74939239018091	Rubicon 18	90cm
H74939239018131	Rubicon 18	135cm
H74939239018151	Rubicon 18	150cm

Packaged 5/box

Order Number	Description	Length (cm)
H74939240035061	Rubicon 35	65cm
H74939240035091	Rubicon 35	90cm
H74939240035131	Rubicon 35	135cm
H74939240035151	Rubicon 35	150cm

Packaged 5/box

6 Accessories

Medikit Supersheath™ Introducer Sheath

Medikit Supersheath™ RO Introducer Sheath

Medikit Supersheath™ XL Introducer Sheath

Encore™ 26 Inflation Device

Gateway™ Plus Y-Adaptor

Vascular Dilators

Femoral Entry Needles

Torque Vises

Connecting Tubes

Floswitch™ HP High Pressure Flow Control Device



Medikit Supersheath™

Introducer Sheath / Einführschleuse / Introducteur / Introductor

- The latest generation of the Medikit Supersheath™ with smooth transitions to facilitate insertion of the sheath. Silicone-coated valve and smooth inner surface promote device passage. Rotating suture wing allows for quick sheath repositioning, innovative dilator/hub twist-lock feature for security and safety during insertion

Each set supplied with 0.035" (0.89 mm) J-tip guidewire (45 cm length with 7 cm and 11 cm introducers; 80 cm length with 25 cm introducers)

Color-coded hub to facilitate French size identification

4F (1.3 mm) Sheath is compatible with a 0.035" (0.89 mm) guidewire. All other sizes are compatible with 0.038" (0.97 mm) guidewires

- Die jüngste Generation der Medikit Supersheath™ mit glatten Übergängen erleichtert die Einführung der Schleuse. Das silikonbeschichtete Ventil und die glatte Innenfläche erleichtern die Passage des Systems. Der drehbare Nahtflügel erlaubt eine schnelle Neupositionierung der Einführschleuse, die innovative Dilator/Twist-Lock-Funktion für den Schaft gewährleistet die Sicherheit während der Einführung

Jedes Set wird mit einem Führungsdraht von 0,035 Zoll (0,889 mm) mit J-Spitze geliefert (45 cm Länge bei 7 cm und 11 cm Einführschleusen; 80 cm Länge bei 25 cm Einführschleusen)

Der farbcodierte Schaftansatz erlaubt eine einfache Identifizierung der Größe

Die 4F (1,3 mm) Einführschleuse ist kompatibel mit einem Führungsdraht 0,035" (0,889 mm). Alle anderen Größen sind kompatibel mit Führungsdrähten 0,038" (0,965 mm)

- La dernière génération de Medikit Supersheath™ avec des transitions lisses pour faciliter et simplifier l'insertion de l'introducteur. La valve à revêtement en silicone et la surface interne lisse facilitent le passage du dispositif. L'ailette de suture pivotante permet le repositionnement rapide de l'introducteur, verrouillage innovant du dilateur/embase pour plus de sécurité au cours de l'insertion

Chaque kit est fourni avec un guide à extrémité en J de 0,035" (0,889 mm) (45 cm de longueur avec les introducteurs de 7 cm et de 11 cm ; 80 cm de longueur avec les introducteurs de 25 cm)

Embase à code couleur pour une identification simple des tailles Fr L'introducteur 4F (1,3 mm) est compatible avec un guide de 0,035" (0,889 mm). Toutes les autres tailles sont compatibles avec les guides de 0,038" (0,965 mm)

- La última generación de Medikit Supersheath™ con suaves transiciones para facilitar la inserción del introductor. La válvula recubierta de silicona y la superficie interna suave proporcionan el fácil tránsito del dispositivo. La rotación de la aleta de sutura ofrece un rápido reposicionamiento del introductor y una innovadora función del cono del dilador con cierre por torsión para una mayor seguridad durante la inserción

Cada set se suministra con una guía con punta en J de 0,035 pulgadas (0,889 mm) (45 cm de largo con introductores de 7 cm y 11 cm; 80 cm de largo con introductores de 25 cm)

Válvula con codificación por colores para garantizar una fácil identificación del tamaño French

El introductor de 4F (1,3 mm) es compatible con una guía de 0,035 pulgadas (0,889 mm). Los tamaños restantes son compatibles con guías de 0,038 pulgadas (0,965 mm)

Medikit Supersheath™ Set*

Sheath Length 7 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115704B1	4	1.3	Red
M00115705B1	5	1.7	Grey
M00115706B1	6	2.0	Green
M00115707B1	7	2.3	Orange
M00115708B1	8	2.7	Blue

Packaged 10/box

Sheath Length 11 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115710B1	4	1.3	Red
M00115711B1	5	1.7	Grey
M00115712B1	6	2.0	Green
M00115713B1	7	2.3	Orange
M00115714B1	8	2.7	Blue
M00115715B1	9	3.0	Black

Packaged 10/box

* Each set supplied with 0.035" (0.89 mm) J-tip guidewire (45 cm length, with 7 cm and 11 cm introducers; 80 cm length, with 25 cm introducer)

Medikit Supersheath™

Sheath Length 11 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115720B1	4	1.3	Red
M00115721B1	5	1.7	Grey
M00115722B1	6	2.0	Green
M00115723B1	7	2.3	Orange
M00115724B1	8	2.7	Blue
M00115725B1	9	3.0	Black

Packaged 10/box

Sheath Length 25 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115730B1	5	1.7	Grey
M00115731B1	6	2.0	Green
M00115732B1	7	2.3	Orange
M00115733B1	8	2.7	Blue
M00115734B1	9	3.0	Black

Packaged 10/box

Medikit Supersheath™ RO

Introducer Sheath / Einführschleuse / Introducteur / Introducutor

- Combining the features of the Medikit SuperSheath™ with a distal radiopaque marker to facilitate visualization of the distal tip 4F (1.3 mm) Sheath is compatible with a 0.035" (0.89 mm) guidewire. All other sizes are compatible with 0.038" (0.97 mm) guidewires

- Kombiniert die Funktionen des Medikit SuperSheath™ mit einem distalen, röntgendichten Marker zur leichteren Visualisierung der distalen Spitze

4F (1,3 mm) Einführschleuse, kompatibel mit einem Führungsdraht 0,035" (0,889 mm). Alle anderen Größen sind kompatibel mit Führungsdrähten 0,038" (0,965 mm)

- Associe les caractéristiques du Medikit SuperSheath™ à un marqueur radio-opaque distal pour faciliter la visualisation de l'extrémité distale

L'introducteur 4F (1,3 mm) est compatible avec un guide de 0,035" (0,889 mm). Toutes les autres tailles sont compatibles avec les guides de 0,038" (0,965 mm)

- Combinación de las características del Medikit SuperSheath™ con un marcador distal radiopaco para facilitar la visualización de la punta distal

El introductor de 4F (1,3 mm) es compatible con una guía de 0,035 pulgadas (0,889 mm). Los tamaños restantes son compatibles con guías de 0,038 pulgadas (0,965 mm)

Medikit Supersheath™ RO

Sheath Length 7 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115956B1	5	1.7	Grey
M00115957B1	6	2.0	Green
M00115958B1	7	2.3	Orange
M00115959B1	8	2.7	Blue
M00115955B1	9	3.0	Black

Packaged 10/box

Medikit Supersheath™ RO

Sheath Length 25 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115970B1	4	1.3	Red
M00115971B1	5	1.7	Grey
M00115965B1	6	2.0	Green
M00115966B1	7	2.3	Orange
M00115967B1	8	2.7	Blue
M00115968B1	9	3.0	Black

Packaged 10/box

Sheath Length 11 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115960B1	4	1.3	Red
M00115961B1	5	1.7	Grey
M00115962B1	6	2.0	Green
M00115963B1	7	2.3	Orange
M00115964B1	8	2.7	Blue
M00115969B1	9	3.0	Black

Packaged 10/box



Medikit Supersheath™ XL

Introducer Sheath / Einführschleuse / Introducteur / Introducutor

- Slide-free one-way Silicone hemostasis valve designed to allow smooth passage of devices and leak-free performance; smooth dilator/sheath transition to promote atraumatic insertion; locking hub design to protect against movement during percutaneous insertion

Each set includes short 0.035" (0.89 mm) guidewires

All Medikit Supersheath™ XL sheaths are compatible with 0.038" (0.97 mm) guidewires

- Silikonbeschichtetes Hämostaseventil für gute Passageeigenschaften ohne Leckagen; glatter Übergang zwischen Dilator und Einführschleuse zur leichteren atraumatischen Einführung; Locking-Hub-Design schützt vor Verrutschen während der perkutanen Einführung

Zu jedem Satz gehören kurze Führungsdrähte 0,035" (0,889 mm).

Alle Medikit Supersheath™ XL-Einführschleusen sind kompatibel mit Führungsdrähten 0,038" (0,965 mm)

- Valve hémostatique en silicone à sens unique, sans glissement, conçue pour faciliter le passage des dispositifs et permettre des performances sans fuite ; transition dilateur/introducteur lisse pour favoriser l'insertion atraumatique ; embase verrouillable pour empêcher les mouvements pendant l'insertion percutanée

Chaque kit comprend des guides courts de 0,035" (0,889 mm)

Tous les introducteurs Medikit Supersheath™ XL sont compatibles avec des guides de 0,038" (0,965 mm)

- Válvula hemostática de silicona diseñada para permitir el paso progresivo de dispositivos y el rendimiento sin pérdidas; suave transición del dilatador/vaina para proporcionar la inserción atraumática; diseño de válvula de bloqueo para proteger frente al movimiento durante la inserción percutánea

Cada set incluye guías cortas de 0,035 pulgadas (0,889 mm)

Todos los introductores Medikit Supersheath™ XL son compatibles con guías de 0,038 pulgadas (0,965 mm)

Medikit Supersheath™ XL Set*

Sheath Length 11 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115716B1	10	3.3	Grey
M00115717B1	11	3.7	White
M00115718B1	12	4.0	Green
M00115741B1	14	4.7	Orange

Packaged 10/box

* Each set includes short 0.035" (0.89 mm) guidewire

Medikit Supersheath™ XL

Sheath Length 11 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115726B1	10	3.3	Grey
M00115727B1	11	3.7	White
M00115728B1	12	4.0	Green
M00115740B1	14	4.7	Orange

Packaged 10/box

Sheath Length 25 cm

Order Number	Size (F)	Size (mm)	Color Code
M00115735B1	10	3.3	Grey
M00115736B1	11	3.7	White
M00115737B1	12	4.0	Green
M00115739B1	14	4.7	Orange

Packaged 10/box



Encore™ 26

Inflation Device / Inflationssystem / Inflateur / Dispositivo de inflado

- Specifically designed for high pressure applications
Working range to 26 ATM (2634 kPa) to meet the requirements of high pressure applications
High pressure stopcock designed to prevent leakage during high pressure inflation throughout the working range
 - Speziell für Hochdruckanwendungen
Arbeitsbereich bis 26 ATM (2634 kPa) entsprechend den Anforderungen für Hochdruckanwendungen
Das Hochdruck-Absperrventil verhindert eine Leckage bei Inflation mit hohem Druck im gesamten Arbeitsbereich
- Conçu spécialement pour les applications haute pression
Plage de travail jusqu'à 26 ATM (2 634 kPa) pour répondre aux exigences des applications haute pression
Le robinet d'arrêt haute pression prévient les fuites pendant le gonflage à haute pression sur toute la plage de travail
 - Diseñado específicamente para aplicaciones de alta presión
Rango de trabajo de 26 ATM (2634 kPa) para cumplir con los requisitos de aplicaciones de alta presión
La llave de paso de alta presión previene el goteo durante el inflado de alta presión en el rango de trabajo

Order Number	Description	Units per box
M001151050	20cc Encore™ inflation device	1
M001151062	20cc Encore™ inflation device	5

Gateway™ Plus

Y-Adapter / Y-Adapter / Adaptateur en Y / Adaptador Y

- Hemostasis valve designed to provide tight seal around devices and stabilize guidewire position during over-the-wire injection
Large 0.118" (2.997 mm) I.D. facilitates device passage
Rotating connector positions the injection port in an optimal location throughout the procedure
 - Die Konstruktion des Hämostaseventils gewährleistet eine sichere Abdichtung des Systems und stabilisiert die Führungsdrahtposition während der Injektion
Großer Innendurchmesser von 0,118" (2,997 mm) erleichtert die Systempassage
Durch den drehbaren Verbinder lässt sich der Injektionsport während des gesamten Eingriffs an der optimalen Stelle positionieren
- Le design de la valve hémostatique assure une étanchéité parfaite autour des dispositifs et stabilise la position du guide pendant l'injection sur guide
Le diamètre interne de 0,118" (2,997 mm) facilite le passage des dispositifs
Le connecteur pivotant permet de toujours positionner l'orifice d'injection de façon optimale tout au long de la procédure
 - El diseño de la válvula hemostática proporciona un sellado estrecho alrededor de los dispositivos y estabiliza la posición de la guía durante la inyección sobre la guía.
El gran I.D. de 0,118 pulgadas (2.997 mm) facilita el paso del dispositivo
La rotación del conector sitúa el puerto de inyección en una ubicación óptima durante el procedimiento

Order Number	Description	Units per box
M001153223	Gateway™ Plus Y-Adapter	10

Vascular Dilators

Vascular Dilators / Gefäß-Dilatatoren / Dilateurs vasculaires / Dilatadores vasculares


- A full line of high quality dilators for vascular and fascial applications. Smooth end-rounding and gradual taper ease insertion through tissue. PTFE construction provides strength to resist buckling during advancement
- Une gamme complète de dilateurs de haute qualité pour les applications vasculaires et fasciales. L'effilage lisse progressif et arrondi facilite l'insertion à travers les tissus. La construction en PTFE offre la résistance nécessaire contre le gauchissement pendant la progression
- Komplette Produktlinie hochwertiger Dilatatoren für Gefäßanwendungen und fasziale Anwendungen. Glatte, abgerundete Enden und die konische Form erleichtern die Einführung durch das Gewebe. PTFE-Konstruktion mit ausreichender Festigkeit verhindert das Abknicken während des Vorschubs
- Una línea completa de dilatadores de alta calidad para aplicaciones vasculares y faciales. La suave punta y gradual conificación facilita la inserción a través del tejido. La construcción en PTFE proporciona resistencia para evitar la deformación durante el avance

Order Number	Description	Size (F)	Size (mm)	Recommended Guidewire (in)	Recommended Guidewire (mm)	Length (cm)
M001481491	TD-4/35	4	1.3	0.035	0.89	20
M001481511	TD-5/35	5	1.7	0.035	0.89	20
M001481551	TD-5/38	5	1.7	0.038	0.97	20
M001481521	TD-6/38	6	2.0	0.038	0.97	20
M001481531	TD-7/38	7	2.3	0.038	0.97	20
M001481541	TD-8/38	8	2.7	0.038	0.97	20

Femoral Entry Needles

Femoral Entry Needles / Femorale Einführnadeln / Aiguilles de ponction fémorale / Agujas de acceso femoral

- Femoral Entry Needles for vascular access
- Femorale Einführnadeln für den Gefäßzugang
- Aiguilles de pose fémorale pour l'accès vasculaire
- Agujas de acceso femoral para entrada vascular




Order Number	Description	Style	Gauge (G)	OD (mm)	ID (mm)	Guidewire (in)	Guidewire (mm)	Length (mm)
M001441631	SEN-18/038	Seldinger	18	1.270	0.84	0.038	0.97	73
M001441691	EN-18/038	Single piece	18	1.270	0.84	0.038	0.97	70
M001441601	SEN-19/035	Seldinger	19	1.067	0.69	0.035	0.89	73
M001441661	EN-19/035	Single piece	19	1.067	0.69	0.035	0.89	70

Packaged 10/box

Torque Vises

Torque Vises / Torque Vises / Torqueurs / Torque

- Adjustable size accomodates 0.025" (0.64 mm) to 0.038" (0.97 mm) guidewires. Ribbed surface promotes handling
- Anpassbare Größe für Führungsdrähte zwischen 0,025" (0,635 mm) und 0,038" (0,965 mm). Gerippte Oberfläche erleichtert die Handhabung
- La taille ajustable est compatible avec les guides de 0,025" (0,635 mm) à 0,038" (0,965 mm). La surface striée facilite la manipulation
- El tamaño ajustable se adapta a guías de 0,025 pulgadas (0,635 mm) a 0,038 pulgadas (0,965 mm). La superficie estriada proporciona un fácil manejo



Order Number	Description	Units per box
M001465501	Torque vises fit 0.025" (0.64 mm) to 0.038" (0.97 mm) wires	12

Connecting Tubes

Connecting Tubes / Anschlussschläuche / Raccords haute pression / Tubos de conexión

Low pressure connecting tube
PVC - female to male fixed luer lock - ID 2.26 mm

Order Number	Description	Length (cm)	Maximum Working Pressure (PSI) (kPa)	
M001203371	LPCT/25/500 PSI/34 BAR	25	500	3448
M001203381	LPCT/51/500 PSI/34 BAR	51	500	3448
M001203391	LPCT/76/500 PSI/34 BAR	76	500	3448
M001203401	LPCT/102/500 PSI/34 BAR	102	500	3448
M001203411	LPCT/122/500 PSI/34 BAR	122	500	3448
M001203421	LPCT/152/500 PSI/34 BAR	152	500	3448

Packaged 50/box

High pressure connecting tube with rotating luer
PVC - female to male rotating luer lock - ID 2.26 mm

Order Number	Description	Length (cm)	Maximum Working Pressure (PSI) (kPa)	
M001203251	HPCT/25/1000 PSI/69 BAR	25	1000	6895
M001203261	HPCT/51/1000 PSI/69 BAR	51	1000	6895
M001203271	HPCT/76/1000 PSI/69 BAR	76	1000	6895
M001203281	HPCT/102/1000 PSI/69 BAR	102	1000	6895
M001203291	HPCT/122/1000 PSI/69 BAR	122	1000	6895
M001203361	HPCT/152/1000 PSI/69 BAR	152	1000	6895

Packaged 50/box

Flexible high pressure connecting tube braided with rotating luer / braided polyurethane
– female to male rotating luer lock - ID 1.80 mm


Order Number	Description	Length (cm)	Maximum Working Pressure (PSI) (kPa)	
M001203311	HPCT/25/1200 PSI/84 BAR	25	1200	8274
M001203321	HPCT/51/1200 PSI/84 BAR	51	1200	8274
M001203331	HPCT/76/1200 PSI/84 BAR	76	1200	8274
M001203341	HPCT/102/1200 PSI/84 BAR	102	1200	8274
M001203351	HPCT/122/1200 PSI/84 BAR	122	1200	8274

Packaged 50/box

Floswitch™ HP

High Pressure Flow Control Device / Hochdruckdurchflussregler / Robinet 1 voie / Dispositivos de control de flujo de alta presión

- Sliding thumb switch allows simple, one-handed operation
Large lumen permits unrestricted flow in angiographic procedures
Design allows operation with or without guidewire in place, accomodates guidewires up to 0.045" (1.143 mm) and permits use as torque vise for guidewire manipulation
Transparent lumen enables bubble detection
High strength material withstands pressures up to 1050 PSI (7239 kPa) to accomodate high pressure injections
- Schieberändelschalter zur einfachen Bedienung mit einer Hand
Unbeschränkter Durchfluss durch großes Lumen für Angiografie-Eingriffe.
Durch die Konstruktion ist ein Einsatz mit oder ohne Führungsdraht möglich. Geeignet für Führungsdrähte bis 0,045" (1,143 mm) Durchmesser; zur Arbeit mit Führungsdrähten auch als Torque Vise einsetzbar
Einfache Erkennung von Bläschen durch transparentes Lumen
Extrem festes Material, beständig für Drücke bis 1.050 psi (7239 kPa) für Injektionen mit hohem Druck
- L'interrupteur à glissière permet une utilisation simple d'une seule main
La grande lumière permet un flux illimité dans les procédures d'angiographie
Le design permet une utilisation avec ou sans guide en place, est compatible avec les guides de jusqu'à 0,045" (1,143 mm) et permet d'utiliser le régulateur comme étau de torsion pour la manipulation du guide
Lumière transparente pour une détection aisée des bulles
Le matériau haute résistance résiste à des pressions de jusqu'à 1 050 PSI (7 239 kPa) et est compatible avec les injections haute pression
- El deslizamiento del botón con el pulgar permite una simple operación con una mano
El lumen grande permite un flujo ilimitado en los procedimientos angiográficos
Su diseño permite el funcionamiento con o sin guía colocada, acomoda guías de hasta 0,045 pulgadas (1,143 mm) y permite el uso de un torque para la manipulacion de la guía
El lumen transparente permite la fácil detección de burbujas
El material altamente resistente soporta presiones de hasta 1.050 PSI (7239 kPa) para permitir inyecciones de alta presión



Order Number	Description	Units per box
M001442001	One Way High Pressure Flow Control Device 1050 PSI (7239 kPa)	12
M001442011	One Way High Pressure Flow Control Device 1050 PSI (7239 kPa)	24

Floswitch™ is a registered trademark of Becton Dickinson


7 Stents

7.1 Balloon Expandable Stenting

***Express™ Vascular LD** Peripheral Stent System*

***Express™ Vascular SD** Monorail™ Premounted Stent System*

PROMUS Element™ Plus Drug Eluting Stent for Below The Knee treatment



7.1 Balloon Expandable Stenting

Express™ Vascular LD

Peripheral Stent System / Peripheres Stentsystem / Système de stent périphérique / Sistema de stent periférico

Short instrument - 75 cm catheter

Order Number	Stent Diameter at NP* (mm)	Stent Length (mm)	Max. Diameter (mm)	Balloon Diameter (mm)	Balloon Length (mm)	Balloon Nominal Pressure (ATM)	Balloon Nominal Pressure (kPa)	Balloon Rated Burst Pressure (ATM)	Balloon Rated Burst Pressure (kPa)	Minimum Introd. Sheath (F)	Minimum Introd. Sheath (mm)	Crossing Profile (mm/in)
H74938162520750	5	17	9	5	20	8	811	12	1216	6	2.16	1.93/0.076
H74938162530750	5	27	9	5	30	8	811	12	1216	6	2.16	1.96/0.077
H74938162540750	5	37	9	5	40	8	811	12	1216	6	2.16	1.93/0.076
H74938162560750	5	57	9	5	60	8	811	12	1216	6	2.16	1.96/0.077
H74938162620750	6	17	9	6	20	8	811	12	1216	6	2.16	1.96/0.077
H74938162630750	6	27	9	6	30	8	811	12	1216	6	2.16	2.03/0.080
H74938162640750	6	37	9	6	40	8	811	12	1216	6	2.16	1.93/0.076
H74938162660750	6	57	9	6	60	8	811	12	1216	6	2.16	1.98/0.078
H74938162720750	7	17	9	7	20	8	811	12	1216	6	2.16	1.98/0.078
H74938162730750	7	27	9	7	30	8	811	12	1216	6	2.16	2.06/0.081
H74938162740750	7	37	9	7	40	8	811	12	1216	6	2.16	2.03/0.080
H74938162760750	7	57	9	7	60	8	811	12	1216	6	2.16	2.00/0.079
H74938162820750	8	17	9	8	20	8	811	12	1216	6	2.16	2.06/0.081
H74938162830750	8	27	9	8	30	8	811	12	1216	6	2.16	2.08/0.082
H74938162840750	8	37	9	8	40	8	811	12	1216	6	2.16	2.08/0.082
H74938162860750	8	57	9	8	60	8	811	12	1216	7	2.51	2.18/0.086
H74938162920750	9	25	11	9	30	8	811	12	1216	7	2.51	2.34/0.092
H74938162940750	9	37	11	9	40	8	811	12	1216	7	2.51	2.34/0.092
H74938162960750	9	57	11	9	60	8	811	12	1216	7	2.51	2.34/0.092
H74938162102070	10	25	11	10	30	10	1013	12	1216	7	2.51	2.39/0.094
H74938162104070	10	37	11	10	40	10	1013	12	1216	7	2.51	2.34/0.092
H74938162106070	10	57	11	10	60	10	1013	12	1216	7	2.51	2.31/0.091

*Nominal Pressure

Long instrument - 135 cm catheter

Order Number	Stent Diameter at NP* (mm)	Stent Length (mm)	Max. Diameter (mm)	Balloon Diameter (mm)	Balloon Length (mm)	Balloon Nominal Pressure (ATM)	Balloon Nominal Pressure (kPa)	Balloon Rated Burst Pressure (ATM)	Balloon Rated Burst Pressure (kPa)	Minimum Introd. Sheath (F)	Minimum Introd. Sheath (mm)	Crossing Profile (mm/in)
H74938162520130	5	17	9	5	20	8	811	12	1216	6	2.16	1.93/0.076
H74938162530130	5	27	9	5	30	8	811	12	1216	6	2.16	1.96/0.077
H74938162540130	5	37	9	5	40	8	811	12	1216	6	2.16	1.93/0.076
H74938162560130	5	57	9	5	60	8	811	12	1216	6	2.16	1.96/0.077
H74938162620130	6	17	9	6	20	8	811	12	1216	6	2.16	1.96/0.077
H74938162630130	6	27	9	6	30	8	811	12	1216	6	2.16	2.03/0.080
H74938162640130	6	37	9	6	40	8	811	12	1216	6	2.16	1.93/0.076
H74938162660130	6	57	9	6	60	8	811	12	1216	6	2.16	1.98/0.078
H74938162720130	7	17	9	7	20	8	811	12	1216	6	2.16	1.98/0.078
H74938162730130	7	27	9	7	30	8	811	12	1216	6	2.16	2.06/0.081
H74938162740130	7	37	9	7	40	8	811	12	1216	6	2.16	2.03/0.080
H74938162760130	7	57	9	7	60	8	811	12	1216	6	2.16	2.00/0.079
H74938162820130	8	17	9	8	20	8	811	12	1216	6	2.16	2.06/0.081
H74938162830130	8	27	9	8	30	8	811	12	1216	6	2.16	2.08/0.082
H74938162840130	8	37	9	8	40	8	811	12	1216	6	2.16	2.08/0.082
H74938162860130	8	57	9	8	60	8	811	12	1216	7	2.51	2.18/0.086
H74938162920130	9	25	11	9	30	8	811	12	1216	7	2.51	2.34/0.092
H74938162940130	9	37	11	9	40	8	811	12	1216	7	2.51	2.34/0.092
H74938162960130	9	57	11	9	60	8	811	12	1216	7	2.51	2.34/0.092
H74938162120130	10	25	11	10	30	10	1013	12	1216	7	2.51	2.39/0.094
H74938162140130	10	37	11	10	40	10	1013	12	1216	7	2.51	2.34/0.092
H74938162160130	10	57	11	10	60	10	1013	12	1216	7	2.51	2.31/0.091

Recommended guidewire 0.035" (0.89 mm)

*Nominal Pressure

Ordering Information - 64
Peripheral Interventions

Boston
Scientific



7.1 Balloon Expandable Stenting

Express™ Vascular SD

Monorail™ Premounted Stent System / Vormontiertes Monorail™ Stentsystem / Système de stent prémonté Monorail™
Sistema de stent periférico

Order Number	Stent Diam. at NP* (mm)	Stent Length (mm)	Cath. Usable Length (mm)	Balloon Length (mm)	Balloon Nominal Pressure (ATM) (kPa)	Balloon Rated Burst Pressure (ATM) (kPa)	Minimum Introducer Sheath (F) (mm)	Recommended Guiding Catheter OD (F) ID (in)	Recommended Guiding Catheter OD (mm) ID (mm)	Crossing Profile (in) (mm)
H74937917415900	4.0	15	90	16	10 1013	14 1419	5 1.67	6/0.064	2.00/1.63	0.055 1.40
H74937918419900	4.0	19	90	20	10 1013	14 1419	5 1.67	6/0.064	2.00/1.63	0.055 1.40
H74937917515900	5.0	15	90	16	10 1013	14 1419	5 1.67	6/0.067	2.00/1.70	0.059 1.50
H74937918519900	5.0	19	90	20	10 1013	14 1419	5 1.67	6/0.067	2.00/1.70	0.060 1.52
H74937917614900	6.0	14	90	15	10 1013	14 1419	5 1.67	6/0.067	2.00/1.70	0.063 1.60
H74937918618900	6.0	18	90	19	10 1013	14 1419	5 1.67	6/0.067	2.00/1.70	0.063 1.60
H74937917715900	7.0	15	90	16	10 1013	14 1419	6 2.00	7/0.078	2.33/1.98	0.070 1.78
H74937918719900	7.0	19	90	20	10 1013	14 1419	6 2.00	7/0.078	2.33/1.98	0.070 1.78
H74937917415150	4.0	15	150	16	10 1013	14 1419	5 1.67	6/0.064	2.00/1.63	0.055 1.40
H74937918419150	4.0	19	150	20	10 1013	14 1419	5 1.67	6/0.064	2.00/1.63	0.055 1.40
H74937917515150	5.0	15	150	16	10 1013	14 1419	5 1.67	6/0.067	2.00/1.70	0.059 1.50
H74937918519150	5.0	19	150	20	10 1013	14 1419	5 1.67	6/0.067	2.00/1.70	0.060 1.52
H74937917614150	6.0	14	150	15	10 1013	14 1419	5 1.67	6/0.067	2.00/1.70	0.063 1.60
H74937918618150	6.0	18	150	19	10 1013	14 1419	5 1.67	6/0.067	2.00/1.70	0.063 1.60
H74937917715150	7.0	15	150	16	10 1013	14 1419	6 2.00	7/0.078	2.33/1.98	0.070 1.78
H74937918719150	7.0	19	150	20	10 1013	14 1419	6 2.00	7/0.078	2.33/1.98	0.070 1.78

Maximum OD compatible Guidewire 0.018" (0.46 mm)
*Nominal Pressure

7.1 Balloon Expandable Stenting

PROMUS Element™ Plus DES BTK **Over-The-Wire™**

Drug Eluting Stent for Below The Knee treatment

Order Number	Stent Diam. (mm)	Stent Length (mm)	Cath. Usable Length (cm)	Balloon Nominal Pressure (ATM) (kPa)	Balloon Rated Burst Pressure (ATM) (kPa)	Guide Catheter compatibility (F) (mm)	Recommended Guide Sheath (F) (mm)	Catheter Outer Diameter (F) (mm)	Guidewire compatibility (in) (mm)	Crossing Profile (in) (mm)							
H74939227221210	2.25	12	144	11	1117	18	1827	6	1.68	5	1.93	3.4-2.4(*)	1.20-0.85(*)	0.014	0.36	0.040	1.02
H74939227251210	2.50	12	144	11	1117	18	1827	6	1.68	5	1.93	3.4-2.4(*)	1.20-0.85(*)	0.014	0.36	0.040	1.02
H74939227301210	3.00	12	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227351210	3.50	12	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227401210	4.00	12	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227222010	2.25	20	144	11	1117	18	1827	6	1.68	5	1.93	3.4-2.4(*)	1.20-0.85(*)	0.014	0.36	0.040	1.02
H74939227252010	2.50	20	144	11	1117	18	1827	6	1.68	5	1.93	3.4-2.4(*)	1.20-0.85(*)	0.014	0.36	0.040	1.02
H74939227302010	3.00	20	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227352010	3.50	20	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227402010	4.00	20	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227222810	2.25	28	144	11	1117	18	1827	6	1.68	5	1.93	3.4-2.4(*)	1.20-0.85(*)	0.014	0.36	0.040	1.02
H74939227252810	2.50	28	144	11	1117	18	1827	6	1.68	5	1.93	3.4-2.4(*)	1.20-0.85(*)	0.014	0.36	0.040	1.02
H74939227302810	3.00	28	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227352810	3.50	28	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227402810	4.00	28	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227253810	2.50	38	144	11	1117	18	1827	6	1.68	5	1.93	3.4-2.4(*)	1.20-0.85(*)	0.014	0.36	0.040	1.02
H74939227303810	3.00	38	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.4(*)	1.20-0.85(*)	0.014	0.36	0.040	1.02
H74939227353810	3.50	38	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02
H74939227403810	4.00	38	144	11	1117	16	1620	6	1.68	5	1.93	3.4-2.7(**)	1.20-0.95(**)	0.014	0.36	0.040	1.02

(*) ≤ 1.20 mm (3.4 F) proximally and ≤ 0.85 mm (2.4 F) distally

(**) ≤ 1.20 mm (3.4 F) proximally and ≤ 0.95 mm (2.7 F) distally

PROMUS Element™ Plus DES BTK **Monorail™**

Drug Eluting Stent for Below The Knee treatment

Order Number	Stent Diam. (mm)	Stent Length (mm)	Cath. Usable Length (cm)	Balloon Nominal Pressure (ATM) (kPa)	Balloon Rated Burst Pressure (ATM) (kPa)	Guide Catheter compatibility (F) (mm)	Recommended Guide Sheath (F) (mm)	Catheter Outer Diameter (F) (mm)	Guidewire compatibility (in) (mm)	Crossing Profile (in) (mm)							
H74939215221210	2.25	12	144	11	1117	18	1827	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215251210	2.50	12	144	11	1117	18	1827	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215301210	3.00	12	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215351210	3.50	12	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215401210	4.00	12	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215222010	2.25	20	144	11	1117	18	1827	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215252010	2.50	20	144	11	1117	18	1827	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215302010	3.00	20	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215352010	3.50	20	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215402010	4.00	20	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215222810	2.25	28	144	11	1117	18	1827	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215252810	2.50	28	144	11	1117	18	1827	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215302810	3.00	28	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215352810	3.50	28	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215402810	4.00	28	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215253810	2.50	38	144	11	1117	18	1827	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215303810	3.00	38	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215353810	3.50	38	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02
H74939215403810	4.00	38	144	11	1117	16	1620	5	1.42	4	1.5	2.3-2.7(*)	0.80-0.95(*)	0.014	0.36	0.040	1.02

(*) ≤ 0.80 mm (2.3 F) proximally and ≤ 0.95 mm (2.7 F) distally

Ordering Information - 66
Peripheral Interventions

Boston
Scientific



7 Stents

7.2 Self-Expandable Stenting

WALLSTENT-Uni™ Endoprosthesis Multi-Purpose Vascular Stent with Unistep™ Plus Delivery System


WALLSTENT™ RP Iliac Endoprosthesis with Unistep™ Plus Delivery System

Epic™ Self-Expanding Nitinol Vascular Stent

Innova™ Self-expanding Stent System

7F (2.37 mm) Placehit WALLSTENT™ Endoprosthesis for Transhepatic Use

WALLGRAFT™ Endoprosthesis with Unistep™ Plus Delivery System



7.2 Self-Expandable Stenting

WALLSTENT-Uni™ Endoprosthesis

Multi-Purpose Vascular Stent with Unistep™ Plus Delivery System / Vielseitig einsetzbarer vaskulärer Stent mit Unistep™ Plus Applikationssystem / Stent vasculaire polyvalent avec système de mise en place Unistep™ Plus / Stent vascular multipropósito con sistema de liberación Unistep™ Plus

Short instrument

Order Number	Fully Open Stent		Size of WALLSTENT-Uni™ implanted in Vessel (mm)						Delivery System			
	Diam. (mm)	Length (mm)	Vessel Diam.	Stent Length	Vessel Diam.	Stent Length	Vessel Diam.	Stent Length	Catheter Size (F)	Recomm. Introd. (F)	Sheath (mm)	
M001731000	5	20	4	26	3	30	2	33	6	≤ 2.1	6	2.00
M001731020	5	40	4	54	3	63	2	68	6	≤ 2.1	6	2.00
M001731040	5	55	4	78	3	89	2	97	6	≤ 2.1	6	2.00
M001731060	5	80	4	116	3	134	2	145	6	≤ 2.1	6	2.00
M001731080	6	20	5	25	4	29	3	32	6	≤ 2.1	6	2.00
M001731100	6	45	5	52	4	60	3	65	6	≤ 2.1	6	2.00
M001731120	6	60	5	74	4	85	3	93	6	≤ 2.1	6	2.00
M001731140	6	90	5	111	4	128	3	139	6	≤ 2.1	6	2.00
M001731160	7	20	6	25	5	28	4	30	6	≤ 2.1	6	2.00
M001731180	7	40	6	50	5	57	4	63	6	≤ 2.1	6	2.00
M001731200	7	60	6	72	5	82	4	90	6	≤ 2.1	6	2.00
M001731220	7	90	6	108	5	123	4	134	6	≤ 2.1	6	2.00
M001731240	8	20	7	29	6	36	5	41	6	≤ 2.1	6	2.00
M001731260	8	40	7	49	6	56	5	61	6	≤ 2.1	6	2.00
M001731280	8	60	7	70	6	79	5	86	6	≤ 2.1	6	2.00
M001731300	8	80	7	105	6	119	5	130	6	≤ 2.1	6	2.00
M001731320	10	20	9	27	8	33	7	38	6	≤ 2.1	6	2.00
M001731340	10	42	9	48	8	54	7	58	7	≤ 2.4	7	2.33
M001731360	10	68	9	69	8	77	7	83	7	≤ 2.4	7	2.33
M001731380	10	94	9	103	8	115	7	124	7	≤ 2.4	7	2.33
M001731400	12	20	11	26	10	31	9	36	8	≤ 2.9	9	3.00
M001731420	12	40	11	47	10	51	9	55	8	≤ 2.9	9	3.00
M001731440	12	60	11	66	10	73	9	79	8	≤ 2.9	9	3.00
M001731460	12	90	11	100	10	110	9	119	8	≤ 2.9	9	3.00
M001731480	14	20	13	27	12	33	11	38	9	≤ 3.2	10	3.33
M001731490	14	40	13	46	12	50	11	54	9	≤ 3.2	10	3.33
M001731500	14	60	13	65	12	72	11	77	9	≤ 3.2	10	3.33
M001731510	14	90	13	98	12	107	11	115	9	≤ 3.2	10	3.33
M001731520	16	20	15	23	14	28	13	32	9	≤ 3.2	10	3.33
M001731530	16	40	15	45	14	49	13	52	9	≤ 3.2	10	3.33
M001731540	16	60	15	64	14	70	13	75	9	≤ 3.2	10	3.33
M001731550	16	90	15	97	14	105	13	112	9	≤ 3.2	10	3.33
M001731560	18	40	17	45	16	48	15	51	10	≤ 3.5	11	3.67
M001731570	18	60	17	64	16	69	15	73	10	≤ 3.5	11	3.67
M001731580	18	90	17	95	16	103	15	110	10	≤ 3.5	11	3.67
M001731590	20	40	19	40	18	44	17	47	10	≤ 3.5	11	3.67
M001731600	20	55	19	57	18	63	17	68	10	≤ 3.5	11	3.67
M001731610	20	80	19	86	18	94	17	102	10	≤ 3.5	11	3.67
M001731620	22	35	21	35	20	40	19	44	10	≤ 3.5	11	3.67
M001731630	22	45	21	50	20	57	19	62	10	≤ 3.5	11	3.67
M001731640	22	70	21	75	20	85	19	93	10	≤ 3.5	11	3.67
M001731650	24	35	23	35	22	39	21	43	11	≤ 3.9	12	4.00
M001731660	24	45	23	50	22	56	21	61	11	≤ 3.9	12	4.00
M001731670	24	70	23	75	22	84	21	92	11	≤ 3.9	12	4.00

Total instrument length 100 cm
Usable instrument length 75 cm
Recommended guidewire max. 0.035" (0.89 mm), 180 cm

7.2 Self-Expandable Stenting

Long instrument

Order Number	Fully Open Stent		Size of WALLSTENT-Uni™ implanted in Vessel (mm)						Delivery System			
	Diam. (mm)	Length (mm)	Vessel Diam.	Stent Length	Vessel Diam.	Stent Length	Vessel Diam.	Stent Length	Catheter Size (F)	Recomm. Introd. (F)	Sheath (mm)	
M001731010	5	20	4	26	3	30	2	33	6	≤ 2.1	6	2.00
M001731030	5	40	4	54	3	63	2	68	6	≤ 2.1	6	2.00
M001731050	5	55	4	78	3	89	2	97	6	≤ 2.1	6	2.00
M001731070	5	80	4	116	3	134	2	145	6	≤ 2.1	6	2.00
M001731090	6	20	5	25	4	29	3	32	6	≤ 2.1	6	2.00
M001731110	6	45	5	52	4	60	3	65	6	≤ 2.1	6	2.00
M001731130	6	60	5	74	4	85	3	93	6	≤ 2.1	6	2.00
M001731150	6	90	5	111	4	128	3	139	6	≤ 2.1	6	2.00
M001731170	7	20	6	25	5	28	4	30	6	≤ 2.1	6	2.00
M001731190	7	40	6	50	5	57	4	63	6	≤ 2.1	6	2.00
M001731210	7	60	6	72	5	82	4	90	6	≤ 2.1	6	2.00
M001731230	7	90	6	108	5	123	4	134	6	≤ 2.1	6	2.00
M001731250	8	20	7	29	6	36	5	41	6	≤ 2.1	6	2.00
M001731270	8	40	7	49	6	56	5	61	6	≤ 2.1	6	2.00
M001731290	8	60	7	70	6	79	5	86	6	≤ 2.1	6	2.00
M001731310	8	80	7	105	6	119	5	130	6	≤ 2.1	6	2.00
M001731330	10	20	9	27	8	33	7	38	6	≤ 2.1	6	2.00
M001731350	10	42	9	48	8	54	7	58	7	≤ 2.4	7	2.33
M001731370	10	68	9	69	8	77	7	83	7	≤ 2.4	7	2.33
M001731390	10	94	9	103	8	115	7	124	7	≤ 2.4	7	2.33
M001731410	12	20	11	26	10	31	9	36	8	≤ 2.9	9	3.00
M001731430	12	40	11	47	10	51	9	55	8	≤ 2.9	9	3.00
M001731450	12	60	11	66	10	73	9	79	8	≤ 2.9	9	3.00
M001731470	12	90	11	100	10	110	9	119	8	≤ 2.9	9	3.00

Total instrument length 160 cm
Usable instrument length 135 cm
Recommended guidewire max. 0.035" (0.89 mm), 260-300 cm

7.2 Self-Expandable Stenting

WALLSTENT™ RP

Iliac Endoprosthesis with Unistep™ Plus Delivery System / Endoprosthesis mit Unistep™ Plus-Trägersystem
Endoprothèse iliaque avec système de mise en place Unistep™ Plus / Endoprótesis ilíaca con sistema de liberación Unistep™ Plus

Short instrument

Order Number	Fully Open* Stent		Size of WALLSTENT-RP™ implanted in Vessel (mm)			
	Diam. (mm)	Length (mm)	Vessel Diam.	Stent Length	Vessel** Diam.	Stent Length
M001712000	6	24	5	33	4	41
M001712020	6	36	5	50	4	61
M001712040	6	46	5	67	4	81
M001712060	6	59	5	83	4	101
M001712080	7	23	6	31	5	38
M001712100	7	34	6	46	5	57
M001712120	7	55	6	77	5	95
M001712140	7	67	6	93	5	114
M001712160	8	20	7	29	6	36
M001712180	8	38	7	57	6	72
M001712200	8	47	7	72	6	89
M001712220	8	66	7	100	6	124
M001712240	9	18	8	27	7	34
M001712260	9	35	8	53	7	67
M001712280	9	52	8	80	7	101
M001712300	9	61	8	83	7	117
M001712320	10	20	9	27	8	33
M001712340	10	39	9	54	8	66
M001712360	10	49	9	67	8	83
M001712380	10	69	9	93	8	115

Total instrument length 100 cm
Usable instrument length 75 cm
Recommended guidewire 0.035" (0.89 mm), 240 cm
Recommended introducer sheath 6F (2.00 mm)

*As Labelled on box
**Recommended target vessel diameter

Long instrument

Order Number	Fully Open* Stent		Size of WALLSTENT-RP™ implanted in Vessel (mm)			
	Diam. (mm)	Length (mm)	Vessel Diam.	Stent Length	Vessel** Diam.	Stent Length
M001712010	6	24	5	33	4	41
M001712030	6	36	5	50	4	61
M001712050	6	46	5	67	4	81
M001712070	6	59	5	83	4	101
M001712090	7	23	6	31	5	38
M001712110	7	34	6	46	5	57
M001712130	7	55	6	77	5	95
M001712150	7	67	6	93	5	114
M001712170	8	20	7	29	6	36
M001712190	8	38	7	57	6	72
M001712210	8	47	7	72	6	89
M001712230	8	66	7	100	6	124
M001712250	9	18	8	27	7	34
M001712270	9	35	8	53	7	67
M001712290	9	52	8	80	7	101
M001712310	9	61	8	83	7	117
M001712330	10	20	9	27	8	33
M001712350	10	39	9	54	8	66
M001712370	10	49	9	67	8	83
M001712390	10	69	9	93	8	115

Total instrument length 160 cm
Usable instrument length 135 cm
Recommended guidewire 0.035" (0.89 mm), 340 cm
Recommended introducer sheath 6F (2.00 mm)

*As Labelled on box
**Recommended target vessel diameter

7.2 Self-Expandable Stenting

Epic™

Self-Expanding Nitinol Vascular Stent / Selbstexpandierender Nitinol-Stent / Stent vasculaire auto expansible en nitinol / Stent vascular de nitinol autoexpandible

Order Number		Stent Diameter (mm)	Reference Vessel Diam. (mm)	Stent Length (mm)	Sheath Size (F / mm)	Guidewire Compatibility OD (mm)/ID (mm)
Short Instrument*	Long Instrument**					
H74939054052070	H74939054052020	5	4	20	6F / 2.0	0.035/ 0.89
H74939054053070	H74939054053020	5	4	29	6F / 2.0	0.035/ 0.89
H74939054054070	H74939054054020	5	4	42	6F / 2.0	0.035/ 0.89
H74939054055070	H74939054055020	5	4	52	6F / 2.0	0.035/ 0.89
H74939054056070	H74939054056020	5	4	62	6F / 2.0	0.035/ 0.89
H74939054057070	H74939054057020	5	4	72	6F / 2.0	0.035/ 0.89
H74939054058070	H74939054058020	5	4	82	6F / 2.0	0.035/ 0.89
H74939054051070	H74939054051020	5	4	101	6F / 2.0	0.035/ 0.89
H74939054051270	H74939054051220	5	4	122	6F / 2.0	0.035/ 0.89
H74939054062070	H74939054062020	6	4.0 - 5.0	19	6F / 2.0	0.035/ 0.89
H74939054063070	H74939054063020	6	4.0 - 5.0	29	6F / 2.0	0.035/ 0.89
H74939054064070	H74939054064020	6	4.0 - 5.0	41	6F / 2.0	0.035/ 0.89
H74939054065070	H74939054065020	6	4.0 - 5.0	51	6F / 2.0	0.035/ 0.89
H74939054066070	H74939054066020	6	4.0 - 5.0	61	6F / 2.0	0.035/ 0.89
H74939054067070	H74939054067020	6	4.0 - 5.0	71	6F / 2.0	0.035/ 0.89
H74939054068070	H74939054068020	6	4.0 - 5.0	80	6F / 2.0	0.035/ 0.89
H74939054061070	H74939054061020	6	4.0 - 5.0	100	6F / 2.0	0.035/ 0.89
H74939054061270	H74939054061220	6	4.0 - 5.0	120	6F / 2.0	0.035/ 0.89
H74939054072070	H74939054072020	7	5.0 - 6.0	19	6F / 2.0	0.035/ 0.89
H74939054073070	H74939054073020	7	5.0 - 6.0	29	6F / 2.0	0.035/ 0.89
H74939054074070	H74939054074020	7	5.0 - 6.0	41	6F / 2.0	0.035/ 0.89
H74939054075070	H74939054075020	7	5.0 - 6.0	50	6F / 2.0	0.035/ 0.89
H74939054076070	H74939054076020	7	5.0 - 6.0	60	6F / 2.0	0.035/ 0.89
H74939054077070	H74939054077020	7	5.0 - 6.0	70	6F / 2.0	0.035/ 0.89
H74939054078070	H74939054078020	7	5.0 - 6.0	79	6F / 2.0	0.035/ 0.89
H74939054071070	H74939054071020	7	5.0 - 6.0	99	6F / 2.0	0.035/ 0.89
H74939054071270	H74939054071220	7	5.0 - 6.0	118	6F / 2.0	0.035/ 0.89
H74939054082070	H74939054082020	8	6.0 - 7.0	20	6F / 2.0	0.035/ 0.89
H74939054083070	H74939054083020	8	6.0 - 7.0	31	6F / 2.0	0.035/ 0.89
H74939054084070	H74939054084020	8	6.0 - 7.0	41	6F / 2.0	0.035/ 0.89
H74939054085070	H74939054085020	8	6.0 - 7.0	51	6F / 2.0	0.035/ 0.89
H74939054086070	H74939054086020	8	6.0 - 7.0	61	6F / 2.0	0.035/ 0.89
H74939054087070	H74939054087020	8	6.0 - 7.0	72	6F / 2.0	0.035/ 0.89
H74939054088070	H74939054088020	8	6.0 - 7.0	82	6F / 2.0	0.035/ 0.89
H74939054081070	H74939054081020	8	6.0 - 7.0	100	6F / 2.0	0.035/ 0.89
H74939054081270	H74939054081220	8	6.0 - 7.0	120	6F / 2.0	0.035/ 0.89
H74939054092070	H74939054092020	9	7.0 - 8.0	21	6F / 2.0	0.035/ 0.89
H74939054093070	H74939054093020	9	7.0 - 8.0	30	6F / 2.0	0.035/ 0.89
H74939054094070	H74939054094020	9	7.0 - 8.0	40	6F / 2.0	0.035/ 0.89
H74939054095070	H74939054095020	9	7.0 - 8.0	51	6F / 2.0	0.035/ 0.89
H74939054096070	H74939054096020	9	7.0 - 8.0	61	6F / 2.0	0.035/ 0.89
H74939054097070	H74939054097020	9	7.0 - 8.0	71	6F / 2.0	0.035/ 0.89
H74939054098070	H74939054098020	9	7.0 - 8.0	81	6F / 2.0	0.035/ 0.89
H74939054091070	H74939054091020	9	7.0 - 8.0	99	6F / 2.0	0.035/ 0.89
H74939054102070	H74939054102020	10	8.0 - 9.0	20	6F / 2.0	0.035/ 0.89
H74939054103070	H74939054103020	10	8.0 - 9.0	30	6F / 2.0	0.035/ 0.89
H74939054104070	H74939054104020	10	8.0 - 9.0	40	6F / 2.0	0.035/ 0.89
H74939054105070	H74939054105020	10	8.0 - 9.0	50	6F / 2.0	0.035/ 0.89

* Short instrument length 75 cm


** Long instrument length 120 cm



7.2 Self-Expandable Stenting

Order Number		Stent Diameter (mm)	Reference Vessel Diam. (mm)	Stent Length (mm)	Sheath Size (F / mm)	Guidewire Compatibility OD (mm)/ID (mm)
Short Instrument*	Long Instrument**					
H74939054106070	H74939054106020	10	8.0 - 9.0	60	6F / 2.0	0.035/ 0.89
H74939054107070	H74939054107020	10	8.0 - 9.0	71	6F / 2.0	0.035/ 0.89
H74939054108070	H74939054108020	10	8.0 - 9.0	80	6F / 2.0	0.035/ 0.89
H74939054101070	H74939054101020	10	8.0 - 9.0	98	6F / 2.0	0.035/ 0.89
H74939054123070	H74939054123020	12	9.0 - 11.0	32	6F / 2.0	0.035/ 0.89
H74939054124070	H74939054124020	12	9.0 - 11.0	42	6F / 2.0	0.035/ 0.89
H74939054125070	H74939054125020	12	9.0 - 11.0	52	6F / 2.0	0.035/ 0.89
H74939054126070	H74939054126020	12	9.0 - 11.0	61	6F / 2.0	0.035/ 0.89
H74939054143070	H74939054143020	14	11.0 - 13.0	31	6F / 2.0	0.035/ 0.89
H74939054144070	H74939054144020	14	11.0 - 13.0	41	6F / 2.0	0.035/ 0.89
H74939054145070	H74939054145020	14	11.0 - 13.0	50	6F / 2.0	0.035/ 0.89
H74939054146070	H74939054146020	14	11.0 - 13.0	60	6F / 2.0	0.035/ 0.89

* Short instrument length 75 cm
 ** Long instrument length 120 cm



7.2 Self-Expandable Stenting

Innova™

Self-expanding Stent System / Selbstexpandierende Stent-System / Système de Stent Auto-expandible
Sistema de Stent Autoexpandible

- The **Innova™** Self-expanding Stent System introduces a new level of performance in SFA and proximal popliteal stenting, with an optimal combination of **superior flexibility, high radial force** as well as significantly **reduced fracture rates**, and **excellent deliverability**.* The Innova Stent is 6F (2.1 mm) compatible and ranges from 5 mm to 8 mm in diameter and 20 mm to 200 mm in length
- Das selbstexpandierende Stent-System eröffnet eine neue Leistungsdimension im AFS- und proximalen Poplitea-Stenting. Es kombiniert **exzellente Flexibilität und hohe Radialkraft** mit signifikant reduzierten Bruchraten und **hervorragender Platzierbarkeit**.* Der Innova Stent ist 6F (2,1 mm) kompatibel und mit Durchmessern von 5 bis 8 mm sowie in Längen von 20 bis 200 mm erhältlich
- Le système de stent auto-expandible offre un nouveau niveau de performance dans le stenting de l'artère fémorale superficielle et de l'artère poplitée proximale, combinant de manière optimale **une flexibilité et une force radiale supérieures**, des taux de fracture **significativement réduits et une mise en place aisée**.* Le stent Innova est compatible 6F (2,1 mm) et existe dans un diamètre de 5 mm à 8 mm et une longueur de 20 mm à 200 mm
- El sistema de stent autoexpandible de introduce un nuevo nivel de rendimiento en los stents SFA y poplíteos proximales, con una combinación óptima de **flexibilidad superior, gran fuerza radial**, al igual que una significativa **reducción en las tasas de fracturas y una excelente capacidad de liberación**.* El stent Innova es compatible con 6F (2,1 mm) y varía entre 5 mm y 8 mm de diámetro y entre 20 mm y 200 mm de longitud

* Bench test data on file

Order Number		Stent Diameter (mm)	Reference Vessel Diam. (mm)	Stent Length (mm)	Sheath Size		Guidewire Compatibility	
Short Instrument*	Long Instrument**				(F)	(mm)	OD (mm)	ID (mm)
H749 39180050270	H749 39181050230	5	4.0	20	6	2.0	0.035	0.89
H749 39180054070	H749 39181054030	5	4.0	40	6	2.0	0.035	0.89
H749 39180056070	H749 39181056030	5	4.0	60	6	2.0	0.035	0.89
H749 39180058070	H749 39181058030	5	4.0	80	6	2.0	0.035	0.89
H749 39180051070	H749 39181051030	5	4.0	100	6	2.0	0.035	0.89
H749 39180051270	H749 39181051230	5	4.0	120	6	2.0	0.035	0.89
H749 39180051570	H749 39181051530	5	4.0	150	6	2.0	0.035	0.89
H749 39180051870	H749 39181051830	5	4.0	180	6	2.0	0.035	0.89
H749 39180052070	H749 39181052030	5	4.0	200	6	2.0	0.035	0.89
H749 39180060270	H749 39181060230	6	4.0 - 5.0	20	6	2.0	0.035	0.89
H749 39180064070	H749 39181064030	6	4.0 - 5.0	40	6	2.0	0.035	0.89
H749 39180066070	H749 39181066030	6	4.0 - 5.0	60	6	2.0	0.035	0.89
H749 39180068070	H749 39181068030	6	4.0 - 5.0	80	6	2.0	0.035	0.89
H749 39180061070	H749 39181061030	6	4.0 - 5.0	100	6	2.0	0.035	0.89
H749 39180061270	H749 39181061230	6	4.0 - 5.0	120	6	2.0	0.035	0.89
H749 39180061570	H749 39181061530	6	4.0 - 5.0	150	6	2.0	0.035	0.89
H749 39180061870	H749 39181061830	6	4.0 - 5.0	180	6	2.0	0.035	0.89
H749 39180062070	H749 39181062030	6	4.0 - 5.0	200	6	2.0	0.035	0.89
H749 39180070270	H749 39181070230	7	5.0 - 6.0	20	6	2.0	0.035	0.89
H749 39180074070	H749 39181074030	7	5.0 - 6.0	40	6	2.0	0.035	0.89
H749 39180076070	H749 39181076030	7	5.0 - 6.0	60	6	2.0	0.035	0.89
H749 39180078070	H749 39181078030	7	5.0 - 6.0	80	6	2.0	0.035	0.89
H749 39180071070	H749 39181071030	7	5.0 - 6.0	100	6	2.0	0.035	0.89
H749 39180071270	H749 39181071230	7	5.0 - 6.0	120	6	2.0	0.035	0.89
H749 39180071570	H749 39181071530	7	5.0 - 6.0	150	6	2.0	0.035	0.89
H749 39180071870	H749 39181071830	7	5.0 - 6.0	180	6	2.0	0.035	0.89
H749 39180072070	H749 39181072030	7	5.0 - 6.0	200	6	2.0	0.035	0.89
H749 39180080270	H749 39181080230	8	6.0 - 7.0	20	6	2.0	0.035	0.89
H749 39180084070	H749 39181084030	8	6.0 - 7.0	40	6	2.0	0.035	0.89
H749 39180086070	H749 39181086030	8	6.0 - 7.0	60	6	2.0	0.035	0.89
H749 39180088070	H749 39181088030	8	6.0 - 7.0	80	6	2.0	0.035	0.89
H749 39180081070	H749 39181081030	8	6.0 - 7.0	100	6	2.0	0.035	0.89
H749 39180081270	H749 39181081230	8	6.0 - 7.0	120	6	2.0	0.035	0.89
H749 39180081570	H749 39181081530	8	6.0 - 7.0	150	6	2.0	0.035	0.89
H749 39180081870	H749 39181081830	8	6.0 - 7.0	180	6	2.0	0.035	0.89
H749 39180082070	H749 39181082030	8	6.0 - 7.0	200	6	2.0	0.035	0.89

* Short instrument length 75 cm

** Long instrument length 130 cm

7.2 Self-Expandable Stenting

7F (2.37 mm) Placehit WALLSTENT™

Delivery Device with WALLSTENT™ / Applikationssystem / Dispositif de mise en place avec WALLSTENT™
Dispositivo de liberación con WALLSTENT™

Order Number	Fully Open Stent (mm)		Size of Placehit WALLSTENT implanted in Vessel (mm)					
	Diam.	Length	Diam.	Length	Diam.	Length	Diam.	Length
H965SCH650700	8.0	50	6.5	66	7.0	62	7.5	58
H965SCH650710	8.0	70	6.5	93	7.0	87	7.5	81
H965SCH650730	10.0	50	8.5	62	9.0	59	9.5	56
H965SCH650740	10.0	70	8.5	86	9.0	82	9.5	78
H965SCH650750	10.0	90	8.5	111	9.0	106	9.5	101

Total instrument length: 90.5 +/- 1.7 cm
Usable instrument length: 62.7 +/- 1.7 cm
Recommended guidewire: 0.035" (0.89mm) or 0.038" (0.97mm)
Recommended introducer sheath: 7F (2.4 mm)
System profile maximum outer diameter: 2.37 mm

WALLGRAFT™ Covered Stent

Endoprosthesis with Unistep™ Plus Delivery System / Endoprothese mit Unistep™ Plus Applikationssystem
Endoprothèse avec système de mise en place Unistep™ Plus / Endoprótesis con sistema de liberación Unistep™ Plus


Order Number	Fully Open Stent (mm)		Size of WALLGRAFT implanted in Vessel (mm)				Shaft Size		Introducer size	
	Diam.	Length	Diam.	Length	Diam.	Length	(F)	(mm)	(F)	(mm)
H965520250	6	20	-	-	5	28	9	3.1	10	3.3
H965520010	6	30	-	-	5	41	9	3.1	10	3.3
H965520020	6	50	-	-	5	69	9	3.1	10	3.3
H965520030	6	70	-	-	5	97	9	3.1	10	3.3
H965520260	7	20	5	30	6	27	9	3.1	10	3.3
H965520050	7	30	5	44	6	40	9	3.1	10	3.3
H965520060	7	50	5	74	6	67	9	3.1	10	3.3
H965520070	7	70	5	104	6	94	9	3.1	10	3.3
H965520270	8	20	6	29	7	26	9	3.1	10	3.3
H965520090	8	30	6	43	7	39	9	3.1	10	3.3
H965520100	8	50	6	72	7	66	9	3.1	10	3.3
H965520110	8	70	6	101	7	92	9	3.1	10	3.3
H965520280	9	20	7	28	8	26	10	3.4	11	3.7
H965520130	9	30	7	43	8	39	10	3.4	11	3.7
H965520140	9	50	7	71	8	65	10	3.4	11	3.7
H965520150	9	70	7	100	8	91	10	3.4	11	3.7
H965520290	10	20	8	27	9	25	10	3.4	11	3.7
H965520170	10	30	8	41	9	38	10	3.4	11	3.7
H965520180	10	50	8	68	9	63	10	3.4	11	3.7
H965520190	10	70	8	96	9	88	10	3.4	11	3.7
H965520210	12	30	9	43	10	40	11	3.7	12	4.1
H965520220	12	50	9	72	10	67	11	3.7	12	4.1
H965520230	12	70	9	100	10	94	11	3.7	12	4.1
H965520320	14	50	11	70	12	66	12	4.1	12	4.1
H965520330	14	70	11	98	12	92	12	4.1	12	4.1

Total instrument length 116 cm
Usable instrument length 90 cm
Recommended Guidewire: <= 0.035 in (0.89 mm)

7 Stents

7.3 Carotid Stenting

Adapt™ Carotid Stent and **Monorail™** Delivery System
Carotid WALLSTENT™ Monorail™ Stenting Device




Adapt™

Carotid Stent and Monorail™ Delivery System / Carotid-Stent und Monorail™ Applikationssystem
Stent carotidien et système de mise en place Monorail™ / Stent carotídeo y sistema de liberación Monorail™

- The Adapt Monorail Carotid Stent System provides the increased support and lesion coverage of closed cell design with the same flexibility of an open cell architecture*
 - Das Adapt Monorail™ Carotis-Stentsystem bietet die Flexibilität einer offenzelligen Architektur mit verbesserter Läsionsabstützung und - abdeckung des geschlossenzelligen Designs*
- Le système de stent carotidien Monorail™ ADAPT associe la flexibilité d'une architecture à cellules ouvertes au support et à la couverture artérielle qu'offre un design à cellules fermées*
 - El sistema de stent carotideo Monorail™ Adapt combina la flexibilidad de la arquitectura de celda abierta con el incremento de soporte y la cobertura de la lesion del diseño de celda cerrada*

*Data on file at Boston Scientific

Order Number	Labeled Stent		Size of ADAPT implanted in Vessel				Delivery Catheter			
	Diam. (mm) Dynamic Tapering	Length (mm)	Vessel Diam. (mm)	Stent Length (mm)	Vessel Diam. (mm)	Stent Length (mm)	Outer Diameter (F / mm)	Working Length (mm)	Guide Sheath (F / in / mm)	Guide Catheter (F / in / mm)
M001552000	4-9	21	4	24	9	21	5 / 0.068 / 1.75	135	6 / 0.087 / 2.21	8 / 0.087 / 2.21
M001552020	4-9	32	4	37	9	32	5 / 0.068 / 1.75	135	6 / 0.087 / 2.21	8 / 0.087 / 2.21
M001552040	4-9	40	4	46	9	40	5 / 0.068 / 1.75	135	6 / 0.087 / 2.21	8 / 0.087 / 2.21



Carotid WALLSTENT™ Monorail™

Carotid Stenting Device / Carotis-Stent-System / Dispositif de stent carotidien / Dispositivo stent carotídeo

Order Number	Fully Open Diameter (mm)	Size of Carotid WALLSTENT implanted in Vessel				System OD (F/mm)	Compatibility: Minimum ID	
		Vessel Diam. (mm)	Stent Length (mm)	Vessel Diam. (mm)	Stent Length (mm)		Guiding Catheter (F/in/mm)	Guiding Sheath (F/in/mm)
H965SCH647010	6	5	30	4	36	5.0/1.75	7/0.073/1.85	5/0.073/1.85
H965SCH647070	8	7	30	6	36	5.0/1.75	7/0.073/1.85	5/0.073/1.85
H965SCH647080	8	7	40	6	48	5.0/1.75	7/0.073/1.85	5/0.073/1.85
H965SCH647090	8	7	50	6	62	5.0/1.75	7/0.073/1.85	5/0.073/1.85
H965SCH647120	10	9	30	8	36	5.9/2.1	8/0.086/2.18	6/0.086/2.18
H965SCH647130	10	9	40	8	49	5.9/2.1	8/0.086/2.18	6/0.086/2.18
H965SCH647140	10	9	50	8	59	5.9/2.1	8/0.086/2.18	6/0.086/2.18

Recommended guidewire: 0.014 in (0.36mm)



8 Embolic Protection

8.1 FilterWire EZ™ Embolic Protection System

FilterWire EZ™ Embolic Protection System



8.1 FilterWire EZ™ Embolic Protection System

FilterWire EZ™

Embollic Protection System / Embolie-Protektionssystem / Système de protection embolique / Sistema de protección embólica

Order Number	Description	Crossing Profile F/inches/mm	Retrieval Profile F/inches/mm	Vessel Compatibility (mm)
H74920105 1900	FilterWire EZ 190 cm	3.2/0.042/1.07	4.3/0.056/1.42	3.5 - 5.5 mm
H74920105 3000	FilterWire EZ 300 cm	3.2/0.042/1.07	4.3/0.056/1.42	3.5 - 5.5 mm
H74950100 1500	Bent Tip Retrieval Sheath	4.0/0.053/1.33	4.0/0.053/1.33	N/A
H74922150 010	Trooper™ Patriot™ (AddWire™) Extension Wire 0.014" (0.36 mm) 145 cm	--	--	--



9 Peripheral Embolisations

9.1 Guidewires & Microcatheters

FATHOM™-16 Steerable Guidewires

Transend™ Steerable Guidewires

Renegade™ Fiber Braided Microcatheter Family

Renegade™ HI-FLO™ Fathom™ Pre-loaded Microcatheter System

FasTracker™ 325 Microcatheters



9.1 Guidewires & Microcatheters

FATHOM™-16

Steerable Guidewires / Steuerbare Führungsdrähte / Guides orientables / Guías dirigibles

Order Number	Total Length (cm)	Shaft O.D. (in / mm)	Tip Length (cm)	Shapeable Tip (cm)
M001509000	140	0.016 / 0.41	10	2
M001509010	140	0.016 / 0.41	20	2
M001509100	180	0.016 / 0.41	10	2
M001509110	180	0.016 / 0.41	20	2

Transend™

Steerable Guidewires / Steuerbare Führungsdrähte / Guides orientables / Guías dirigibles

Transend™ - 14

Order Number	Total Length (cm)	Shaft O.D. (in / mm)	Reshapeable Tip Length (cm)
M001468100	135	0.014 / 0.36	2
M001468110	165	0.014 / 0.36	2
M001468010	190	0.014 / 0.36	2

Transend™ - 18

M001468120	135	0.018 / 0.46	2
M001468130	165	0.018 / 0.46	2

Renegade™

Microcatheters / Mikrokatheter / Microcathéters / Microcatéteres

Order Number	Total Length (cm)	Distal Tip Length (cm)	Tip Design	Proximal O.D. (F / mm)	Distal O.D. (F / mm)	Distal I.D. (in / mm)
M001181270	105	20	Angled	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181280	105	30	Angled	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181330	130	20	Angled	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181340	130	30	Angled	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181390	150	20	Angled	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181400	150	30	Angled	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181250	105	20	Straight	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181260	105	30	Straight	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181310	130	20	Straight	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181320	130	30	Straight	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181370	150	20	Straight	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53
M001181380	150	30	Straight	3.0 / 1.0	2.4 / 0.80	0.021 / 0.53



9.1 Guidewires & Microcatheters

Renegade™

Fiber Braided Microcatheters / Faserverstärkte Mikrokatheter / Microcathéters tressés fibrés / Microcatéteres de fibra trenzada

Order Number	Total Length (cm)	Distal Tip Length (cm)	Proximal O.D. (F / mm)	Distal O.D. (F / mm)	Distal I.D. (in / mm)
M001182520	130	20	3.0 / 1.0	2.5 / 0.83	0.021 / 0.53
M001182510	150	10	3.0 / 1.0	2.5 / 0.83	0.021 / 0.53
M001182530	150	20	3.0 / 1.0	2.5 / 0.83	0.021 / 0.53

Renegade™ HI-FLO™

Microcatheters / Mikrokatheter / Microcathéters / Microcatéteres

Order Number	Total Length (cm)	Distal Tip Length (cm)	Proximal O.D. (F / mm)	Distal O.D. (F / mm)	Distal I.D. (in / mm)
M001182840	105	10	3.0 / 1.0	2.8 / 0.93	0.027 / 0.69
M001182850	105	20	3.0 / 1.0	2.8 / 0.93	0.027 / 0.69
M001182860	115	10	3.0 / 1.0	2.8 / 0.93	0.027 / 0.69
M001182870	115	20	3.0 / 1.0	2.8 / 0.93	0.027 / 0.69
M001182880	135	10	3.0 / 1.0	2.8 / 0.93	0.027 / 0.69
M001182890	135	20	3.0 / 1.0	2.8 / 0.93	0.027 / 0.69
M001182900	150	10	3.0 / 1.0	2.8 / 0.93	0.027 / 0.69
M001182910	150	20	3.0 / 1.0	2.8 / 0.93	0.027 / 0.69

Renegade™ HI-FLO™ Transend™ Kits

Microcatheter and Steerable Guidewire / Mikrokatheter und steuerbarer Führungsdraht

Microcathéter et guide orientable / Microcatéter y guía dirigitel

Order Number	Total Length (cm)	Distal Tip Length (cm)	Proximal O.D. (F / mm)	Distal O.D. (F / mm)	Guidewire	Guidewire Length (cm)
M001182980	105	10	3.0 / 1.0	2.8 / 0.93	Transend -18	135
M001182990	105	20	3.0 / 1.0	2.8 / 0.93	Transend -18	135
M001183000	115	10	3.0 / 1.0	2.8 / 0.93	Transend -18	135
M001183010	115	20	3.0 / 1.0	2.8 / 0.93	Transend -18	135
M001183020	135	10	3.0 / 1.0	2.8 / 0.93	Transend -18	165
M001183030	135	20	3.0 / 1.0	2.8 / 0.93	Transend -18	165

Each kit includes: Renegade HI-FLO™ microcatheter in specified size; Transend™ .018" [0.46 mm] guidewire; torque device and guidewire introducer.

Renegade™ HI-FLO™ FATHOM™ Kits

Microcatheter and Steerable Guidewire / Mikrokatheter und steuerbarer Führungsdraht

Microcathéter et guide orientable / Microcatéter y guía dirigitel

Order Number	Total Length (cm)	Distal Tip Length (cm)	Proximal O.D. (F / mm)	Distal O.D. (F / mm)	Guidewire	Guidewire Length / Tip (cm)
M001184580	105	10	3.0 / 1.0	2.8 / 0.93	FATHOM -16	140 / 10
M001184590	105	20	3.0 / 1.0	2.8 / 0.93	FATHOM -16	140 / 20
M001184600	115	10	3.0 / 1.0	2.8 / 0.93	FATHOM -16	140 / 10
M001184610	115	20	3.0 / 1.0	2.8 / 0.93	FATHOM -16	140 / 20
M001184620	135	10	3.0 / 1.0	2.8 / 0.93	FATHOM -16	180 / 10
M001184630	135	20	3.0 / 1.0	2.8 / 0.93	FATHOM -16	180 / 20
M001184640	150	10	3.0 / 1.0	2.8 / 0.93	FATHOM -16	180 / 10
M001184650	150	20	3.0 / 1.0	2.8 / 0.93	FATHOM -16	180 / 20

Each kit includes: Renegade HI-FLO™ microcatheter in specified size; FATHOM™-16 0.016" (0.41 mm) guidewire; torque device and guidewire introducer.



9.1 Guidewires & Microcatheters

Renegade™ HI-FLO™ FATHOM™

Pre-loaded Microcatheter System / Pre-loaded Microcatheter System / Pre-loaded Microcatheter System
Pre-loaded Microcatheter System

Order Number	Total Length (cm)	Distal Tip Length (cm)	Proximal O.D. (F / mm)	Distal O.D. (F / mm)	Guidewire	Guidewire Length / Tip (cm)
M001184500	140	10	3.0 / 1.0	2.8 / 0.93	FATHOM-16	105/10
M001184510	140	10	3.0 / 1.0	2.8 / 0.93	FATHOM-16	105/20
M001184520	140	10	3.0 / 1.0	2.8 / 0.93	FATHOM-16	115/10
M001184530	140	10	3.0 / 1.0	2.8 / 0.93	FATHOM-16	115/20
M001184540	180	10	3.0 / 1.0	2.8 / 0.93	FATHOM-16	135/10
M001184550	180	10	3.0 / 1.0	2.8 / 0.93	FATHOM-16	135/20
M001184560	180	10	3.0 / 1.0	2.8 / 0.93	FATHOM-16	150/10
M001184570	180	10	3.0 / 1.0	2.8 / 0.93	FATHOM-16	150/20

FasTracker™ 325

Infusion Catheters / Infusionskatheter / Cathéters de perfusion / Catéteres de infusión

Order Number	Total Length (cm)	Distal Tip Length (cm)	Proximal O.D. (F / mm)	Distal O.D. (F / mm)	Distal I.D. (in / mm)
M0011611010	105	12	3.0 / 1.0	2.7 / 0.91	0.024/ 0.61
M0011611040	135	12	3.0 / 1.0	2.7 / 0.91	0.024/ 0.61
M0011611030	140	8	3.0 / 1.0	2.7 / 0.91	0.024/ 0.61

FasTracker™ 325

Microcatheters and Preloaded Microcatheter System / Vorgeladene Mikrokatheter-Systeme / Système de microcathéter pré-chargé / Sistemas de microcatéter precargado

Order Number	Total Length (cm)	Distal Tip Length (cm)	Proximal O.D. (F / mm)	Distal O.D. (F / mm)	Guidewire	Guidewire Length (cm)
M0015501010	105	12	3.0 / 1.0	2.8 / 0.93	Taper -14 Flex Tip	155
M0015501030	125	12	3.0 / 1.0	2.8 / 0.93	Taper -14 Flex Tip	155
M0015501040	135	12	3.0 / 1.0	2.8 / 0.93	Taper -14 Flex Tip	165
M0015502030	125	12	3.0 / 1.0	2.8 / 0.93	Taper -16	155
M0015502040	135	12	3.0 / 1.0	2.8 / 0.93	Taper -16	165
M0035503010	105	12	3.0 / 1.0	2.8 / 0.93	Transend -18	135
M0035503040	135	12	3.0 / 1.0	2.8 / 0.93	Transend -18	165

Each system includes: FasTracker 325 with Hydrolene™ hydrophilic; coaxially inserted taper steerable guidewire or Transend™ .018" [0.46 mm] microguidewire; torque device; guidewire introducer.



9 Peripheral Embolisations

9.2 Embolic Agents & Occlusion Balloons

Interlock™ *Fibred IDC™ Occlusion System Family*

IDC Interlocking Detachable Coils

Complex Helical -18 Fibred Platinum Coils

Pushable Coils

Occlusion Balloons

Contour™ SE *Microspheres*

Contour™ *Embolisation Particles*



9.2 Embolic Agents & Occlusion Balloons

Interlock™ -18 Fibered IDC™ Occlusion System

Gefaserte mechanisch ablösbare Spiralen / Coils détachables à verrouillage fibré / Coils fibrados de liberación mecánica

Order Number	Max. Coil Diameter (mm)	Restrained Length (cm)	Shape
M001361500	3	6	2D Helical 
M001361510	3	12	2D Helical 
M001361520	4	8	2D Helical 
M001361530	4	15	2D Helical 
M001361540	5	8	2D Helical 
M001361550	5	15	2D Helical 
M001361560	6	10	2D Helical 
M001361570	6	20	2D Helical 
M001361580	8	20	2D Helical 
M001361590	10	20	2D Helical 
M001361600	10	30	2D Helical 
M001361610	12	20	2D Helical 
M001361620	12	30	2D Helical 
M001361630	14	20	2D Helical 
M001361640	14	30	2D Helical 
M001361740	2 / 3	2.3	VortX Diamond 
M001361750	2 / 4	4.1	VortX Diamond 
M001361760	2 / 5	5.8	VortX Diamond 
M001361770	2 / 6	8	VortX Diamond 

Interlock™ -35 Fibered IDC™ Occlusion System

Order Number	Max. Coil Diameter (mm)	Restrained Length (cm)	Shape
M001363500	3	4	2D
M001363520	4	10	2D
M001363530	4	20	2D
M001363540	6	10	2D
M001363550	6	20	2D
M001363560	6	40	2D
M001363580	8	20	2D
M001363590	8	40	2D
M001363600	10	20	2D
M001363610	10	40	2D
M001363620	12	20	2D
M001363630	12	40	2D
M001363640	15	20	2D
M001363650	15	40	2D
M001363670	18	40	2D
M001363700	4	6	Cube
M001363710	4	15	Cube
M001363720	6	10	Cube
M001363730	6	20	Cube
M001363760	8	20	Cube
M001363770	8	40	Cube
M001363790	10	25	Cube
M001363800	10	40	Cube
M001363810	15	25	Cube
M001363820	15	40	Cube
M001363830	20	40	Cube
M001363910	4	4.5	Diamond
M001363920	6	9	Diamond
M001363930	8	14	Diamond
M001363940	10	21.5	Diamond
M001363950	12	30	Diamond



IDC Interlocking Detachable Coils

Mechanisch ablösbare Spiralen / Coils détachables à verrouillage / Coils de liberación mecánica

Order Number	Max. Coil Diameter (mm)	Restrained Length (cm)	Shape
M0013612040	2	4	2D Helical 
M0013613060	3	6	2D Helical 
M0013613100	3	10	2D Helical 
M0013614040	4	4	2D Helical 
M0013614080	4	8	2D Helical 
M0013614120	4	12	2D Helical 
M0013615080	5	8	2D Helical 
M0013605150	5	15	2D Helical 
M0013606100	6	10	2D Helical 
M0013606200	6	20	2D Helical 
M0013607100	7	10	2D Helical 
M0013607200	7	20	2D Helical 
M0013608100	8	10	2D Helical 
M0013608200	8	20	2D Helical 
M0013609100	9	10	2D Helical 
M0013609200	9	20	2D Helical 
M0033601010	10	10	2D Helical 
M0013601020	10	20	2D Helical 
M0013601210	12	10	2D Helical 
M0013601220	12	20	2D Helical 
M0013601410	14	10	2D Helical 
M0013601420	14	20	2D Helical 

Complex Helical -18

Fibred Platinum Coil / Coil hélicoïdal complexe fibré en platine / Coil de Platino fibrado Complex Helical

Order Number	Product Name	Quantity per Box	Description	Max Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)
M0013120220	Complex Helical -18 (0.46 mm)	1	Complex Helical -18 4/4 mm	4	4	20
M0013120330	Complex Helical -18 (0.46 mm)	1	Complex Helical -18 6/6 mm	6	6	30
M0013120440	Complex Helical -18 (0.46 mm)	1	Complex Helical -18 7/10 mm	7	10	40
M0013120550	Complex Helical -18 (0.46 mm)	1	Complex Helical -18 8/12 mm	8	12	50
M0013120660	Complex Helical -18 (0.46 mm)	1	Complex Helical -18 10/14 mm	10	14	60
M0013120770	Complex Helical -18 (0.46 mm)	1	Complex Helical -18 11/17 mm	11	17	70
M0013120221	Complex Helical -18 (0.46 mm)	5	Complex Helical -18 4/4 mm	4	4	20
M0013120331	Complex Helical -18 (0.46 mm)	5	Complex Helical -18 6/6 mm	6	6	30
M0013120441	Complex Helical -18 (0.46 mm)	5	Complex Helical -18 7/10 mm	7	10	40
M0013120551	Complex Helical -18 (0.46 mm)	5	Complex Helical -18 8/12 mm	8	12	50
M0013120661	Complex Helical -18 (0.46 mm)	5	Complex Helical -18 10/14 mm	10	14	60
M0013120771	Complex Helical -18 (0.46 mm)	5	Complex Helical -18 11/17 mm	11	17	70

9.2 Embolic Agents & Occlusion Balloons

VortX™ Diamond -18

Pushable Coils / Vorschiebbare Spiralen / Coils à pousser / Coils empujables

Order Number	Product Name	Quantity per Box	Description	Max Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)
M0013822030	VortX Diamond -18	1	VortX Diamond -18 3/3.3 mm	3	3.3	23
M0013822040	VortX Diamond -18	1	VortX Diamond -18 4/3.7 mm	4	3.7	41
M0013822050	VortX Diamond -18	1	VortX Diamond -18 5/5.5 mm	5	5.5	58
M0013822060	VortX Diamond -18	1	VortX Diamond -18 6/6.7 mm	6	6.7	80
M0013822031	VortX Diamond -18	5	VortX Diamond -18 3/3.3 mm	3	3.3	23
M0013822041	VortX Diamond -18	5	VortX Diamond -18 4/3.7 mm	4	3.7	41
M0013822051	VortX Diamond -18	5	VortX Diamond -18 5/5.5 mm	5	5.5	58
M0013822061	VortX Diamond -18	5	VortX Diamond -18 6/6.7 mm	6	6.7	80

VortX™ -18

Pushable Coils / Vorschiebbare Spiralen / Coils à pousser / Coils empujables

Order Number	Product Name	Quantity per Box	Description	Max Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)
M0013812030	VortX -18	1	VortX -18 3/2.5 mm	3	2.5	22
M0013812040	VortX -18	1	VortX -18 4/4 mm	4	4	42
M0013812050	VortX -18	1	VortX -18 5/5.5 mm	5	5.5	60
M0013812060	VortX -18	1	VortX -18 6/6.5 mm	6	6.5	85
M0013812031	VortX -18	5	VortX -18 3/2.5 mm	3	2.5	22
M0013812041	VortX -18	5	VortX -18 4/4 mm	4	4	42
M0013812051	VortX -18	5	VortX -18 5/5.5 mm	5	5.5	60
M0013812061	VortX -18	5	VortX -18 6/6.5 mm	6	6.5	85

Straight -18

Pushable Coils / Vorschiebbare Spiralen / Coils à pousser / Coils empujables

Order Number	Product Name	Quantity per Box	Description	Max Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)
M0013120020	Straight -18	1	Straight -18 2 mm	Straight	2	2
M0013120050	Straight -18	1	Straight -18 5 mm	Straight	5	5
M0013120021	Straight -18	5	Straight -18 2 mm	Straight	2	2
M0013120051	Straight -18	5	Straight -18 5 mm	Straight	5	5

Figure 8 -18

Pushable Coils / Vorschiebbare Spiralen / Coils à pousser / Coils empujables

Order Number	Product Name	Quantity per Box	Description	Max Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)
M0013120210	Figure 8-18	1	Figure 8 -18 2/5 mm	2	5	10
M0013120211	Figure 8-18	5	Figure 8 -18 2/5 mm	2	5	10

9.2 Embolic Agents & Occlusion Balloons

Multi-Loop -18

Pushable Coils / Vorschiebbare Spiralen / Coils à pousser / Coils empujables

Order Number	Product Name	Quantity per Box	Description	Max Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)
M0013120430	Figure 8-18	1	Multi-Loop -18 4/7 mm	4	7	30
M0013120431	Figure 8-18	5	Multi-Loop -18 4/7 mm	4	7	30

2D Helical -35

Pushable Coils / Vorschiebbare Spiralen / Coils à pousser / Coils empujables

Order Number	Product Name	Quantity per Box	Description	Max Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)
M0013723020	2D -35	1	2D Helical -35 Coil 3/2.6 mm	3	2.6	20
M0013723040	2D -35	1	2D Helical -35 Coil 3/5.2 mm	3	5.2	40
M0013724030	2D -35	1	2D Helical -35 Coil 4/2.9 mm	4	2.9	30
M0013725030	2D -35	1	2D Helical -35 Coil 5/2.4 mm	5	2.4	30
M0013725050	2D -35	1	2D Helical -35 Coil 5/4.0 mm	5	4.0	50
M0013726040	2D -35	1	2D Helical -35 Coil 6/2.6 mm	6	2.6	40
M0013727040	2D -35	1	2D Helical -35 Coil 7/2.3 mm	7	2.3	40
M0013729060	2D -35	1	2D Helical -35 Coil 9/2.7 mm	9	2.7	60
M0013723021	2D -35	5	2D Helical -35 Coil 3/2.6 mm	3	2.6	20
M0013723041	2D -35	5	2D Helical -35 Coil 3/5.2 mm	3	5.2	40
M0013724031	2D -35	5	2D Helical -35 Coil 4/2.9 mm	4	2.9	30
M0013725031	2D -35	5	2D Helical -35 Coil 5/2.4 mm	5	2.4	30
M0013725051	2D -35	5	2D Helical -35 Coil 5/4.0 mm	5	4.0	50
M0013726041	2D -35	5	2D Helical -35 Coil 6/2.6 mm	6	2.6	40
M0013727041	2D -35	5	2D Helical -35 Coil 7/2.3 mm	7	2.3	40
M0013729061	2D -35	5	2D Helical -35 Coil 9/2.7 mm	9	2.7	60

VortX™ -35

Pushable Coils / Vorschiebbare Spiralen / Coils à pousser / Coils empujables

Order Number	Product Name	Quantity per Box	Description	Max Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)
M0013732040	VortX -35	1	VortX -35 4/4 mm	4	4	30
M0013733050	VortX -35	1	VortX -35 5/4.5 mm	5	4.5	35
M0013733060	VortX -35	1	VortX -35 6/5 mm	6	5.0	53
M0013733070	VortX -35	1	VortX -35 7/5.5 mm	7	5.5	67
M0013732041	VortX -35	5	VortX -35 4/4 mm	4	4	30
M0013733051	VortX -35	5	VortX -35 5/4.5 mm	5	4.5	35
M0013733061	VortX -35	5	VortX -35 6/5 mm	6	5.0	53
M0013733071	VortX -35	5	VortX -35 7/5.5 mm	7	5.5	67

Coil Pusher-16

Order Number	Total Length (mm)	Outer Diameter (in / mm)	Distal Top Outer Diameter (in / mm)
M0014012160	177	0.016 / 0.41	0.018 / 0.46

Ordering Information - 88
Peripheral Interventions

Boston
Scientific



9.2 Embolic Agents & Occlusion Balloons

Occlusion Balloons - Standard

Standard Occlusion Balloon Catheters / Standard Okklusions-Ballonkatheter / Cathéters d'occlusion à ballonnet standard
Catéteres de balón de oclusión estándar

Order Number	Catheter Size (F / mm)	Lumens	Usable Length (cm)	Recommended Guidewire (in / mm)	Inflated Balloon Diameter (mm)
M001171020	7 / 2.4	2	65	0.038 / 0.97	11.5
M001171030	7 / 2.4	2	100	0.038 / 0.97	11.5

Occlusion Balloons - Berenstein

Occlusion Balloon Catheters (Large Lumen)* / Okklusions-Ballonkatheter (großes Lumen)*
Cathéters d'occlusion à ballonnet (grande lumière)* / Catéter de balón de oclusión tipo Berenstein (lumen grande)*

Order Number	Catheter Size (F / mm)	Lumens	Usable Length (cm)	Recommended Guidewire (in / mm)	Inflated Balloon Diameter (mm)
M001173010	6 / 1.95	2	80	0.038 / 0.97	11.5

Contour™ SE Microspheres

Embolisationssphären / Particules d'embolisation / Partículas de embolización

Order Number	Size (µm)	Volume (ml)	Syringes / Box
M001761132	100-300	1.0	5
M001761172	100-300	2.0	5
M001761232	300-500	1.0	5
M001761272	300-500	2.0	5
M001761332	500-700	1.0	5
M001761372	500-700	2.0	5
M001761432	700-900	1.0	5
M001761472	700-900	2.0	5
M001761532	900-1200	1.0	5
M001761572	900-1200	2.0	5

Contour SE references will be discontinued in 2014

Contour™ Embolisation Particles

PVA Embolisationspartikel / Particules d'embolisation / Partículas de embolización

Order Number	Vials per Box	Particle Size Range (µm)
M0017600121	2	45-150
M0017600151	5	45-150
M0017600221	2	150-250
M0017600251	5	150-250
M0017600321	2	250-355
M0017600351	5	250-355
M0017600421	2	355-500
M0017600451	5	355-500
M0017600621	2	500-710
M0017600651	5	500-710
M0017600821	2	710-1000
M0017600851	5	710-1000
M0017601151	5	1000-1180

*Berenstein™ Occlusion Balloon Catheter is approved
for use with microcatheters and embolics
*Der Berenstein Okklusions Ballon ist Katheter zugelassen
für die Nutzung mit Mikrokathetern und embolische Produkten

*El catéter Balón de Oclusión Berenstein es aprobado
para su utilización con Microcatéteres y material embolizante

*Le catheter ballon occlusif Berenstein est approuvé
pour l'utilisation avec des microcatheters et des agents d'embolisation

10 Non-Vascular Interventions

10.1 All-Purpose Drainage Catheters

Flexima™ All-Purpose Drainage Catheters
Percuflex™ All-Purpose Drainage Catheters



10.1 All-Purpose Drainage Catheters

Flexima™

Hydrophilic All-Purpose Drainage Catheters* / Hydrophiler Universal-Drainage-Katheter*
Cathéters de drainage universels hydrophiles* / Catéteres hidrofílicos de drenaje multipropósito*

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001271280	APD	6 / 2.0	20	Reg
M001271290	APD	8 / 2.66	20	Reg
M001271300	APD	10 / 3.33	20	Reg
M001271360	APD	12 / 4.0	25	Reg
M001271370	APD	14 / 4.66	25	Reg
M001271260	APD	6 / 2.0	20	Firm
M001271270	APD	8 / 2.66	20	Firm
M001271250	APD	10 / 3.33	20	Firm
M001271330	APD Loop	6 / 2.0	25	Reg
M001271340	APD Loop	8 / 2.66	25	Reg
M001271350	APD Loop	10 / 3.33	25	Reg
M001271380	APD Loop	12 / 4.0	25	Reg
M001271390	APD Loop	14 / 4.66	25	Reg
M001271310	APD Loop	6 / 2.0	25	Firm
M001271320	APD Loop	8 / 2.66	25	Firm
M001271240	APD Loop	10 / 3.33	25	Firm



As shown: M001271300



As shown: M001271350

Flexima™

Hydrophilic All-Purpose Drainage Catheters* - Kits / Hydrophiler Universal-Drainage-Katheter* - Kits
Cathéters de drainage universels hydrophiles* - Kits / Kits de catéteres hidrofílicos de drenaje multipropósito

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001271960	APD Loop	8 / 2.66	25	Reg
M001271970	APD Loop	10 / 3.33	25	Reg
M001271980	APD Loop	8 / 2.66	25	Firm
M001271990	APD Loop	10 / 3.33	25	Firm

Each Kit Contains: 21 Gauge Diagnostic Needle [O.D. 0.81 mm] ; 0.018" [0.46 mm] Platinum Tip Guidewire; AccuStick™ II Introducer; 0.038" [0.97 mm] J-tip Guidewire; Dilators;

Flexima APD or APDL Catheter; Percufix(TM) Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube

Flexima™

Hydrophilic All-Purpose Drainage Catheters* - with TempTip™ Hydrophilic Tip** / Hydrophiler Universal-Drainage-Katheter* - mit hydrophiler TempTip™-Spitze** / Cathéters de drainage universels hydrophiles* - à extrémité hydrophile TempTip™ / Catéteres hidrofílicos de drenaje multipropósito* con punta hidrofílica**

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001281290	APD TT	8 / 2.66	20	Reg
M001281360	APD TT	12 / 4.0	25	Reg
M001281340	APD Loop TT	8 / 2.66	25	Reg
M001281350	APD Loop TT	10 / 3.33	25	Reg
M001281380	APD Loop TT	12 / 4.0	25	Reg
M001281390	APD Loop TT	14 / 4.66	25	Reg



As shown: M001282390

*All catheters packaged with both metal and flexible stiffening cannulae.

**Catheter distal tip is constructed of dissolvable material, leaving a large opening to facilitate drainage.

Ordering Information - 91
Peripheral Interventions

Boston
Scientific



10.1 All-Purpose Drainage Catheters

Flexima™

Hydrophilic All-Purpose Drainage Catheters* - Tight Loop with TempTip™ Hydrophilic Tip** / Hydrophiler Universal-Drainage-Katheter* - kompakter Loop mit hydrophiler TempTip™-Spitze** / Cathéters de drainage universels hydrophiles* - boucle serrée, à extrémité hydrophile TempTip™** / Catéteres hidrofilicos de drenaje multipropósito* con punta hidrofilica de bucle cerrado TempTip™**

- Small catheter diameter and tight locking loop design facilitates placement in difficult to access anatomical areas
- Le cathéter de petit diamètre et la boucle serrée facilitent la mise en place dans les zones anatomiques difficiles d'accès
- Der kleine Katheterdurchmesser und der kompakte Locking-Loop erleichtern die Platzierung auch beim Zugang in anatomisch schwierige Bereiche
- El diámetro de catéter pequeño y el diseño de bucle de bloqueo cerrado facilita la colocación en áreas anatómica de difícil acceso

Order Number	Description	F / mm	Length (cm)
M001275080	APD TT	7 / 2.33	20
M001275080	APD Loop TT	7 / 2.33	20



M001275080

Flexima™

Hydrophilic All-Purpose Drainage Catheters* - Tight Loop with TempTip™ Hydrophilic Tip** Kits / Hydrophiler-Universal-Drainage-Katheter* Kit - kompakter Loop mit hydrophiler TempTip™-Spitze** / Cathéters de drainage universels hydrophiles* - à boucle serrée, à extrémité hydrophile TempTip™** - Kits / Catéteres de drenaje hidrofilico multipropósito* con kits de punta hidrofilica de bucle cerrado TempTip™**

Order Number	Description	F / mm	Length (cm)
M001275070	APD TT K	7 / 2.33	20

Each Kit Contains: 21 Gauge Diagnostic Needle [0.018" (0.46 mm)] Platinum Tip Guidewire; AccuStick™ II Introducer; 0.038" [0.97 mm] J-tip Guidewire; Dilators; Flexima APD or APDL Catheter; Percufix™ Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube

* All catheters packaged with both metal and flexible stiffening cannulae.

** Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement. The hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.



10 Non-Vascular Interventions

10.2 vanSonnenberg Drainage Catheters

Flexima™ vanSonnenberg Drainage Catheters
Percuflex™ vanSonnenberg Drainage Catheters



10.2 vanSonnenberg Drainage Catheters

Percuflex™ vanSonnenberg VS™

Chest Drainage Catheters* / Thorax-Drainage-Katheter* / Cathéters de drainage thoracique* / Catéteres de drenaje torácico*

- Features large, skived drainage holes and large inner lumen for exceptional flow. Designed for percutaneous drainage of pleural empyemas, lung abscesses and other fluid collections
- Große schräge Drainageöffnungen und großes Innenlumen, dadurch günstige Durchflusseigenschaften. Zur perkutanen Drainage von pleuralen Empyema, Lungenabszessen und anderen Flüssigkeitsansammlungen
- Dotés de grands orifices de drainage parés et d'une grande lumière interne pour un flux exceptionnel. Conçus pour le drainage percutané des empyèmes pleuraux, des abcès pulmonaires et des accumulations d'autres fluides
- Ofrece grandes aberturas de drenaje en capas finas y gran lumen interno para un flujo excepcional. Diseñado para un drenaje percutáneo de empiemas pleurales, abscesos pulmonares y toma de muestras de otros fluidos

Order Number	Description	Diameter F / mm	Length (cm)
M001245130	VCDG Loop	12 / 4.0	25
M001245140	VCDG Loop	14 / 4.66	25



As shown: M001245130

Percuflex™ VS Sump™

Catheters* / Katheter* / Cathéters* / Catéteres*

Order Number	Description	Diameter F / mm	Length (cm)
M001203000	VSC	12 / 4.0	30
M001203010	VSC	14 / 4.66	30
M001203020	VSC	12 / 4.0	30-2
M001203030	VSC	14 / 4.66	30-2
M001203050	VSC	12 / 4.0	35
M001203060	VSC	14 / 4.66	35

Percuflex™ VS Sump™

Catheters* - Non Coated Basic Kits** / Katheter* - nicht beschichtete Basis-Kits**

Cathéters* - Kits de base sans revêtement** / Catéteres* - Kits básicos sin recubrimiento**

Order Number	Description	Diameter F / mm	Length (cm)
M001201000	VSKC	12 / 4.0	30
M001201010	VSKC	14 / 4.66	30
M001203030	VSKC Loop	14 / 4.66	35

Percuflex™ VS Sump™

Catheters* - Non Coated Catheters* with TempTip™ Hydrophilic Tip*** / Katheter* - nicht beschichtete Katheter* mit hydrophiler TempTip™-Spitze*** / Cathéters* - Cathéters sans revêtement* à extrémité hydrophile TempTip™***

Catéteres* - Catéteres sin recubrimiento* con punta hidrofílica*** TempTip™

Order Number	Description	Diameter F / mm	Length (cm)
M001213060	VSC Loop TT	12 / 4.0	35

* All catheters packaged with both metal and flexible stiffening cannulae.

** Each kit includes: Percuflex™ VS Sump™ catheter; stiffening cannula and trocar stylet; microbial airfilter; suction adaptor; .038" [0.97 mm] 15 mm "J" guidewire; 3-way stopcock and two adhesive pads.

*** Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement. the hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.



Flexima™ VS Sump™

Drainage Catheters* / Drainage-Katheter* / Cathéters de drainage* / Catéteres de drenaje*

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001273000	VSC	12 / 4.0	30	Reg
M001273010	VSC	14 / 4.66	30	Reg
M001273130	VSC	14 / 4.66	30	Firm
M001273030	VSC Loop	12 / 4.0	35	Reg
M001273040	VSC Loop	14 / 4.66	35	Reg
M001273150	VSC Loop	12 / 4.0	35	Firm
M001273160	VSC Loop	14 / 4.66	35	Firm



As shown: M001273030

*All catheters packaged with both metal and flexible stiffening cannulae.

10 Non-Vascular Interventions

10.3 Biliary Drainage Catheters

Flexima Biliary Drainage Catheter Systems

Carey-Coons[™] Biliary Endoprosthesis Kits - Soft

AmplatzAnchor[™] Percutaneous Drainage Systems with Locking Malecot



10.3 Biliary Drainage Catheters

Flexima™

Biliary Drainage Catheter Systems* / Gallengangdrainage-Kathetersysteme* / Système de cathéter de drainage biliaire*
Sistemas de catéteres de drenaje biliar*

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001271600	Biliary	10 / 3.33	35	Soft
M001271620	Biliary	14 / 4.66	35	Soft
M001271560	Biliary	8 / 2.66	35	Reg
M001271570	Biliary	10 / 3.33	35	Reg
M001271580	Biliary	12 / 4.0	35	Reg
M001271590	Biliary	14 / 4.66	35	Reg
M001271540	Biliary	8 / 2.66	35	Firm
M001271550	Biliary	10 / 3.33	35	Firm



As shown:
M001271560
M001271570
M001271550

Flexima™

Biliary Drainage Catheter Kits* / Gallengangdrainage-Katheter* - Kits / Cathéters de drainage biliaire* - Kits
Sistemas de catéteres de drenaje biliar*

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001271650	Biliary	8 / 2.66	35	Reg
M001271660	Biliary	10 / 3.33	35	Reg
M001271630	Biliary	8 / 2.66	35	Firm
M001271640	Biliary	10 / 3.33	35	Firm

Each Kit Contains: 21 Gauge Diagnostic Needle [O.D. 0.81 mm]; 0.018" [0.46 mm] Platinum Tip Guidewire; AccuStick™ II Introducer; 0.038" [0.97 mm] Guidewire; Flexima Biliary Catheter; Percufix™ Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.

Flexima™

Biliary Drainage Catheters* with TempTip™ Hydrophilic Tip** / Gallengangdrainage-Katheter* - mit hydrophiler TempTip™ -Spitze**
Cathéters de drainage biliaire* - à extrémité hydrophile TempTip™** / Catéteres de drenaje biliar* con punta hidrofílica** TempTip™

Order Number	Description	Diameter F / mm	Length (cm)	Catheter Tip Shape
M001281560	Biliary TT	8 / 2.66	35	Reg
M001281570	Biliary TT	10 / 3.33	35	Reg
M001281580	Biliary TT	12 / 4.0	35	Reg
M001281590	Biliary TT	14 / 4.66	35	Reg

* All catheters packaged with both metal and flexible stiffening cannulae.

*** Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement. the hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.



10.3 Biliary Drainage Catheters

Flexima™

Hydrophilic Biliary Drainage Catheters* - with RO Markers / Hydrophiler-Gallengangdrainage-Katheter* - mit RO-Marker
Cathéters de drainage biliaire hydrophiles* - avec marqueurs RO / Catéteres hidrofilicos de drenaje biliar* con marcadores RO

Firm Durometer for Initial Placement

Order Number	Diameter F / mm	Working Length (cm)
M001272640	8 / 2.66	35
M001272650	10 / 3.33	35

Soft Durometer for Enhanced Patient Comfort

Order Number	Diameter F / mm	Working Length (cm)
M001272660	10 / 3.33	35
M001272670	12 / 4.0	35

Regular Durometer for Enhanced Column Strength
and Patient Comfort

Order Number	Diameter F / mm	Working Length (cm)
M001272600	8 / 2.66	35
M001272610	10 / 3.33	35
M001272620	12 / 4.0	35
M001272630	14 / 4.66	35



As shown: M001272650

Flexima™

Hydrophilic Biliary Drainage Catheters* - with TempTip™ Hydrophilic Tip** and RO Markers / Hydrophiler-Gallengangdrainage-Katheter* - mit hydrophiler TempTip™-Spitze** und RO-Marker / Cathéters de drainage biliaire hydrophiles* - à extrémité hydrophile TempTip™** et avec marqueurs RO / Catéteres hidrofilicos de drenaje biliar* con puntas hidrofilicas TempTip™** y marcadores RO

TempTip™ Hydrophilic Dissolving Tip for Enhanced
Drainage and Ease of Catheter Exchange

Order Number	Diameter F / mm	Working Length (cm)
M001282600	8 / 2.66	35
M001282610	10 / 3.33	35



As shown: M001282610

Carey-Coons™

Biliary Endoprosthesis Kits - Soft / Gallengang-Endoprothese Kits - weich / Kits d'endoprothèse biliaire Souples
Kits y software de endoprótesis biliar

- Designed for totally internal biliary drainage through the biliary ducts into the duodenum.
- Für die komplett interne Gallengangdrainage über die Gallengänge im Duodenum.
- Conçus pour le drainage biliaire entièrement interne par les voies biliaires vers le duodénum.
- Diseñado para un drenaje biliar interno total a través de los conductos biliares hasta el duodeno.

Order Number	Description	Diameter F / mm	Length (cm)
M001208000	CSS	12 / 4.0	15
M001208010	CSS	14 / 4.66	15

Each Kit Contains: Dilator; Flexible Stiffening Cannula; Stabilizer; Carey-Coons Soft Stent; External Drain.



As shown: M011208010

* All catheters packaged with both metal and flexible stiffening cannulae.

** Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement. The hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.

Carey-Coons™

Biliary Endoprosthesis Kits - Soft with TempTip™ Hydrophilic Tip** / Gallengang-Endoprothese Kits - weich mit hydrophiler TempTip™-Spitze** / Kits d'endoprothèse biliaire - Souples à extrémité hydrophile TempTip™**
Kits y software de endoprótesis biliar con punta hidrofílica** TempTip™

Order Number	Description	Diameter F / mm	Length (cm)
M001218000	CSSTT	12 / 4.0	15

Each kit includes: Carey-Coons™ soft stent biliary endoprosthesis (without or with TempTip™); 12F [4.0 mm] dilator in 12F [4.0 mm] kit; 14F [4.66 mm] dilator in 14F [4.66 mm] kit; 7F [2.33 mm] initial drainage catheter and introducer system consisting of coaxial straightener and stabiliser.

Amplatz

Anchor™ Percutaneous Drainage Systems with Locking Malecot / Systeme mit Locking Malecot
Systèmes Anchor™ avec Malécot verrouillable / Sistemas con bloqueo Malecot Amplatz Anchor™

- System with Malecot for percutaneous nephrostomy, biliary, and abscess drainage
- Système avec Malécot pour la néphrostomie percutanée, le drainage biliaire et le drainage des abcès
- System mit Malecot zur perkutanen Nephrostomie, Gallengang- und Abszessdrainage
- Sistema con Malecot para nefrostomía percutánea, biliar y drenaje de abscesos

Order Number	Description	Diameter F / mm	Length (cm)
M001221000	AA	8 / 2.66	30
M001221010	AA	10 / 3.33	30
M001221020	AA	12 / 4.0	30




** Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement. The hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.

10 Non-Vascular Interventions

10.4 Urinary Drainage Catheters

Flexima™ Urinary Drainage Catheters
Percuflex™ Urinary Drainage Catheters



10.4 Urinary Drainage Catheters

Flexima™

Nephrostomy Catheters* / Nephrostomie-Katheter* - Kits* / Cathéters de néphrostomie* / Catéteres de nefrostomía*

- For percutaneous drainage of urinary collections
- Pour le drainage percutané des accumulations urinaires
- Perkutane Drainage von Urinansammlungen
- Para un drenaje percutáneo de las muestras de orina

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001271830	Neph	10 / 3.33	25	Soft
M001271840	Neph	12 / 4.0	25	Soft
M001271850	Neph	14 / 4.66	25	Soft
M001271790	Neph	8 / 2.66	25	Reg
M001271800	Neph	10 / 3.33	25	Reg
M001271810	Neph	12 / 4.0	25	Reg
M001271820	Neph	14 / 4.66	25	Reg
M001271770	Neph	8 / 2.66	25	Firm
M001271780	Neph	10 / 3.33	25	Firm



As shown: M001271800

Flexima™

Nephrostomy Catheters* - Kits / Nephrostomie-Katheter* - Kits / Cathéters de néphrostomie* - Kits / Kits de catéteres de nefrostomía*

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001271860	Neph	8 / 2.66	25	Reg
M001271870	Neph	10 / 3.33	25	Reg
M001271880	Neph	12 / 4.0	25	Reg
M001271900	Neph	8 / 2.66	25	Firm

Each kit includes: Flexima™ hydrophilic nephrostomy drainage catheter system; Accustick™ II introducer (coaxial 4F [1.33mm] dilator / 6F [2.0mm] sheath) with locking stiffening cannula; 21 gauge [0.018 inch] diagnostic needle, 6" [15.24cm]; .018" [0.46mm] extra stiff guidewire with platinum floppy tip; .038" [0.97mm] heavy duty guidewire, 145cm; Percufix™ catheter

cuff; sheet of moisture vapour permeable dressing, two cable ties; connecting tube with stopcock for drainage bag and fascial dilators (8F [2.66mm] in 8F [2.66mm] kit / 8F [2.66mm] and 10F [3.33mm] in 10F [3.33mm] kit / 8F [2.66mm], 10F [3.33mm] and 12F [4.0mm] in 12F [4.0mm] kit / 8F [2.66mm], 10F [3.33mm], 12F [4.0mm] and 14F [4.66mm] in 14F [4.66mm] kit).

Flexima™

Nephrostomy Catheters* - with TempTip™ Hydrophilic Tip** / Nephrostomie-Katheter* - mit hydrophiler TempTip™-Spitze**
Cathéters de néphrostomie* - à extrémité hydrophile TempTip™** / Catéteres de nefrostomía* con punta hidrofílica** TempTip™

- Dissolvable tip at distal end of catheter designed to create a large opening to facilitate drainage
- L'extrémité distale dissoluble du cathéter crée une grande ouverture pour faciliter le drainage
- Auflösbare Spitze am distalen Katheterende, die eine große Öffnung hinterlässt und so die Drainage erleichtert
- Punta soluble en el extremo distal del catéter diseñada para crear una gran abertura que facilite el drenaje

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001281790	Neph TT	8 / 2.66	25	Reg
M001281800	Neph TT	10 / 3.33	25	Reg
M001281810	Neph TT	12 / 4.0	25	Reg
M001281820	Neph TT	14 / 4.66	25	Reg



As shown: M001281790

* All catheters packaged with both metal and flexible stiffening cannulae.

** Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement. The hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.



10.4 Urinary Drainage Catheters

Flexima™

Nephrostomy Catheters* - with TempTip™ Hydrophilic Tip - Kits / Nephrostomie-Katheter* - mit hydrophiler TempTip™-Spitze - Kits / Cathéters de néphrostomie* - à extrémité hydrophile TempTip™ - Kits / Kits de catéteres de nefrostomía* con puntas hidrofílicas TempTip™

Order Number	Description	Diameter F / mm	Length (cm)	Durometer
M001281660	Neph TT	10 / 3.33	25	Reg

Each kit includes: Flexima™ hydrophilic nephrostomy drainage catheter system; Accustick™ II introducer (coaxial 4F [1.33mm] dilator / 6F [2.0mm] sheath) with locking stiffening cannula; 21 gauge [0.D. 0.81mm] diagnostic needle, 6" [15.24cm]; .018" [0.46mm] extra stiff guidewire with platinum floppy tip; .038" [0.97mm] heavy duty guidewire, 145cm; Percufix™ catheter cuff; sheet of moisture vapour permeable dressing, two cable ties; connecting tube with stopcock for drainage bag and fascial dilators (8F [2.66mm] in 8F [2.66mm] kit / 8F [2.66mm] and 10F [3.33mm] in 10F [3.33mm] kit / 8F [2.66mm], 10F [3.33mm] and 12F [4.0mm] in 12F [4.0mm] kit / 8F [2.66mm], 10F [3.33mm], 12F [4.0mm] and 14F [4.66mm] in 14F [4.66mm] kit).

Percuflex VTC™

Nephrostomy Catheters* / Nephrostomie-Katheter* / Cathéters de néphrostomie* / Catéteres de nefrostomía*

- A choice for cases where additional catheter rigidity is needed.
- Auswahl für Fälle, in denen eine höhere Katheterwiderstandskraft erforderlich ist.
- Pour les cas où une rigidité supplémentaire du cathéter est nécessaire.
- Una elección para casos donde es necesaria una rigidez adicional del catéter.

Non Coated Catheters

Order Number	Description	Catheter Diameter	Length (cm)	Durometer
M001225140	VTC	8F / 2.67mm	25	Firm



As shown: M001225140

Glidex Coated Catheters

Order Number	Description	Catheter Diameter	Length (cm)
M001245260	VTG	8F / 2.67mm	25
M001245270	VTG	10F / 3.33mm	25
M001245280	VTG	12F / 4.00mm	25
M001245290	VTG	14F / 4.67mm	25

Each Kit Contains: 21 Gauge Diagnostic Needle [0.D. 0.81mm] ; .018" [0.46 mm] Platinum Tip Guidewire; AccuStick II Introducer; .038" [0.97 mm] J-Tip Guidewire; Dilators; VTC Nephrostomy Catheter; Percufix Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.

Glidex Coated Kits**

Order Number	Description	Catheter Diameter
M001245300	VTG	8F / 2.67mm
M001245320	VTG	12F / 4.00mm

Each Kit Contains: 21 Gauge Diagnostic Needle; .018" Platinum Tip Guidewire; AccuStick II Introducer; .038" J-Tip Guidewire; Dilators; VTC Nephrostomy Catheter; Percufix Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.

* All catheters packaged with both metal and flexible stiffening cannulae.

** Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement. The hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.

10 Non-Vascular Interventions

10.5 Ureteral Drainage Catheters

Flexima™

Percuflex™

Amplatz Anchor™

Malecot Nephrostomy Catheters



10.5 Ureteral Drainage Catheters

FleximaTM

Ureteral Stents* / Harnleiterstents* / Stents urétéraux* / Stents uretrales*

Order Number	Description	F / mm	Length between Pigtails (cm)	Durometer
M001274000	US	8 / 2.66	20	Reg
M001274010	US	8 / 2.66	22	Reg
M001274020	US	8 / 2.66	24	Reg
M001274030	US	8 / 2.66	26	Reg
M001274040	US	8 / 2.66	28	Reg
M001274050	US	10 / 3.33	20	Reg
M001274060	US	10 / 3.33	22	Reg
M001274070	US	10 / 3.33	24	Reg
M001274080	US	10 / 3.33	26	Reg
M001274090	US	10 / 3.33	28	Reg



As shown: M001274050

FleximaTM

Ureteral Stents* - Kits / Harnleiterstents* - Kits / Stents urétéraux* - Kits / Kits de stents uretrales*

Order Number	Description	F / mm	Length between Pigtails (cm)	Durometer
M001274110	USK	6 / 2.0	22	Reg
M001274120	USK	6 / 2.0	24	Reg
M001274130	USK	6 / 2.0	26	Reg
M001274140	USK	6 / 2.0	28	Reg
M001274150	USK	8 / 2.66	20	Reg
M001274160	USK	8 / 2.66	22	Reg
M001274170	USK	8 / 2.66	24	Reg
M001274180	USK	8 / 2.66	26	Reg
M001274190	USK	8 / 2.66	28	Reg

Each kit includes: .38" [0.97mm] guidewire, stabiliser and ureteral stent.

*All stents packaged with flexible stiffening cannula.



10.5 Ureteral Drainage Catheters

Percuflex™

Non-Coated Nephroureteral Stents - with TempTip™ Hydrophilic Tip* / Nicht beschichtete nephroureterale Stents - mit hydrophiler TempTip™-Spitze* / Stents néphro-urétéraux sans revêtement - à extrémité hydrophile TempTip™*
Stents nefrouretrales sin recubrimiento con punta hidrofílica* TempTip™

Order Number	Description	Catheter size	Length (cm)
M001235490	USTT	8	20
M001235500	USTT	8	22
M001235510	USTT	8	24
M001235530	USTT	8	28
M001235570	USTT	10	20
M001235580	USTT	10	22
M001235590	USTT	10	24
M001235600	USTT	10	26
M001235610	USTT	10	28

Percuflex™

Glidex™ Coated Nephroureteral Stents / Glidex™ beschichtete nephroureterale Stents /
Stents néphro-urétéraux à revêtement Glidex™ / Stents nefrouretrales con recubrimiento Glidex™

Order Number	Description	Catheter size	Length (cm)
M001245490	USG	8	20
M001245500	USG	8	22
M001245510	USG	8	24
M001245520	USG	8	26
M001245530	USG	8	28
M001245570	USG	10	20
M001245580	USG	10	22
M001245590	USG	10	24
M001245600	USG	10	26
M001245610	USG	10	28

Percuflex™

Glidex™ Coated Nephroureteral Stents - Kits / Glidex™ beschichtete nephroureterale Stents - Kits
Stents néphro-urétéraux à revêtement Glidex™ - Kits / Kits de stents nefrouretrales con recubrimiento Glidex™

Order Number	Description	Catheter size	Length (cm)
M001245720	USGK	6	20
M001245730	USGK	6	22
M001245740	USGK	6	24
M001245750	USGK	6	26
M001245760	USGK	6	28
M001245620	USGK	8	20
M001245630	USGK	8	22
M001245640	USGK	8	24
M001245650	USGK	8	26

Each Kit Contains: 0.038" Guidewire; Stabilizer; Ureteral Stent.

* Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement. The hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.



10.5 Ureteral Drainage Catheters

Percuflex™

Nephroureteral Stent Systems (IE) Internal / External / Nephroureterale Stent-Systeme (IE) intern/ extern
Systèmes de stent néphro-urétéral (IE) Interne/ Externe / Sistemas de stents nefrouretrales internos/externos (IE)

Order Number	Description	F / mm	Length between Pigtails (cm)
M001221360	IE	8 / 2.66	22
M001221370	IE	8 / 2.66	24
M001221380	IE	8 / 2.66	26
M001221390	IE	8 / 2.66	28
M001221400	IE	10 / 3.33	22
M001221410	IE	10 / 3.33	24
M001221420	IE	10 / 3.33	26
M001221430	IE	10 / 3.33	28



As shown: M001221400

Percuflex™

Nephroureteral Stent Systems - with TempTip™ Hydrophilic Tip* / Nephroureterale Stent-Systeme - mit hydrophiler TempTip™-Spitze*
Stents néphro-urétéraux - à extrémité hydrophile TempTip™* / Sistemas de stents nefrouretrales con punta hidrofílica* TempTip™

Order Number	Description	F / mm	Length between Pigtails (cm)
M001231360	IETT	8 / 2.66	22
M001231370	IETT	8 / 2.66	24
M001231380	IETT	8 / 2.66	26
M001231390	IETT	8 / 2.66	28
M001231400	IETT	10 / 3.33	22
M001231410	IETT	10 / 3.33	24
M001231420	IETT	10 / 3.33	26
M001231430	IETT	10 / 3.33	28

Amplatz Anchor™

Percutaneous Systems with Locking Malecot / Perkutane Systeme mit Locking Malecot / Systèmes Amplatz avec Malécot
Verrouillable Sistemas con bloqueo Malecot

Order Number	Description	F / mm	Length (cm)
M001221000	AA	8 / 2.66	30
M001221010	AA	10 / 3.33	30
M001221020	AA	12 / 4.0	30



* Dissolvable material molded into the distal tip of the catheter, designed to follow a guidewire facilitating percutaneous placement.
The hydrophilic tip is designed to dissolve, leaving a large opening at the distal end of the catheter to facilitate drainage.



Malecot Nephrostomy Catheters

Nephrostomy Catheters / Nephrostomie-Katheter / Cathéters de Néphrostomie à Malécot
Catéteres de Nefrostomía Malecot

Order Number	Description	F / mm	Length (cm)
M001224100	MN	14 / 4.66	35
M001224110	MN	16 / 5.33	35
M001224120	MN	20 / 6.66	35
M001224130	MN	24 / 8.0	35



As shown: M001224110

Stenting Malecot

Nephrostomy Catheters / Nephrostomie-Katheter / Cathéters de Néphrostomie à Malécot
Catéteres de Nefrostomía Malecot para colocación de stents

Order Number	Description	F / mm	Length (cm)
M001224200	SMN	20-8 / 6.66 - 9.33	52
M001224210	SMN	24-8 / 8.0 - 9.33	52



As shown: M001224200


10 Non-Vascular Interventions

10.6 Introducers and Accessories

AccuStick™

Percufix™

*Connecting Tube for Drainage Bag
Fascial Dilators*



AccuStick™

II Introducer Systems with RO Markers - Systems / Einführschleusen-Systeme mit RO-Marker

Systèmes d'introducteur AccuStick™ II avec marqueurs RO - Systèmes / Sistemas introductores con marcadores RO - Sistemas

- For gaining initial drainage access and facilitating placement of 0.038" [0.97 mm] working guidewire alongside 0.018" [0.46 mm] guidewire
- Pour créer un accès de drainage initial et faciliter la mise en place du guide de travail de 0,038" [0,97 mm] le long du cathéter de 0,018" [0,46 mm]
- Zur ersten Anlage eines Drainagezugangs und zur Platzierung des Führungsdrahts von 0,038" (0,97 mm) neben dem Führungsdraht 0,018" (0,46 mm)
- Para obtener un acceso inicial de drenaje y facilitar la colocación de la guía de trabajo de 0,038 pulgadas [0,97 mm] a lo largo de la guía de 0,018 pulgadas [0,46 mm]

Order Number	Description
M001207020	AccuStick II with .038" [0.97 mm] J-Tip Wire
M001207030	AccuStick II with .038" [0.97 mm] Straight Wire
M001207040	AccuStick II without .038" [0.97 mm] Wire
M001207050	AccuStick II without Guidewires
M001207100	AccuStick II with .038" [0.97 mm] J-Tip and .018" [0.46 mm] Nitinol Guidewire



As shown: M001207100

Each system includes: AccuStick™ II introducer (coaxial 4F [1.33 mm] dilator / 6F [2.0 mm] sheath with RO marker) with locking stiffening cannula and diagnostic needle with stylet, 21 gauge [O.D. 0.81 mm].

AccuStick™

II Introducer Systems with RO Markers - Components / Einführschleusen-Systeme mit RO-Marker - Komponenten / Systèmes d'introducteur AccuStick™ II avec marqueurs RO - Composants / Sistemas introductores con marcadores RO - Componentes

Order Number	Description	Number per Box
M001201081	.018" [0.46 mm], 60 cm, 7.5 cm Floppy Platinum Tip Wire	5
M001206041	AccuStick II Diagnostic Needles	10

Percufix™

Catheter Cuff Kit / Kathettermanschetten-Kit / Kit de colerette pour cathéters / Kit de manguitos de catéte

Order Number	Description	Number per Box
M001207001	PCC	10

Each kit includes: Percufix™ catheter cuff; two cable ties and 6" [15.24 cm] by 8" [20.32 cm] sheet of transparent moisture vapour permeable dressing.



As shown: M001207101

Connecting Tube for Drainage Bag

Connecting Tube for Drainage Bag / Anschlussschlauch für Drainagebeutel / Tube de connexion pour poches de drainage
Tubo de conexión para la bolsa de drenaje

Order Number	Description	Number per Box
M001229101	CTube-2: connecting tube, 30 cm length with removable stopcock	5



As shown: M001229101

Fascial Dilators

Fascial Dilators / Fasziale Dilatoren / Dilateurs fasciaux / Diladores faciales

Order Number	Description	Diameter F / mm	Length (cm)
M001482501	FD	10 /3.33	20
M001482521	FD	14 /4.67	25



As shown: M001482521

11 Radiofrequency Ablation

11.1 RF 3000™ Generator

RF 3000™ Generator



RF 3000™

Generator / Generator / Générateur / Generator

Order Number	Description
M001262210	200W Radiofrequency Generator 230V
M001262200*	200W Radiofrequency Generator 115V

* Not available for sale in the European Economic Area

11 Radiofrequency Ablation

11.2 Needle Electrodes and Accessories

LeVeen™ Needle Electrode

LeVeen™ CoAccess Electrode System

LeVeen™ SuperSlim™ Needle Electrode

Soloist™ Single Needle Electrode



11.2 Needle Electrodes and Accessories

LeVeen™

Needle Electrode / Nadelelektrode / Électrode / Electrodo de aguja

Order Number	Array Diameter (cm)	Cannula Length (cm)
M001262170	5.0	25
M001262160	5.0	15
M001262310	4.0	25
M001262130	4.0	15
M001262150	3.5	25
M001262030	3.5	15
M001262020	3.5	12
M001262050	3.0	15
M001262040	3.0	12

LeVeen™ CoAccess™

Electrode System / Elektrodensystem / Système d'électrodes / Sistema de electrodos

Order Number	Array Diameter (cm)	Cannula Length (cm)
M001262240	4.0	15
M001262230	3.5	15
M001262220	3.0	15

LeVeen™ Superslim™

Needle Electrode / Nadelelektrode / Électrode / Electrodo de aguja

Order Number	Array Diameter (cm)	Cannula Length (cm)
M001262290	3.0	25
M001262280	3.0	15
M001262270	2.0	25
M001262260	2.0	15

Soloist™ Single

Needle Electrode / Single Needle Elektrode / Électrode à une aiguille / Electrodo de aguja

Order Number	Ablation Diameter (cm)	Cannula Length (cm)
M001262500	1.0	18

RF Accessories

RF Accessories / RFA-Zubehör / Accessoires RF / Accesorios RF


Order Number	Description
M001262191	Grounding Pads (Box of 50)
M001262250	LeVeen™ CoAccess™ Introducer Set



12 Miscellaneous

12.1 Thrombolysis

Katzen[™] Infusion Wire



Katzen™

Infusion Wire / Infusionsdraht / Guide d'infusion / Guía de infusión

Order Number	Description	OD (in)	OD (mm)	Overall Length (cm)	Infusion Length (cm)	Infusion Side Holes
M001461900	Katzen™ Infusion Wire	0.035	0.89	145	3	6
M001461920	Katzen™ Infusion Wire	0.035	0.89	145	6	12
M001461940	Katzen™ Infusion Wire	0.035	0.89	145	9	18
M001461960	Katzen™ Infusion Wire	0.035	0.89	145	12	24
M001461981*	Katzen™ Core Wire	0.016-0.005	0.41 - 0.13	152	-	-

Maximum Pressure 500 PSI (3448 kPa)

*Packaged 3/box



Boston Scientific International SA
Parc Val Saint Quentin - Bâtiment H
78960 Voisins-le-Bretonneux - France
+33 1 39 30 97 00

